INSTRUCTIONS FOR USE

IMPORTANT INFORMATION

Please read prior to use

Contours VPS®

Volar Plating System (Bone Plate)

Orthofix Inc.
1720 Bray Central Drive
McKinney, Texas 75069 U.S.A.
800-266-3349
469-742-8801

European Representative

Orthofix Srl.

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37012 Bussolengo (VR) Italia

Tel: 39 045 671 9000

Fax: 39 045 671 9380

English

INDICATIONS FOR USE

The Orthofix Contours VPS, Volar Plating System, is intended for the volar fixation of fractures and osteotomies involving the distal radius.

DESCRIPTION

The Contours VPS is an anatomically contoured, Delta-shaped plate intended for volar applications to the distal radius. The benefits to utilizing a volar approach for fixation of distal radius fractures and osteotomies are that the dorsal tissues are preserved and it may reduce the chances of tendon ruptures occurring.

The Orthofix Contours VPS Plate features an angled head with two rows of holes for placement of screws distally on the radius to secure bone fragments and provide stabilization. If a bone graft is required, access to the site is made possible by the triangular opening in the plate that spans from the plate head to the distal portion of the shaft. The shaft of the Plate is placed proximally on the distal radius and offers a variety of K-wire and screw placement options for secure fixation to the external surface of the bone.

The device provides the stabilization and fixation necessary in the treatment of distal radius fractures and osteotomies to allow early functional use of the hand.

CONTRAINDICATIONS

- Cases in which there is active infection local to the operative site.
- Any patient having inadequate tissue coverage over the operative site or where there is inadequate bone stock, bone quality, or anatomical definition.
- Patients with metabolic bone disorders such as severe osteopenia, osteoporosis, osteomyelitis, poorly controlled diabetes mellitus and Paget’s disease.
- Cases in which the Contours VPS Plate would cross joint spaces or open epiphyseal growth plates such as skeletally immature patients.
- Any time implant utilization would interfere with anatomical structures or physiological functioning.
- Suspected or documented metal allergy or intolerance.
- Any patient unwilling or incapable of following post-operative care instructions.

WARNINGS

- Bone plates and screws must not be reused. Reuse could result in injury or reoperation due to breakage or infection. Do not attempt to resterilize implants that come in contact with body fluids.
- MRI (Magnetic Resonance Imaging) image quality may be affected by the presence of the Contours VPS Plate and Bone Screws.
- The Volar Plating System has not been evaluated for safety and compatibility in the MR environment.
- The Volar Plating System has not been tested for heating or migration in the MR environment.
- Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been determined to affect bone healing and could potentially have an adverse effect on the bone regenerate during fracture healing.

PRECAUTIONS

- It is essential that proper operative technique be followed for implantation.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Scratches may reduce the functional strength of the component.
- Examine all Contours VPS components carefully prior to use to assure there is no damage to any of the components. If you believe a component to be faulty, damaged, or suspect DO NOT USE.
- The Contours VPS Plate and Bone Screws should not be combined with components from another manufacturer or system. Different metal types should not be used together.

POSSIBLE ADVERSE EFFECTS

- Bending or breaking of the plate or loosening, bending, or breaking of the bone screws in situ.
- A secondary surgical procedure may be necessary to replace or remove the bone plate and screws.
- Inadequate tissue surrounding the implant may cause localized irritation and/or pain.
- Scar formation possibly causing pain and/or neurological compromise around nerves.
- Non-union, delayed union, or mal-union could occur.
- Bone loss or decrease in bone density, possibly caused by stress shielding.
- Intractable pain.
- Complications associated with metal sensitivity.
- Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke, excessive bleeding, phlebitis, wound necrosis, wound infection, or damage to blood vessels.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons, or device failures that require surgical re-intervention to remove or replace the plate or bone screws. Proper implantation and patient selection will greatly affect the results.

MATERIALS

The Contours VPS Plate and Bone Screws are available in titanium alloy. Please refer to the individual product labels to determine the material for the components.

It is important that only components of similar materials are used together to eliminate the risk of galvanic corrosion.

NON-STERILE

Unless otherwise noted, Orthofix Volar Plating System components are provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Product integrity and performance are assured only if packaging is undamaged

CLEANING

Prior to use, NON-STERILE product must be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water or with compatible detergent. After cleaning, the device and/or system components should be thoroughly rinsed in sterile distilled water and dried using clean non-woven fabric.

STERILITY

The Contours VPS Plate and Bone Screws are supplied NON-STERILE and require sterilization prior to use. The recommended, validated sterilization cycle is:

Method	Cycle	Temperature	Exposure Time
Steam	Pre-Vacuum	132° - 135°C	Minimum 10 minutes
	(minimum 4 pulses)	(270° -275°F)	
Steam	Vacuum	270° F (132° C)	10 minutes

The Contours VPS Plate and Bone Screws are intended for SINGLE USE ONLY.

Validation and routine monitoring should be performed as per ANSI/AAMI ST79:2006 and A1:2008, A2:2009, “Comprehensive guide to steam sterilization and sterility assurance in health care facilities” and ANSI/AAMI ST19, “Sterilization of health care products- Biological indicators- Part 3: Biological indicators for moist heat sterilization.” Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10-6.

SERVICE

Contact Orthofix at 800-266-3349 or 469-742-8801 for return goods authorization or, if outside the U.S.A., contact your local distributor. Please note that alterations, repair, or modifications of these devices may jeopardize patient safety.

CAUTION

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

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