<p><b>CONTOURS VPS-3®</b> Volar Plating System (bone plate) <b>3</b></p>	EN	<p><b>CONTOURS VPS-3®</b> System för volarplatta (distal radiusplatta) <b>37</b></p>	SV
<p><b>CONTOURS VPS-3®</b> Sistema di placche volari (placche ossee) <b>5</b></p>	IT	<p><b>CONTOURS VPS-3®</b> Σύστημα παραλαμιαίας στερέωσης (πλάκα οστού) <b>40</b></p>	EL
<p><b>CONTOURS VPS-3®</b> Plaque palmaire (plaque osseuse) <b>9</b></p>	FR	<p><b>CONTOURS VPS-3®</b> 掌側プレートシステム (ボーンプレート) <b>44</b></p>	JP
<p><b>CONTOURS VPS-3®</b> Volares Plattensystem (Knochenplatte) <b>13</b></p>	DE	<p><b>CONTOURS VPS-3®</b> 掌側钢板系统 (骨板) <b>47</b></p>	CN
<p><b>CONTOURS VPS-3®</b> Sistema de Placa Volar (placa ósea) <b>17</b></p>	ES	<p><b>CONTOURS VPS-3®</b> Systém volárního dlahování (kostní dlaho) <b>50</b></p>	CS
<p><b>CONTOURS VPS-3®</b> Volart skinnesystem (knogleskinne) <b>21</b></p>	DA	<p><b>CONTOURS VPS-3®</b> Volar Plating System (plytka kostna) <b>53</b></p>	PL
<p><b>CONTOURS VPS-3®</b> Volaarinen levyjärjestelmä (luuhun kiinnitettävä levy) <b>24</b></p>	FI	<p><b>CONTOURS VPS-3®</b> Sistem volarne plošče (kostna plošča) <b>57</b></p>	SL
<p><b>CONTOURS VPS-3®</b> Volarplatesystem (beinplate) <b>27</b></p>	NO	<p><b>CONTOURS VPS-3®</b> Volar Plating System (골접합용 플레이트) <b>60</b></p>	KR
<p><b>CONTOURS VPS-3®</b> Volair plaatsysteem (botplaat) <b>30</b></p>	NL		
<p><b>CONTOURS VPS-3®</b> Sistema de Placa Volar (placa óssea) <b>34</b></p>	PT		

# IMPORTANT INFORMATION

Please read prior to use

RX Only

CE 0086

## INSTRUCTIONS FOR USE

Device System Name:

# CONTOURS VPS-3®

Volar Plating System (bone plate)

Orthofix Inc.  
3451 Plano Pkwy  
Lewisville, Texas 75056 U.S.A.  
800-266-3349  
214-937-3000

European Representative  
Orthofix Srl  
Via delle Nazioni 9  
37012 Bussolengo (VR) Italia  
39 045 671 9000 Fax 39 045 671 9380

**Indications for use:** The Orthofix Contours VPS-3, Volar Plating System, is intended for the volar fixation of fractures and osteotomies involving the distal radius.

**Description:** The Contours VPS-3 is an anatomically contoured, Delta-shaped plate intended for volar applications to the distal radius. The benefits to utilizing a volar approach for fixation of distal radius fractures and osteotomies are that the dorsal tissues are preserved and it may reduce the chances of tendon ruptures occurring.

The Orthofix Contours VPS-3 Plate features an angled head with two rows of holes for placement of screws distally on the radius to secure bone fragments and provide stabilization. If a bone graft is required, access to the site is made possible by the triangular opening in the plate that spans from the plate head to the distal portion of the shaft. The shaft of the Plate is placed proximally on the distal radius and offers a variety of K-wire and screw placement options for secure fixation to the external surface of the bone. The device provides the stabilization and fixation necessary in the treatment of distal radius fractures and osteotomies to allow early functional use of the hand.

### Contraindications:

- Cases in which there is active infection local to the operative site.
- Any patient having inadequate tissue coverage over the operative site or where there is inadequate bone stock, bone quality, or anatomical definition.
- Patients with metabolic bone disorders such as severe osteopenia, osteoporosis, osteomyelitis, poorly controlled diabetes mellitus and Paget's disease.
- Cases in which the Contours VPS-3 Plate would cross joint spaces or open epiphyseal growth plates such as in skeletally immature patients.
- Any time implant utilization would interfere with anatomical structures or physiological performance.

- Suspected or documented metal allergy or intolerance.
- Any patient unwilling or incapable of following post-operative care instructions.

## Warnings:

- Bone plates and screws must not be reused. Reuse could result in injury or reoperation due to breakage or infection. Do not attempt to resterilize implants that come in contact with body fluids.
- The Contours VPS-3 Plate has not been evaluated for safety and compatibility in the MR environment. Contours VPS-3 Plate has not been tested for heating or migration in the MR environment.
- Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been determined to affect bone healing and could potentially have an adverse effect on the bone regenerate during fracture healing.

## Precautions:

- It is essential that proper operative technique be followed for implantation. Please refer to the appropriate Operative Technique Guide.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Scratches may reduce the functional strength of the component.
- Examine all Contours VPS-3 components carefully prior to use to assure there is no damage to any of the components. If you believe a component to be faulty, damaged, or suspect **DO NOT USE**.
- The Contours VPS-3 Plate and Bone Screws should not be combined with components from another manufacturer or system.
- Different metal types should not be used together. Combining implants of dissimilar metals may result in galvanic corrosion.
- The Steri-Tray is for sterilization of Contours VPS-3 components **ONLY**.

## Possible Adverse Effects:

- Bending or breaking of the plate or loosening, bending, or breaking of the bone screws in situ.
- A secondary surgical procedure may be necessary to replace or remove the bone plate and screws.
- Inadequate tissue surrounding the implant may cause localized irritation and/or pain.
- Scar formation possibly causing pain and/or neurological compromise around nerves.
- Non-union, delayed union, or mal-union could occur.
- Bone loss or decrease in bone density, possibly caused by stress shielding.
- Intractable pain.
- Complications associated with metal sensitivity.
- Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke, excessive bleeding, phlebitis, wound necrosis, wound infection, or damage to blood vessels.

**Important:** A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons, or device failures that require surgical re-intervention to remove or replace the plate or bone screws. Proper implantation and patient selection will greatly affect the results.

**Materials:** The Contours VPS-3 Plate and Bone Screws are currently available in titanium alloy. Please refer to the individual product labels for component materials.

It is important that only components of similar materials are used together to eliminate the risk of galvanic corrosion.

**Cleaning:** Prior to use, NON-STERILE product should be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water or with compatible detergent. After cleaning, the device and/or system components should be thoroughly rinsed in sterile distilled water and dried using clean non-woven fabric.

**Sterility:** Unless otherwise noted, Contours VPS-3 components are provided NON-STERILE and require sterilization prior to use. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and

sterilization procedures. Product integrity and performance are assured only if packaging is undamaged.

The recommended, validated sterilization cycle achieving a sterility assurance level (SAL) of  $10^{-6}$  is:

METHOD	ITEM	CYCLE	TEMPERATURE	EXPOSURE TIME	DRYING TIME
Steam	Steri-Tray; single wrapped (FDA cleared wrap recommended)	Pre-Vacuum (minimum 4 pulses)	132°C /270°F	4 minutes	30 minutes (place a towel underneath Steri-Tray)

The Steri-Tray is for sterilization of Contours VPS-3 components **ONLY**.

The Contours VPS-3 implants are intended for **SINGLE USE ONLY**.

Validation and routine monitoring should be performed as per AAMI recommended practice AMSI/AAMI ST79:2006 and A1:2008, A2:2009, "Comprehensive guide to steam sterilization and sterility assurance in health care facilities" and ANSI/AAMI ST 19, "Sterilization of health care products- Biological indicators- Part 3: Biological indicators for moist heat sterilization".

**Service:** Contact Orthofix at 800-266-3349 or 214-937-3000 for return goods authorization or, if outside the U.S.A., contact your local distributor. Please note that alterations, repair, or modifications of these devices may jeopardize patient safety.

**Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.**

For more information contact your local representative or call 800-266-3349.

<b>RX Only</b> Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
 See Instructions for Use	 Lot Number
 Single Use Only Do Not Reuse	 Catalogue Number
 Titanium and its Alloys	 Manufacturer
 Provided Non-Sterile	 Authorised Representative
 CE Mark	

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CV-1302-PL-US © Orthofix Holdings, Inc.