IM6058 Rev.01/10.2021  (0424980)

€ <sup>0051</sup>

**IT**
**EN**
**FR**
**DE**
**ES**
**DA**
**NL**
**SV**
**PL**
**AR**
**RU**
**NO**

**IT**

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### EN

## INSTRUCTIONS FOR USE – STERILE KIT FOR PERISTALTIC PUMP

### DESCRIPTION

The following instructions are applicable to the device called “STERILE KIT FOR PERISTALTIC PUMP” manufactured by EMOTEC. This device, supplied sterile and pyrogen-free, are designed for single use and on a single patient. If the surgical protocol entails the use in connections to other devices or accessories manufactured by other companies, the user must also comply with the "Instructions for use" of the latter.

### INDICATIONS

The“STERILE KIT FOR PERISTALTIC PUMP”is a line consisting of: piercing needle with the air-filter fitted with a hydrophobic antibacterial, PVC tubing, Silicone tubing, water cap, caps for protection and spacers clip. It is used solely for irrigation of the surgical site with sterile physiological solution. The working principle is the peristaltic actuation of the portion of the silicone tube.

### CONTRAINDICATIONS

None, provided that intende use and provided information are respected.

### STORAGE

Store product at room temperature and protected from light, in its original protective packaging.

### WARNINGS AND PRECAUTIONS

- Sterility is guaranteed only if packaging in intact. Follow handling procedure for sterile packaging.
- Before removing the device, control the integrity of the pack; discard any kit that is damaged or not perfectly sealed.
- Use the device immediately after opening the package and removing the protective caps posed at the ends.
- EMOTEC srl declines all responsibility for the non-integrity of the device in the event of transport damage.
- EMOTEC srl will replace the device only when its original defect has been ascertained.
- This device must be used by qualified professionals only.
- The user is responsible for its correct use and for any changes to the following standard operating mode.
- All system connections must be tightened by hand, avoiding the use of excessive force to avoid breaking or damaging the fittings. Test by hand the degree of locking of each connection.
- Before connecting to the patient and during use, check that there are no leaks, with solution leakage and / or air bubble formation. In the event of visible leaks or defects, check the connections and, if necessary, immediately replace the device.
- THIS DEVICE IS FOR SINGLE USE ONLY, DO NOT REUSE OR RESTERILIZE.
- THE REUSE OF THE DEVICE MAY CAUSE CROSS INFECTIONS TO THE PATIENT.
- DISPOSE OF THE USED PRODUCT OBSERVING ALL THE PRECAUTIONS FOR POTENTIALLY CONTAMINATED PRODUCTS.

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