TrueLok™ Elevate

US	<u>TRUELOK™ SYSTEMS</u> <u>2</u>	ZH	<u>TRUELOK™ SYSTEMS</u> <u>101</u>
EN	<u>TRUELOK™ SYSTEMS</u> <u>10</u>	CS	<u>TRUELOK™ SYSTEMS</u> <u>109</u>
IT	<u>TRUELOK™ SYSTEMS</u> <u>17</u>	PL	<u>TRUELOK™ SYSTEMS</u> <u>116</u>
FR	<u>TRUELOK™ SYSTEMS</u> <u>24</u>	SL	<u>TRUELOK™ SYSTEMS</u> <u>123</u>
DE	<u>TRUELOK™ SYSTEMS</u> <u>31</u>	KO	<u>TRUELOK™ SYSTEMS</u> <u>130</u>
ES	<u>TRUELOK™ SYSTEMS</u> <u>38</u>	AR	<u>TRUELOK™ SYSTEMS</u> <u>137</u>
DA	<u>TRUELOK™ SYSTEMS</u> <u>45</u>	BG	<u>TRUELOK™ SYSTEMS</u> <u>145</u>
FI	<u>TRUELOK™ SYSTEMS</u> <u>52</u>	LT	<u>TRUELOK™ SYSTEMS</u> <u>152</u>
NO	<u>TRUELOK™ SYSTEMS</u> <u>59</u>	MS	<u>TRUELOK™ SYSTEMS</u> <u>159</u>
NL	<u>TRUELOK™ SYSTEMS</u> <u>66</u>	RO	<u>TRUELOK™ SYSTEMS</u> <u>166</u>
PT	<u>TRUELOK™ SYSTEMS</u> <u>73</u>	TR	<u>TRUELOK™ SYSTEMS</u> <u>173</u>
SV	<u>TRUELOK™ SYSTEMS</u> <u>80</u>	HR	<u>TRUELOK™ SYSTEMS</u> <u>180</u>
EL	<u>TRUELOK™ SYSTEMS</u> <u>87</u>	LV	<u>TRUELOK™ SYSTEMS</u> <u>187</u>
JA	<u>TRUELOK™ SYSTEMS</u> <u>94</u>		

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

0123

PQTL E D 11/24 (0433546)

TrueLok™ Elevate

US

TRUELOK™ SYSTEMS

ORTHOFIX Srl

Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PQTLE D 11/24 (0433546)

TrueLok™ Elevate

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online.

Important information - please read prior to use.

See also instruction leaflet PQSCR for implantable devices and related instruments and PQRMD for reusable medical devices.

This Instruction For Use (IFU) is for US market

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MEDICAL DEVICE INFORMATION

DESCRIPTION

The TrueLok™ Elevate System (TrueLok Elevate hereinafter) consists of a device that is connected to the bone by means of half pins. The device consists of two bars (the fixed bar is connected to the bone; the sliding bar is connected to the bone segment) and an adjustment knob that allow transverse distraction and compression movement of the bone segment. For the fixation of the fixing bar, Half Pins from D4mm to D6mm of the Orthofix portfolio can be used. Application and removal of the TrueLok Elevate can be performed with Orthofix general orthopedic instrumentation. The TrueLok™ Elevate system can be used in conjunction with TrueLok™, the TL-HEX™ TrueLok Hexapod System and the TrueLok™ EVO.

INDICATIONS FOR USE

The TrueLok Elevate is intended for treatment of non-union or pseudoarthrosis of long bones and correction of bony or soft tissue defects or deformities. The TrueLok Elevate is indicated for adult and pediatric (greater than 2 through 21 years of age) patients.

CONTRAINDICATIONS

DO NOT USE the TrueLok Elevate if a surgical candidate exhibits or is predisposed to any of the following contraindications:

- Mental or physiological conditions who are unwilling or incapable of following postoperative care instructions as it could result in a treatment failure in the intended population.

INTENDED PATIENTS

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. The TrueLok Elevate is intended to be used in adult and pediatric (greater than 2 through 21 years of age) patients.

INTENDED USERS

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal). During treatment, the TrueLok Elevate is intended to be used by patient/caregiver as the system requires to perform gradual rotation of the Knob.

NOTES FOR USE - IMPLANT REMOVAL

Once the treatment is complete, the implants must be removed. HCP should consider premature removal in case of adverse events.

DISCLAIMER

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

MATERIAL

The implants are made from implant grade stainless steel AISI 316LVM, conforming to ASTM F138 and ISO-5832-1. This material meets the biocompatibility requirements for implants.

The frame components are made from the following materials:

- AISI 303 conforming to ASTM F899/ISO 7153-1
- AISI 304 conforming to ASTM F899/ISO 7153-1
- AISI 302 conforming to ASTM F899/ISO 7153-1
- AISI 316LVM conforming to ASTM F138/ISO 5832-1
- AISI 630 conforming to ISO 7153-1/ASTM F899
- AISI 420 conforming to ISO 7153-1/ASTM F899
- Iglidur X
- EN-AW 6082 T6 conforming to UNI EN 573
- PPSU

The instruments' parts coming in contact with the patient are made from:

- AISI 303 conforming to ASTM F899/ISO 7153-1
- AISI 304 conforming to ASTM F899/ISO 7153-1
- AISI 316LVM conforming to ASTM F138/ISO 5832-1

WARNINGS

1. The template must be positioned within the medullary canal during fluoroscopy check.
2. During initial corticotomy do not cut all the sides of the bone segment.

PRECAUTIONS

1. Pay attention that the template is positioned at least 40mm from the joint to allow the positioning of bicortical half pins in the static bar.
2. Utilizing fluoroscopy, ensure that the template is positioned within the medullary canal in order to assure an unicortical transport segment.
3. When drilling angulate the drill guide slightly (approximately 3–4°) to form a trapezoidal shape of the transport segment to prevent its impingement.
4. Before half pins insertion, do not remove simultaneously all the wires that stabilize the template, but only one at a time.
5. Check with fluoroscopy the correct insertion of the half pins in the first cortex.
6. Before marking, to avoid instability of half pins inserted in the bone segment, perform provisional tightening of both nuts of the external fixator.
7. During and after loosening of the nuts on transport segment half pin fixation bolts, pay attention not to push down the external fixator to avoid half pin instability.
8. During the insertion of bi-cortical half pins, pay attention that the universal half pin fixation bolts are in contact with the external fixator. You can help holding the bolts attached to the fixator.
9. Check with fluoroscopy the correct insertion of the half pins in the second cortex.
10. Ensure to untighten the nuts of the external fixator before completing the corticotomy.
11. Before concluding the surgical operation, check that the transverse bone transport occurs correctly.
12. Before concluding the surgical operation, check that the nuts of the external fixator are securely tightened.
13. The surgeon must evaluate the integrity of the construct at follow-up visits.
14. The bone segment must be checked periodically during treatment, making any necessary adjustments to the fixation.
15. The surgeon shall instruct the patient/caregiver about the correct adjustment to be performed during the treatment.

All Orthofix devices should be used together with their corresponding Orthofix implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the TruLok Elevate when used in conjunction with devices of other manufacturers or with other Orthofix devices if not specifically indicated in the Operative Technique.

POSSIBLE ADVERSE EVENTS

- Non-union or delayed union, malunion
- Superficial infection
- Deep infection
- Loss of fixation
- Bending, breakage or migration of the device
- Additional surgery for soft tissue defects
- Reoperation to replace a component or entire frame configuration
- Bone fracture during or after treatment
- Bone loss or reduced bone density
- Damage to surrounding tissues due to surgical trauma
- Possible tension affecting the soft tissues and/or the frame during callus manipulation
- Wound healing complications
- Tissue necrosis
- Joint contracture, instability or loss of range of motion
- Arthritic changes
- Pain, discomfort or abnormal sensations due to the presence of the device
- Complex Regional Pain Syndrome
- Residual or new deformities
- Persistence or recurrence of the initial condition requiring treatment
- Premature bone callus consolidation during distraction
- Stiffness at surgery site
- Compartment syndrome
- Events caused by intrinsic risks associated with anesthesia and surgery

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP.

MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

The TruLok™ Elevate external fixator is MR unsafe. Keep away from MRI examination room.

Instrumentations, since they do not enter in an MR environment, have not been evaluated for safety in the MR environment. They have not been tested for heating or unwanted movement in the MR environment. The safety of instrumentations in the MR environment is unknown. Performing an MR exam on a person who has these medical devices may result in injury or device malfunction.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE

The TruLok Elevate allows the achievement of correction of bony or soft tissue defects or deformities.

SPECIFIC INFORMATION ON THE DEVICE

To distract the bone segment, the knob of the device shall be pulled and rotated counterclockwise (following the direction of rotation indicated by the arrow marked on the device). Every quarter of a turn the knob returns pushed in contact with the bar. To activate it again, the knob shall be pulled and rotated as previously described. During each turn (counterclockwise) the number of dots (from 1 to 4) aligned with the Orthofix logo will increase. Contrarily, to compress the bone segment, the knob of the device shall be pulled and rotated clockwise. Every quarter of a turn, the device moves the bone segment of 0.25mm from the previous position.

RISKS DUE TO THE RE-USE OF “SINGLE USE” DEVICE

Implantable Device*

The “SINGLE USE” implantable device* of Orthofix is identified through symbol “⊗” reported on the product label. After the removal from the patient, the implantable device* has to be discarded. The re-use of an implantable device* introduces contamination risks for users and patients. The re-use of an implantable device* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

Non Implantable Device

The “SINGLE USE” non implantable device of Orthofix is identified through symbol “⊗” reported on the label or indicated in the “Instructions For Use” supplied with the products. The re-use of a “SINGLE USE” non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled “FOR SINGLE USE ONLY” can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Aluminum based devices are damaged by alkaline (pH>7) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.
- Cleaning and disinfection solutions with a pH 7–10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three time with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using a clean, lint-free cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a Washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices in order to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 min
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 min at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 min;
 - d. Final rinsing with deionised water for 3 min;
 - e. Thermal disinfection at least 90°C or 194°F (max 95°C or 203°F) for 5 minutes or until A0 = 30000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110°C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.

12. Wearing protective equipment, unload the washer disinfecter when it completes the cycle.
13. If necessary, drain excessive water and dry by using clean, lint-free cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains, repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use. All functional checks and inspections described below also cover the interfaces with other instruments or components. The failure modes below may be caused by end of life of the product, improper use or improper maintenance. Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use. Specific instructions may be available for some product codes. These instructions are linked to the product code and are available in a dedicated Orthofix website. Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization, Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix, the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AAMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization.

Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

ADDITIONAL INFORMATION

INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical follow-up and of the removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical reintervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician. The distraction of the limb will be activated in home by the patient/caregiver or in clinic theatre by the HCP. Treatment activation for pediatric patients in the home environment may require the assistance of a caregiver.

The HCP shall instruct the patient or caregiver on the correct use of the device:

- Make the adjustments or get help in making the adjustments as needed
- Report if adjustment schedule cannot be met
- Report any adverse or unexpected effects (components damage)

The HCP can provide to the patient or caregiver home environment instructions available in the following document: "TL Elevate Patient's Informative Guide - TL-2406-PL".

NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

Symbol	Description	
	Medical Device	
 	Consult instructions for use or consult electronic instructions for use	Caution: Consult instructions for use for important cautionary information
	Single Use. Do not re-use. Orthofix note: discard appropriately after use (treatment) on patient	
	Do not re-sterilize	
	STERILE. Sterilized by irradiation	
	NON-STERILE	
	Double sterile barrier system	
	Unique Device Identifier	
 	Catalogue number	Batch code
	Use-by date (year-month-day)	
  0123	CE marking in conformity to applicable European Medical Devices Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	Symbol for MR unsafe. The item is known to pose hazards in all MRI environments	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	

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PQTLE D 11/24 (0433546)

TrueLok™ Elevate

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online.

Important information - please read prior to use.

See also instruction leaflet PQSCR for implantable devices and related instruments and PQRMD for reusable medical devices.

This Instruction For Use (IFU) is NOT for US market.

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MEDICAL DEVICE INFORMATION

DESCRIPTION

The TrueLok™ Elevate System (TrueLok Elevate hereinafter) consists of a device that is connected to the bone by means of half pins. The device consists of two bars (the fixed bar is connected to the bone; the sliding bar is connected to the bone segment) and an adjustment knob that allow transverse distraction and compression movement of the bone segment. For the fixation of the fixing bar, Half Pins from D4mm to D6mm of the Orthofix portfolio can be used. Application and removal of the TrueLok Elevate can be performed with Orthofix general orthopedic instrumentation. The TrueLok™ Elevate system can be used in conjunction with TrueLok™, the TL-HEX™ TrueLok Hexapod System and the TrueLok™ EVO.

INTENDED PURPOSE AND INDICATIONS

INTENDED PURPOSE

The TrueLok Elevate is intended to provide bone fixation.

INDICATIONS FOR USE

The TrueLok Elevate is intended for correction of bony or soft tissue defects or deformities.

CONTRAINDICATIONS

DO NOT USE the TrueLok Elevate if a surgical candidate exhibits or is predisposed to any of the following contraindications:

- Mental or physiological conditions who are unwilling or incapable of following postoperative care instructions as it could result in a treatment failure in the intended population.

INTENDED PATIENTS

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. The TrueLok Elevate is intended to be used in adult patients.

INTENDED USERS

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal). During treatment, the TrueLok Elevate is intended to be used by patient/caregiver as the system requires to perform gradual rotation of the Knob.

NOTES FOR USE - IMPLANT REMOVAL

Once the treatment is complete, the implants must be removed. HCP should consider premature removal in case of adverse events.

DISCLAIMER

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

MATERIAL

The implants are made from implant grade stainless steel AISI 316LVM, conforming to ASTM F138 and ISO-5832-1. This material meets the biocompatibility requirements for implants.

The frame components are made from the following materials:

- AISI 303 conforming to ASTM F899/ISO 7153-1
- AISI 304 conforming to ASTM F899/ISO 7153-1
- AISI 302 conforming to ASTM F899/ISO 7153-1
- AISI 316LVM conforming to ASTM F138/ISO 5832-1
- AISI 630 conforming to ISO 7153-1/ASTM F899
- AISI 420 conforming to ISO 7153-1/ASTM F899
- Iglidur X
- EN-AW 6082 T6 conforming to UNI EN 573
- PPSU

The instruments' parts coming in contact with the patient are made from:

- AISI 303 conforming to ASTM F899/ISO 7153-1
- AISI 304 conforming to ASTM F899/ISO 7153-1
- AISI 316LVM conforming to ASTM F138/ISO 5832-1

WARNINGS

1. The template must be positioned within the medullary canal during fluoroscopy check.
2. During initial corticotomy do not cut all the sides of the bone segment.

PRECAUTIONS

1. Pay attention that the template is positioned at least 40mm from the joint to allow the positioning of bicortical half pins in the static bar.
2. Utilizing fluoroscopy, ensure that the template is positioned within the medullary canal in order to assure an unicortical transport segment.
3. When drilling angulate the drill guide slightly (approximately 3-4°) to form a trapezoidal shape of the transport segment to prevent its impingement.
4. Before half pins insertion, do not remove simultaneously all the wires that stabilize the template, but only one at a time.
5. Check with fluoroscopy the correct insertion of the half pins in the first cortex.
6. Before marking, to avoid instability of half pins inserted in the bone segment, perform provisional tightening of both nuts of the external fixator.
7. During and after loosening of the nuts on transport segment half pin fixation bolts, pay attention not to push down the external fixator to avoid half pin instability.
8. During the insertion of bi-cortical half pins, pay attention that the universal half pin fixation bolts are in contact with the external fixator. You can help holding the bolts attached to the fixator.
9. Check with fluoroscopy the correct insertion of the half pins in the second cortex.
10. Ensure to untighten the nuts of the external fixator before completing the corticotomy.
11. Before concluding the surgical operation, check that the transverse bone transport occurs correctly.
12. Before concluding the surgical operation, check that the nuts of the external fixator are securely tightened.
13. The surgeon must evaluate the integrity of the construct at follow-up visits.
14. The bone segment must be checked periodically during treatment, making any necessary adjustments to the fixation.
15. The surgeon shall instruct the patient/caregiver about the correct adjustment to be performed during the treatment.

All Orthofix devices should be used together with their corresponding Orthofix implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the TruLok Elevate when used in conjunction with devices of other manufacturers or with other Orthofix devices if not specifically indicated in the Operative Technique.

POSSIBLE ADVERSE EVENTS

- Non-union or delayed union, malunion
- Superficial infection
- Deep infection
- Loss of fixation
- Bending, breakage or migration of the device
- Additional surgery for soft tissue defects
- Reoperation to replace a component or entire frame configuration
- Bone fracture during or after treatment
- Bone loss or reduced bone density
- Damage to surrounding tissues due to surgical trauma
- Possible tension affecting the soft tissues and/or the frame during callus manipulation
- Wound healing complications
- Tissue necrosis
- Joint contracture, instability or loss of range of motion
- Arthritic changes
- Pain, discomfort or abnormal sensations due to the presence of the device
- Complex Regional Pain Syndrome
- Residual or new deformities
- Persistence or recurrence of the initial condition requiring treatment
- Premature bone callus consolidation during distraction
- Stiffness at surgery site
- Compartment syndrome
- Events caused by intrinsic risks associated with anesthesia and surgery

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP.

MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

The TruLok™ Elevate external fixator is MR unsafe. Keep away from MRI examination room.

Instrumentations, since they do not enter in an MR environment, have not been evaluated for safety in the MR environment. They have not been tested for heating or unwanted movement in the MR environment. The safety of instrumentations in the MR environment is unknown. Performing an MR exam on a person who has these medical devices may result in injury or device malfunction.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE

The TrueLok Elevate allows the achievement of correction of bony or soft tissue defects or deformities.

SPECIFIC INFORMATION ON THE DEVICE

To distract the bone segment, the knob of the device shall be pulled and rotated counterclockwise (following the direction of rotation indicated by the arrow marked on the device). Every quarter of a turn the knob returns pushed in contact with the bar. To activate it again, the knob shall be pulled and rotated as previously described. During each turn (counterclockwise) the number of dots (from 1 to 4) aligned with the Orthofix logo will increase. Contrarily, to compress the bone segment, the knob of the device shall be pulled and rotated clockwise. Every quarter of a turn, the device moves the bone segment of 0.25mm from the previous position.

RISKS DUE TO THE RE-USE OF “SINGLE USE” DEVICE

Implantable Device*

The “SINGLE USE” implantable device* of Orthofix is identified through symbol “” reported on the product label. After the removal from the patient, the implantable device* has to be discarded.

The re-use of an implantable device* introduces contamination risks for users and patients. The re-use of an implantable device* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

Non Implantable Device

The “SINGLE USE” non implantable device of Orthofix is identified through symbol “” reported on the label or indicated in the “Instructions For Use” supplied with the products. The re-use of a “SINGLE USE” non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled “FOR SINGLE USE ONLY” can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Aluminum based devices are damaged by alkaline (pH>7) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.

5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three time with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using a clean, lint-free cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a Washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices in order to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:

- Pre-cleaning for 4 min
 - Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 min at 55°C;
 - Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 min;
 - Final rinsing with deionised water for 3 min;
 - Thermal disinfection at least 90°C or 194°F (max 95°C or 203°F) for 5 minutes or until A0 = 30000 is reached. The water used for thermal disinfection must be purified.
 - Drying at 110°C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.
- The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
- Select and start a cycle according to the recommendations of the washer manufacturer.
 - On completion of the cycle, ensure that all stages and parameters have been achieved.
 - Wearing protective equipment, unload the washer disinfectant when it completes the cycle.
 - If necessary, drain excessive water and dry by using clean, lint-free cloth.
 - Visually inspect each device for remaining soil and dryness. If soil remains, repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use. All functional checks and inspections described below also cover the interfaces with other instruments or components. The failure modes below may be caused by end of life of the product, improper use or improper maintenance. Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use. Specific instructions may be available for some product codes. These instructions are linked to the product code and are available in a dedicated Orthofix website. Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization, Orthofix recommends using one of the following packaging systems:

- Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used.
- Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix, the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

ADDITIONAL INFORMATION

INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical follow-up and of the removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical reintervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician.

The HCP shall instruct the patient or caregiver on the correct use of the device:

- Make the adjustments or get help in making the adjustments as needed
- Report if adjustment schedule cannot be met
- Report any adverse or unexpected effects (components damage)

NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

Symbol	Description	
	Medical Device	
 	Consult instructions for use or consult electronic instructions for use	Caution: Consult instructions for use for important cautionary information
	Single Use. Do not re-use. Orthofix note: discard appropriately after use (treatment) on patient	
	Do not re-sterilize	
	STERILE. Sterilized by irradiation	
	NON-STERILE	
	Double sterile barrier system	
	Unique Device Identifier	
 	Catalogue number	Batch code
	Use-by date (year-month-day)	
 	CE marking in conformity to applicable European Medical Devices Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	Symbol for MR unsafe. The item is known to pose hazards in all MRI environments	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	

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