Ankle

Hindfoot Nailing System™

EN	ANKLE HINDFOOT NAILING (AHN) SYSTEM 2	JA	足関節後足部ネイル (AHN) システム 128
IT	SISTEMA AHN (ANKLE HINDFOOT NAILING) 12	ZH	踝后足髓内钉系统 139
FR	SYSTÈME D'ENCLOUAGE CHEVILLE/ARRIÈRE-PIED (ANKLE HINDFOOT NAILING, AHN) 23	CS	SYSTÉM ANKLE HINDFOOT NAILING (AHN) 149
DE	SPRUNGGELENK-ARTHRODESE-NAGELSYSTEM 34	PL	SYSTEM DO GWOŹDZIOWANIA DO STAWU SKOKOWEGO I TYŁOSTOPIA (AHN) 159
ES	SISTEMA DE ENCLAVAMIENTO DE TOBILLO Y RETROPIÉ 45	SL	ŽEBELJNI SISTEM ZADNJEGA STOPALA (AHN SISTEM) 170
DA	ANKLE HINDFOOT NAILING-SYSTEM (AHN) 56	KO	ANKLE HINDFOOT NAILING(AHN) SYSTEM 181
FI	NILKKA-KANTAPÄÄ NAULAUS (AHN) -JÄRJESTELMÄ 66	AR	جهاز تسمير الكاحل والعقب (AHN) 191
NO	ANKLE HINDFOOT NAILING-SYSTEM (AHN) 76	BG	СИСТЕМА ЗА ЗАКОВАВАНЕ НА ГЛЕЗЕНА - ЗАДНАТА ЧАСТ НА ХОДИЛОТО (ЗГЗЧХ) 202
NL	AHN-SYSTEEM (ANKLE HINDFOOT NAILING) 86	LT	UŽPAKALINĖS KULKŠNIES VINIŲ (UKV) SISTEMA 214
PT	SISTEMA ANKLE HINDFOOT NAILING (AHN) 97	MS	SISTEM PEMAKUAN KAKI PADA BUKU LALI (AHN) 224
SV	SPIKSYSTEM FÖR FOTLEDEN (AHN) 107	RO	SISTEM DE OSTEOSINTEZĂ A GLEZNEI ȘI RETROPICIORULUI (AHN) 235
EL	ΣΥΣΤΗΜΑ ΗΛΩΣΗΣ ΑΣΤΡΑΓΑΛΟΥ-ΟΠΙΣΘΙΟΥ ΠΟΔΑ (AHN) 117	TR	ANKLE HINDFOOT NAILING (AHN) SYSTEM 245

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

0123

PQAHN D 02/22 (0423779-0424548)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online.

Important information - please read prior to use

See also instruction leaflet PQRMD for reusable medical devices

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RX ONLY

INDICATIONS

The Orthofix Ankle Hindfoot Nailing (AHN) System is intended to facilitate tibiotalarcanal arthrodesis (fusion). Specific indications include:

1. Avascular necrosis of the talus
2. Failed total ankle arthroplasty
3. Trauma (malunited tibial pilon fracture)
4. Severe deformity or instability as a result of talipes equinovarus, cerebral vascular accident, paralysis or other neuromuscular disease
5. Revision ankle arthrodesis
6. Neuroarthropathy
7. Rheumatoid arthritis
8. Osteoarthritis
9. Pseudoarthrosis
10. Post-traumatic arthrosis
11. Previously infected arthrosis
12. Charcot foot
13. Severe endstage degenerative arthritis
14. Severe defects after tumor resection
15. Pantalar arthrodesis

DESCRIPTION

The Orthofix Ankle Hindfoot Nailing System comprises a variety of implantable intramedullary nails with respective locking end caps and locking screws manufactured from titanium alloy (Ti6Al4V) which are intended to facilitate tibiotalarcanal arthrodesis (fusion). The system also includes instrumentation required for implantation and explantation of the implants.

CONTRAINDICATIONS

The Orthofix Ankle Hindfoot Nailing (AHN) System is not designed or sold for any use except as indicated. Use of the Orthofix Ankle Hindfoot Nailing (AHN) System is contraindicated in the following situations:

1. Active or latent infection in the affected area.
2. General medical conditions including: impaired blood supply, pulmonary insufficiency (i.e. ARDS, fat embolism) and insufficient quantity or quality of bone.
3. Patients who are unwilling or incapable of following post-operative care instructions.
4. Suspected or documented metal allergy or intolerance.
5. Severe longitudinal deformity
6. Insufficient plantar heel pad
7. Situations where an isolated ankle or subtalar fusion can be performed
8. Dysvascular limb

GENERAL WARNINGS

These implants are not intended to replace normal healthy bone. Loads produced by weight bearing and activity levels will dictate the longevity of the implant. The patient should understand that stress on an implant can be produced without actual weight bearing. In the absence of solid bony union, the weight of the limb alone, muscular forces associated with moving a limb, or repeated stresses of apparent relatively small magnitude, can result in failure of the implant. Therefore, the patient should follow the postoperative instructions given by the surgeon.

The product is intended for professional use only. Surgeons who supervise the use of the product must have full awareness of orthopedic fixation procedures and should have received adequate training in the use of the product. Prior to surgery, surgeons should be familiar with the devices, instruments and surgical procedure, including the application and removal. Detailed operative technique guidance is available on request; please contact Orthofix or your local distributor.

1. Processing and re-processing of instrumentation should be performed according to the dedicated leaflet PQAHR.
2. Ensure that all components needed for the operation are available and fully functional in the operating theatre before surgery. Perform targeting and functional control before the insertion of the nail as well to avoid nail damage during drilling operation.
3. Implants, nail end caps, locking screws as well as certain parts of the instrumentation (where indicated on the label), are for single use only and must never be reused. If any implant has come into contact with any body fluid it should be considered to have been used.

4. Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been shown to affect bone healing and could potentially have an adverse effect on bone repair during fracture healing.
5. Components of this device are not approved for screw attachment to the posterior elements of the cervical, thoracic, or lumbar spine.
6. Implants can loosen, fracture, corrode, migrate, or cause pain. If an implant remains implanted after complete healing, the implant may cause stress shielding, which may increase the risk of refracture in an active patient. The surgeon should weigh the risks versus benefits when deciding whether to remove the implant. Adequate postoperative management to avoid refracture should follow implant removal.

GENERAL PRECAUTIONS

1. Implants are made of titanium alloy. Do not use with implants of dissimilar metals, since this may cause an electrolytic reaction.
2. Do not use components of the Orthofix Ankle Hindfoot Nailing (AHN) System in conjunction with products of other manufacturers, unless otherwise specified.
3. Care should be used in the handling and storage of the components. The implants should not be scratched, notched or otherwise damaged as these may reduce the functional strength of the component. Contours or bending of an implant should be avoided.
4. All components should be carefully examined prior to use. Product integrity, sterility (in the case of sterile products) and performance are assured only if the packaging is undamaged. DO NOT USE if packaging is compromised or if a component is believed to be faulty, damaged or suspect.
5. It is suggested to verify all the system assembly prior to the implantation.
6. Correct selection of the implant size is extremely important. The appropriate size should be selected for the patient. Failure to use the correct components or their improper positioning may result in loosening, bending, cracking, or fracture of the device or bone or both. The X-Ray Ruler should be used to determine the appropriate nail diameter and length in order to ensure the nail does not protrude from the base of the calcaneus and cause pain.
7. Care should be taken to attach all components correctly and avoid damage to equipment. Ensure that any locations marks correspond, and lock them securely. It is important to check that the required interface between instruments remain firmly fixed during the procedure.
8. Protect the plantar neurovascular structures, both in the dissection and during the procedure, as these structures are at risk. Joint exposure and preparation is required.
9. Targeting devices: Alignment of the Talar Screw should be checked to confirm the position of the internal mechanism. For alignment check ensure that all holes in the nail are properly targeted by the appropriate Targeting Assembly. If both PA and LM screws will be used in calcaneus and LM screws will be used in talus, swing the Distal Targeting Assembly to both PA and LM positions and ensure proper targeting. Ensure all holes are targeted properly before drilling the hole. As with any targeting device, success depends on many factors and cannot be guaranteed in all cases. A free hand technique can be used alternatively. Use fluoroscopy to ensure that drill has properly targeted the nail.
10. The length of locking screws should be determined by using the Drill Guide pressed firmly against the cortex and reading the appropriate screw length directly from the calibrations on the 4.3mm Calibrated Drill, or by using the Long Depth Gauge through the Screw Guide.
11. At least one locking screw should be inserted in each bone segment (tibia, talus, calcaneus). Provide adequate countersink for threaded head screws and low profile screws if inserted in the calcaneus to avoid soft tissue damage.
12. When locking screws are inserted, care should be taken not to penetrate the articulation.
13. Surgeon must use care during oblique screw placement as the oblique mechanism allows up to 10 degrees of angulation. Surgeon must be careful to not hit other calcaneal or talar screws when using this oblique screw.
14. It is important that the locking screws are not inserted close to the fracture line, as effective screw fixation may be compromised, causing implant or fixation failure. In very distal fractures, it is important to ensure that the more proximal of the distal locking screws is at least 1 cm distal to the fracture line.
15. There should never be an empty locking hole between locking screws and the fracture, nor at the fracture site. Care should be taken to use locking screws of correct diameter, length and type.
16. Threaded head screws are intended to be used in the calcaneus.
17. Nails are inserted with or without reaming, depending on nail diameter (1-1.5mm more than nail diameter), patient, fracture type and bone diameter and quality. This decision is up to the operating surgeon.
18. Countersink the nail 5.0mm-10mm to allow for desired compression. This also prevents painful prominence in plantar soft tissue. NOTE: To avoid excessive penetration of the drill bit and subsequent loss of "working" threads, check that the colored mark on the drill bit 6.1 is not hidden by the soft tissue sleeve while countersinking for screw head. Use fluoroscopy to ensure that drill does not hit the ML screws already inserted into the nail.
19. Use of end cap is mandatory to avoid bone ingrowth and to lock calcaneal screws. The length of the end caps should be chosen carefully, in order to avoid protrusion of the nail from the calcaneus.
20. Hammering of the nails may only be performed through the nail attachment rod or with specifically designed insertion tools. Hammering should always be gentle; vigorous hammering should never be necessary. The surgeon should never persist with hammering if the implant is not advancing, but should review the situation, and consider further reaming or, if this is not possible, a smaller implant.
21. Cannulated instruments should be inspected prior to use to confirm that the lumen is free from obstruction. The correctly sized K-wire should be passed through it, to check that it slides easily, before and after each use.
22. It is recommended that all wires that are used to guide cannulated instruments or implants should be single patient use only. Prior to any use, they should be checked and should be discarded if found to be scratched or bent.
23. During the introduction of any instrument or implant over a wire, the wire tip should be screened, using fluoroscopy, as continuously as possible to exclude inadvertently driving the wire further than intended. This is particularly important when the wires are pointing towards potentially dangerous locations. The surgeon should check that there is no bony or other debris built up on the wire or inside the instrument or implant which might cause it to bind on the wire and push it forward.
24. After the procedure, check that fracture reduction and the position of all implants are correct using the image intensifier.
25. Tibiotalar compression, up to 7mm, may be achieved with the internal compression mechanism and should be monitored with fluoroscopy. When necessary, external compression can be applied. Once the Heel Cup is in contact with soft tissue, care should be taken to monitor (using fluoroscopy) the amount of compression applied to avoid collapse of the joint and avoid soft tissue impingement.
26. General advice on weightbearing (unless otherwise specified): full weightbearing can be commenced when there is radiological evidence of bridging callus. Where there is good contact between the two intact segments of bone, so that load sharing can be anticipated, early weightbearing as tolerated should be encouraged. When the bone is comminuted or there is bone loss, so that load sharing will not be possible until callus is formed, weightbearing should be partial only initially. The exact amount of load to be carried depends on the size of implant inserted and the stature of the patient. Hip and knee mobility should always be encouraged within the pain limits. However, the best clinical results are obtained by encouraging early mobility and full weightbearing as tolerated as early as possible and according to the patient's condition.
27. Careful monitoring of the progress of healing must be undertaken in all patients. If callus is slow to develop, other measures may be required to promote its formation, such as dynamisation of the implant, a bone graft, or replacing the implant.
28. Additional equipment may be required for fixation application and removal such as soft tissue retractors, flexible reaming set, cannulated drills, etc.
29. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.

30. Before inserting the Slotted Hammer Adaptor, for proper nail removal unlock the talar screw by inserting and turning the Locking Driver (177301) counterclockwise and, if necessary, release the internal compression mechanism by inserting and turning the Compression Driver (177305) counterclockwise. To remove the nail, insert the Adaptor before removing all the screws. If broken screws are needed to be removed, use the vise grips.
31. It is essential that the proper Operative Technique be followed for Orthofix Ankle Hindfoot Nailing (AHN) System implantation. Refer to the Orthofix Ankle Hindfoot Nailing (AHN) System Operative Technique.

POSSIBLE ADVERSE EFFECTS

1. Loosening, bending or breakage of implanted components.
2. Loss of anatomic positioning with malunion.
3. Scar formation possibly causing pain and/or neurological compromises around nerves.
4. Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke, excessive bleeding, phlebitis, wound or bone necrosis, wound infection or damage to blood vessels or nerves.
5. Non-union or delayed union, which may lead to implant breakage.
6. Metal sensitivity, or allergic reaction to a foreign body.
7. Pain, discomfort, or abnormal sensations due to the presence of the device.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the devices by the surgeon. Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE the Orthofix Ankle Hindfoot Nailing (AHN) System.

MATERIALS

The Orthofix Ankle Hindfoot Nailing (AHN) System is comprised of implantable titanium parts while the instrumentation is made of stainless steel, aluminum alloy and plastic. Those components of the instrumentation that come into contact with patients are made of surgical grade stainless steel or plastic.

MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

The Orthofix Ankle Hindfoot Nailing System has not been evaluated for safety in the MR environment. It has not been tested for heating or unwanted movement in the MR environment. The safety of Orthofix Ankle Hindfoot Nailing System in the MR environment is unknown. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

STERILITY

All implantable components of the Orthofix Ankle Hindfoot Nailing (AHN) System are supplied STERILE. The instrumentation is supplied NON-STERILE (except two guide wires which are gamma sterilized and provided in sterile packaging) and requires cleaning and sterilization prior to use, following the respective recommended procedures. Please review the product label to determine the sterility of each device.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION




Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.


Disassembling of the Targeting Assembly

Targeting Assembly


1. Unscrew the Nail Attachment Rod (177110) from Targeting Assembly simply by turning it anti-clockwise.
2. Unscrew the Heel Cup (177125) by turning it by hand anti-clockwise.

3

3. Remove Locking Cams (177026) and untighten Locking Screw (A).

4

4. Depress button (351350) on Distal Targeting Arm (177120), align the WHITE dot of the Distal Targeting Arm with the BLACK dot on the shaft of the Proximal Targeting Arm (177100) and remove the Proximal Targeting Arm. In newest Distal Targeting Arm the button can be removed.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations


Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning


1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.

For the Targeting Assembly take special care in

- cleaning the threaded hole at the end of the Nail Attachment Rod (177110):

- cleaning the internal threads of the compression nut of the Heel Cup (177125)


- cleaning threaded and non-threaded holes in the Distal Targeting Arm (177120)

- cleaning threaded and non-threaded holes in the Distal Targeting Arm (177120), in latest versions of Distal Targeting Arm, the button can be removed to allow a more accurate cleaning

- cleaning the button and reassemble it taking into account the correct positioning, the CE0123 marking shall face the proximal part of the handle.

- cleaning Proximal Targeting Arm (177100)

- cleaning the internal threads of the Nut on the Targeting Assembly (177100)

- cleaning the threads and threaded holes of all instruments

6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
 - d. Final rinsing with deionized water for 3 minutes;
 - e. Thermal disinfection at least 90°C or 194° F (max 95 °C or 203° F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.

The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfector when it completes the cycle.
13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132° C (270° F)	132° C (270° F)	134° C (273° F)	134° C (273° F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%

- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

IMPLANTABLE DEVICE*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol ⓘ reported on the product label.

After the removal from the patient, the implantable device* must be disposed of properly.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device which is intended:

Any device intended to be totally / partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is also considered an implantable device.

NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol ⓘ reported on the label or are indicated in the "Instructions For Use" supplied with the products.

The re-use of "SINGLE USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

INFORMATION

For further information please contact your local representative or call 800.266.3349 (within US) or +39 0456719000 (Customer Service International, outside US).

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description
 	Consult instructions for Use
	Single use. Do not reuse
	STERILE. Sterilised by irradiation
	NON STERILE
 	Catalogue number
	Expiry date (year-month - day)
 	CE marking in conformity to applicable European Directives/Regulations
 	Date of manufacture
	Do not use if package is opened or damaged
Rx Only	CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

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Orthofix Srl