Rival Bite™

System

EN	Bite Compression Screws™ 2	JA	Rival Bite コンプレッションボーンスクリュー 86
IT	Viti a compressione Rival Bite 9	ZH	Rival Bite 加压螺钉 96
FR	Vis de compression Rival Bite 16	CS	Kompresní šrouby Rival Bite 104
DE	Rival Bite Kompressionsschrauben 23	PL	Śruby kompresyjne Rival Bite 111
ES	Tornillos de compresión Rival Bite 30	SL	Kompresijski vijaki Rival Bite 118
DA	Rival Bite-kompressionsskruer 37	KO	Rival Bite 압박 나사 125
FI	Rival Bite -kompressoruuvit 44	AR	Rival Bite Compression Screws 133
NO	Rival Bite kompresjonsskruer 51	BG	Rival Bite компресиращи винтове 141
NL	Rival Bite Compressieschroeven 58	LT	„Rival Bite“ suspaudimo varžtai 148
PT	Parafusos de compressão Rival Bite 65	MS	Skru Mampatan Rival Bite 155
SV	Rival Bite Kompressionsskruvar 72	RO	Șuruburi cu compresie Rival Bite 162
EL	Βίδες συμπίεσης Rival Bite 79	TR	Rival Bite Kompresyon Vidaları 169

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

0123

PQBITE E 01/22 (0424536)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online.

Important information - please read prior to use
See also instruction leaflet PQRMD for reusable medical devices

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Telephone +39 045 671 9000 - Fax +39 045 671 9380

GENERAL DESCRIPTION

The RIVAL BITE Headed Cannulated and Headless Compression Screw System is comprised of sterile, cannulated, partially- threaded- titanium alloy screws. The RIVAL BITE Screws are self-drilling and self-tapping and are offered in two head designs: headed cannulated and headless compression versions. Both screw types are available in a variety of diameters and lengths. In addition to the range of above mentioned screws the product line also includes washers and guide wires, as part of the system. The corresponding instrumentation necessary for insertion is found in the RIVAL Instrument Case and/or in sterile packs. The unique low profile head design, cutting characteristics and vast screw options of the RIVAL BITE Screws provide extensive versatility in one comprehensive system.

INTENDED USE and INDICATIONS OF USE

The RIVAL BITE Compression Screws are indicated for use in bone reconstruction, osteotomy, arthrodesis, fracture repair, and fracture fixation of bones appropriate for the size of the device. BITE Compression Screws are intended for use over a guide pin or wire. RIVAL washers may be used with the headed screws in cases where the patient has poor bone quality.

Examples of indications of use for which RIVAL BITE Compression screws may be used are:

- Minimally invasive fracture/joint reconstructions
- Multiple-fragment joint fractures
- Simple epiphyseal fractures
 - Fractures of the head of the humerus
 - Fractures of the head of the tibia
 - Cooper fractures of the tibia
 - Fractures of the radius
- Fractures of the wrist, ankle, elbow and shoulder
- Scaphoid fractures and other fractures of the hand
- Metatarsal fractures and other fracture of the foot
- Ligament fixation of the proximal humerus
- Ligament avulsion injuries (Apophysis)
- Fractures of small joint bones
 - Malleolar fractures
 - Navicular fractures
- Fractures of the calcaneus and talus
- Arthrodesis of the ankle joint
- Avulsion fracture
- Fractures of the tarsal region
- Osteotomies

RIVAL BITE Compression Screws are not intended for use in the spine.

CONTRAINDICATIONS

Use of the RIVAL BITE Compression Screws is contraindicated:

- In the presence of an active infection, or suspected latent infection or marked local inflammation in or about the affected area, or in patients who are immunocompromised.
- In patients with conditions such as limited blood supply, insufficient quantity or quality of bone and /or previous infections that may tend to delay healing and increase the possibility of infection and/or rejection of the implant.
- In patients with foreign body sensitivity. Where material sensitivity is suspected, a test should be performed to rule out this possibility prior to implantation.
- Patients with previous sensitivity to titanium.
- Obesity. An overweight or obese patient can produce loads on the implant that can lead to failure of the fixation of the device or to failure of the device itself.
- In conditions that limit the patient's ability and/or willingness to cooperate and follow instructions during the healing process.
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

WARNINGS

- Preoperative procedures, knowledge of surgical technique and proper selection and placement of the implant are important aspects for the successful use of the device by the surgeon.
- Implanted devices should never be re-used or re-sterilized. An explanted metal implant should never be re-implanted.
- Any implant that comes into contact with bone, blood and/or bodily fluids and was not implanted shall be discarded.

- Instruments, guide wires and screws are to be treated as sharps.
- RIVAL BITE Compression Screws, washers and guide wires are intended for single use only.
- Use of an undersized screw in areas of high functional stresses may lead to implant fracture and failure.
- Use of an oversized screw may lead to bone damage and/or breakage.
- Advise the patient to report any unusual changes to the operative site to the surgeon. Closely monitor patients if a change in the fixation site has been detected. The surgeon should evaluate the possibility of subsequent clinical failure and discuss with the patient the need for any measures deemed necessary to aid healing.
- The design of the instruments must not be modified in any way. RIVAL BITE Compression Screws are not intended for screw attachment or fixation in the spine.
- RIVAL cases/trays have not been validated for use with other Orthofix products and/or other medical device manufacturers' products.
- Do not use RIVAL BITE Compression Screws in combination with other Orthofix products and/or other medical device manufacturers' products, unless otherwise specified. The use of products from other manufacturers with RIVAL implants can cause additional risks to the patient due to the possible contamination of the material and misalignment of the implant to the instruments.

PRECAUTIONS

- Federal Law (U.S.A) restricts this device to sale by or on the order of a licensed physician.
- Deformed screw recesses can increase the risk of breakage and impair the accurate fit of the screw.
- Careful handling and storage of the product is required. Scratching or damage to the component can significantly reduce the strength and fatigue resistance of the products.
- Exactly align the screwdriver/screw head connection coaxially to avoid mechanical damage to the implant or the screwdriver.
- When engaging the fixation screw, axial pressure of the screwdriver into the screw head must be adequately applied to ensure that the tip is fully inserted into the screw head.
- In the final phase of screw insertion, the underside of the screw head contacts the bone and a steep rise in resistance is clearly perceptible. Sensitive tightening of the screw must be exercised to reduce the risk of mechanical damage to the screw, screwdriver, or bony hole.
- Users must ensure the cleanliness of the instruments, specifically cannulated instruments to reduce the risk of binding during use. Further, prior to surgical use, user must ensure all devices function as intended.
- Ensure that drilling and cutting tools are sharp.
- It is particularly important to screen with an image intensifier during guide wire insertion and whenever cannulated instruments are advanced over a guide wire. Frequent screening should also be carried out during screw insertion. In all cases, the benefit of fluoroscopy should be weighed against the risk from radiation exposure on an individual patient basis.

POTENTIAL ADVERSE EFFECTS

- Loosening, bending, cracking or fracturing of the device or premature loss of fixation with the bone
- Delayed union, malunion or nonunion of the fracture site
- Increased fibrous tissue response around the fracture site
- Early or late infection, both deep and/or superficial
- Nerve damage
- Metal sensitivity reactions
- Shortening of the affected bone/fracture site
- Deep venous thrombosis
- Avascular necrosis
- Pain and soft tissue irritation at the implantation site
- Hardware migration/dislocation
- Wound breakdown







MRI SAFETY INFORMATION

RIVAL BITE Compression Screws have not been evaluated for safety and compatibility in the MR environment. The implants have not been tested for heating, migration or image artifact in the MR environment. The safety of the RIVAL BITE Compression Screws implants in the MR environment is unknown. Scanning a patient who has these devices may result in patient injury.

STERILE & NON-STERILE PRODUCTS

- Orthofix provides certain devices STERILE while others are provided NON-STERILE. Please review the product label to determine the sterility of each device.
- Devices or kits provided STERILE are labeled as such. Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.
- In case of sterile screw provided in double tube, do not use if the shrink wrap around the tube is not present and intact.
- Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Please refer to section below.
- Product integrity and performance are assured only if packaging is undamaged.

In case of screw provided in double tube, follow the handling instructions below.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7–10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
 - d. Final rinsing with deionized water for 3 minutes;
 - e. Thermal disinfection at least 90°C or 194° F (max 95 °C or 203° F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfector when it completes the cycle.
13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132° C (270° F)	132° C (270° F)	134° C (273° F)	134° C (273° F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

IMPLANTABLE DEVICE*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol "⊗" reported on the product label.

After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

Manufacturer contact

Please contact your local Orthofix sales representative for further details.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description
 	Consult instructions for use or consult electronic instructions for use Caution: Consult instructions for use for important cautionary information
	Single use. Do not reuse
	STERILE. Sterilised by irradiation
	NON STERILE
 	Catalogue number Batch code
	Expiry date (year-month-day)
 	CE marking in conformity to applicable European Directives/Regulations
 	Date of manufacture Manufacturer
	Do not use if package is opened or damaged
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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Orthofix Srl