EN	Guided Growth Plate System Plus 2	JA	ガイドグロースプレートシステムプラス 38
IT	Sistema di placche Guided Growth Plus 5	ZH	导向生长钢板升级版 42
FR	Système de plaque Guided Growth Plus 8	CS	Systém dlah pro řízený růst Plus 46
DE	Das Guided Growth System Plus 11	PL	System wspomagający wzrost kości Plus 49
ES	Sistema de crecimiento guiado con placas Plus 14	SL	Sistem za voden rast kosti Guided Growth Plate System Plus 52
DA	Guided Growth Plate System Plus 17	KO	Guided Growth Plate System Plus 55
FI	Kasvua ohjaava levyjärjestelmä Plus 20	AR	نظام صفائح تقويم النمو بلس 59
NO	Guided Growth Plate System Plus 23	BG	Система за направляван растеж Плюс 62
NL	Plaatsysteem voor gecontroleerde groei Plus 26	LT	Valdomo augimo plokštelių sistema „Plus“ 65
PT	Sistema de Placa de Crescimento Guiado Plus 29	MS	Sistem Plat Pertumbuhan Berpandu Tambahan 68
SV	Stödjande tillväxtplattesystem Plus 32	RO	Sistemul de plăci de creștere ghidată Plus 71
EL	Σύστημα Καθοδηγούμενης Πλάκας Αύξησης Plus 35	TR	Kılavuzlu Büyüme Plakası Sistemi Plus 74

ORTHOFIX®
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

PQ EPP C 01/18

Important information - please read prior to use**Guided Growth Plate System Plus**

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Telephone +39 045 671 9000 - Fax +39 045 671 9380

GENERAL DESCRIPTION

The Guided Growth Plate System Plus consists of different sizes of eight-Plate and quad-Plate plus different cannulated and solid screw options. The device is designed for the gradual correction of pediatric congenital as well as acquired deformities in both the upper and lower extremities, provided that the physis (growth plates) are not fused.

The plates feature a contoured waist and low profile for pediatric usage. There is a center hole in the plates for a temporary guide pin to be implanted to aid application and removal of the plate. The plates are attached to the external surface of the bone over the growth plate by screws. These screws are not locked to the plate, but rather are allowed to swivel and diverge in their position as bone growth occurs. The implant acts like a flexible hinge, permitting growth at the growth plate to gradually straighten the limb.

MATERIAL

The implants are made from Titanium alloy (Ti-6Al-4V).

INDICATIONS FOR USE

The Guided Growth Plate System Plus is intended for redirecting the growth of long bones. It is indicated for gradually correcting angular growth deformities in growing children.

Specific conditions/diseases for which the device will be indicated include: valgus, varus or flexion, extension deformities of the knee (femur and/or tibia); valgus, varus, or plantar flexion deformities of the ankle; valgus or varus deformities of the elbow (humerus), and radial or ulnar deviation, flexion or extension deformities of the wrist (radius).

CONTRAINDICATIONS

The Guided Growth Plate System Plus should not be used for adult deformities: skeletal maturity precludes further growth.

The Guided Growth System Plus should not be used where the growth plate has closed due to trauma or infection since further growth is precluded.

WARNINGS

1. The plates already have a 10 degrees bending and may be bent an additional 10 degrees for more severe anatomical bone contours. Bending beyond 20 degrees offset could result in injury or reoperation due to breakage.
2. Particular care should be taken that bone screws do not enter the joints or damage the growth plates in growing children.
3. Bone plates and screws must not be reused. If any implants has come into contact with any body fluid it should be considered to have been used. In case repositioning of implants is required, new implants should be used.
4. The Guided Growth Plate System Plus has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Guided Growth Plate System Plus in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
5. This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

PRECAUTIONS

1. It is essential that proper operative technique be followed for implantation.
2. Examine all components carefully prior to use. Product integrity, sterility (in the case of sterile products) and performance are assured only if the packaging is undamaged.
3. Do not use if packaging is compromised or if a component is believed to be faulty, damaged or suspect.
4. Do not combine Guided Growth Plate System Plus implantable components with those from other systems, including the previous version of the Guided Growth Plate System.
5. Before inserting the screws, ensure the plate is flush to the bone, bending up to 10 additional degrees if necessary.
6. Advance bone screws in an alternating manner until the screws are fully seated in the plate and the plate is flush with the bone.
7. Do not use the tap with power drill, only manually.
8. Make sure not to over tap and stop once the tap reaches the mechanical stop.
9. Use the correct drill-bit when pre-drilling.
10. After insertion, ensure correct positioning of the implants under image intensification.
11. Re-operation may be necessary to reposition or replace the bone plate and screws.
12. The implants should be removed at conclusion of treatment, but anyway before the screws reach their maximum angle.
13. Choose the appropriate implant size and post-operative care for each patient and application.

POSSIBLE ADVERSE EFFECTS

- Pain and stiffness at the surgery site.
- Damage to the growth plate.
- Rebound or failure to achieve the desired angular correction and possible over-correction.
- Breaking of the plate or loosening, bending, or breaking of the bone screws.
- Complications associated with metal sensitivity.
- Intrinsic risks associated with anesthesia and surgery.

MRI SAFETY INFORMATION

The Guided Growth Plate System Plus has not been evaluated for safety and compatibility in the MR environment.

It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Guided Growth Plate System Plus in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the plate or bone screws. Proper implantation and patient selection will greatly affect the results. Patients must be closely followed since achievement of the correct angle may require plate removal as early as 3 to 12 months after insertion. Rebound overgrowth after plate removal can sometimes require repeat surgery.

STERILITY

Orthofix provides Guided Growth Plates and Screws in sterile or non-sterile configuration. Labelling and product codes permit an easy identification of configurations.

The products supplied NON-STERILE require cleaning and sterilization prior to use, following the respective recommended procedures.

Cleaning, Sterilization and Maintenance

Cleaning is an essential pre-requisite to ensure effective disinfection or sterilization. All instruments must be cleaned prior to use and after each re-use of the instruments. Cleaning may be done following the validated cleaning process described below and using washing instruments and machines validated by the hospital.

WARNINGS

- Aluminium based instruments are damaged by alkaline (pH>7) detergents and solutions.
- Anodised coating is damaged by detergents with free halogen ions or sodium hydroxide.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED.

Limitations and restrictions on reprocessing

- Repeated reprocessing has minimal effect on reusable instruments;
- End of life is normally determined by wear and damage due to use.

Point of use

- It is recommended that instruments are reprocessed as soon as is reasonably practicable following use, putting them inside a sink filled with cold water (< 40°C) for at least 10 minutes and removing gross soiling with a soft cloth or a soft brush.
- Don't use a fixating detergent or hot water as this can cause the fixation of residue which may influence the result of the reprocessing process.

Containment and transportation

- Follow hospital protocols when handling contaminated and bio-hazardous materials;
- All used surgical instruments should be regarded as contaminated. Their handling, collection and transportation should be strictly controlled to minimize any possible risks to patients, personnel and any area of the healthcare facility.

Preparation for decontamination

- Disassemble the devices where applicable. See Orthofix operative techniques for further details.

Cleaning: manual

1. Soak the single components in the cleaning solution. Orthofix recommends the use of a pH-neutral enzymatic cleaning agent. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
2. Brush the single components in the cleaning solution with a surface soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens
3. Ultrasonic clean the single components in a degassed cleaning solution. Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
4. Brush the single components in running tap water
5. Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer
6. Rinse the components in sterile or deionized water
7. Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Cleaning: automated

When the devices to be cleaned have lumens or present complexity, a preliminary manual cleaning may be required.

- 1) Place all the instruments into washer baskets.
- 2) Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
- 3) Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature. Orthofix recommends that cycle steps are at least as follows:
 1. prewash at low temperatures
 2. main wash
 3. rinse with demineralized water
 4. thermal disinfection
 5. drying
- 4) Rinse the components in sterile or deionized water
- 5) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Maintenance, inspection and testing

All instruments and product components should be visually inspected for cleanliness.

If the equipment does not appear visually clean, cleaning steps must be repeated until the products do not result visually clean.

All instruments and product components should be visually inspected for any signs of deterioration that may cause failure in use (such as cracks or damage to surfaces) and functions tested before being sterilized (see detailed Operative Technique Manuals and Instructions for use). Particular attention should be given to:

- Cannulated devices (NB: cannulated drill-bits are single-patient use only).
- Cutting edges: Discard blunt or damaged instruments.
- If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- When instruments form part of an assembly, check assembly with matching components.

Packaging

- Wrap the tray before sterilization with a sterilization wrap approved in compliance with the local requirements (in USA: FDA-cleared) or insert it into a rigid sterilization container in order to prevent contamination after sterilization.
- Do not include additional systems or instruments in the sterilization tray. Sterility cannot be guaranteed if the sterilization tray is overloaded.

Sterilization

- Steam sterilization is recommended. Gas plasma and EtO sterilization should be avoided as they were not validated for Orthofix products;
- Use a validated, properly maintained and calibrated steam sterilizer;
- Do not exceed 140°C (284°F);
- Do not stack trays during sterilization
- Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

STEAM STERILIZER TYPE	GRAVITY	PRE-VACUUM
Preconditioning pulses	N/A	4 pulses
Exposure Temperature	132°C (270°F)	132°C (270°F)
Minimum Holding Time	15 minutes	4 minutes
Drying Time	30 minutes	30 minutes

Storage

Store the sterilized instrument in a dry, clean, and dust-free environment at room temperature.

SERVICE

For any information please contact your local distributor.

CAUTION

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

IMPLANTABLE DEVICE*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol "⊗" reported on the product label.

After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

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EN All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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Orthofix Srl