# Fitbone™

## Intramedullary Lengthening System

**ORTHOFIX Srl**  
**Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy**  
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PQFBC E 02/25 (0433563)

**US EN**

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online

**Important information - please read prior to use**

**See also instruction leaflet PQFBP for FITBONE Control Set and Retraction Control Set**

**See also instruction leaflet PQFBR or PQRMD for FITBONE LOCKING SCREWS, STANDARD and REVISION LOCKING SCREWS and FITBONE reusable instruments**

**This Instruction For Use (IFU) is for US market**

**FITBONE™ TAA**

**Orthofix Srl**  
**Via delle Nazioni 9 - 37012 Bussolengo (VR) Italy**  
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

**MEDICAL DEVICE INFORMATION****DESCRIPTION**

The FITBONE TAA System consists of the intramedullary lengthening nail with a Receiver connected by a bipolar feed line. The external FITBONE Control Set consists of control electronics and a Transmitter. The energy required for the distraction process is transmitted to the Receiver implanted beneath the skin by applying the Transmitter. There is no contact between the implanted intramedullary lengthening nail and the surface of the body.

The FITBONE TAA intramedullary lengthening nail is available in various configurations. You can find the technical data of the different intramedullary lengthening nails FITBONE in the respective data sheets and dimension sheets. These are provided by the manufacturer upon request.

**INTENDED PURPOSE AND INDICATIONS**

FITBONE TAA is an intramedullary lengthening system for limb lengthening of the femur and tibia.  
FITBONE TAA intramedullary lengthening system is indicated for adult and pediatric (greater than 12 through 21 years of age) patients.

**CONTRAINDICATIONS**

1. Patients with any open wounds or areas with poor soft tissue coverage near the operative site.
2. Patients with anatomic deformities which prevent the device from fitting
3. Patients with poor bone quality that would prevent adequate fixation of the device
4. Patients with compromised capacity for healing
5. Patients with metal allergies and sensitivities
6. Patients in which the implant would cross open, healthy epiphyseal growth plates
7. Blood supply limitations, peripheral vascular disease or evidence of inadequate vascularity
8. Insufficient intramedullary space which would lead to cortical weakening or vascular damage during an implantation
9. Patients with a body weight of > 100 kg for TAA11/13
10. Patients with a body weight of > 50 kg for TAA09
11. Differences in leg length of less than 20 mm
12. No free access for proximal insertion of the intramedullary lengthening nail (e.g. coxa valga)
13. No reliable exclusion of bone infection
14. Expected non-compliance, mentally ill patient, or patient with clouded consciousness
15. Pregnancy

**INTENDED PATIENTS**

FITBONE TAA intramedullary lengthening system is indicated for adult and pediatric (greater than 12 through 21 years of age) patients.

**INTENDED USERS**

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal). One component of the system (FITBONE Control Set) is also intended to be used by the patient or caregiver. There are no specific requirements for the patient or caregiver. The HCP shall instruct the patient or caregiver on its correct use.

**NOTES FOR USE****Retraction Function**

In case of unintended over-distraction, a FITBONE Retraction Control Set can be used to activate the nail retraction function. The Retraction Control Set is only available upon request from Orthofix or your distributor. The Retraction Control Set must not be handed over to the patient. Detailed instructions on safe use of the retraction function are available in the FITBONE Retraction Control Set Quick Guide.

## **Implant Removal**

The FITBONE TAA is not a permanent implant and must be removed. Removal of the implant is recommended when, according to the treating surgeon, the regenerated bone can support a sufficient load. In general, a period of 1 to 1 ½ years after implantation is recommended. If the explantation is delayed or not carried out, the FITBONE intramedullary lengthening nail may break. Detailed instructions on implant removal are available in the dedicated FITBONE Operative Techniques: FITBONE TAA Retrograde Femur Removal, FITBONETAA Antegrade Femur Removal and FITBONE TAA Tibia Removal.

## **DISCLAIMER**

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

## **MATERIAL**

The implants are made from implant grade stainless steel, conforming to ASTM F138 and ISO-5832, and have also some parts, in contact with the body, made of other materials (NuSil MED-4870, NuSil MED-1511, NuSil MED-4750, NuSil MED2-4502). All these materials meet the biocompatibility requirements for implants.

## **WARNINGS**

1. Participation of the surgeon in a FITBONE training course is one of the basic prerequisites for autonomous use of the FITBONE TAA System. The indication for treatment, intraoperative procedure and postoperative aftercare are the responsibility of the surgeon.
2. The supplied components constitute a treatment system and must not be substituted under any circumstances without the written consent of the manufacturer.
3. Never put the grid plate under the patient directly without padding. The material can crack and lead to injury if the grid plate does not rest flat over the entire surface.
4. Align the beam path perpendicular to the surface of the plate. The double line must be located precisely in the center of the image for the evaluation. Do not change the position and rotation of the extremity during the measurement.
5. The FITBONE TAA intramedullary lengthening nail is delivered in plasma sterilized condition. Indicators on the sterile packaging serve as proof of plasma sterilization. Check the sterile packaging (both) for integrity and expiration date before opening. Do not use if the sterile packaging is damaged.
6. If the leg is not deformed, make sure that the osteotomy does not result in axial and torsional deviations. The insertion of Schanz screws for torsional control is therefore strongly recommended. These must be placed in such a way that they do not interfere with the insertion of the intramedullary lengthening nail into the medullary canal.
7. Direct use of striking tools to impact a cone or tube can damage the surgical instruments. Use the corresponding Cone-Sinker CS and Tube-Sinker TS (CS=Cone- Sinker, TS=Tube-Sinker) to avoid damage.
8. Axial correction is determined definitely by the choice of the entry point for the FITBONE TAA intramedullary lengthening nail in the bone and the reaming direction up to the osteotomy site.
9. It is strongly recommended that you monitor the entire reaming process with the image converter in two planes to detect any reaming errors in good time.
10. Never use reamers with a flexible shaft as this can lead to unnecessary weakening of the wall or cause the FITBONE TAA intramedullary lengthening nail to jam later.
11. Make sure that the front surfaces remain in corresponding contact at osteotomy level, as otherwise inadequate bone regeneration can be expected.
12. Make sure that the ventral circumference is not weakened at the proximal end of the ream (risk of fracture). In the diaphyseal region, corrective possibilities are limited, and the risk of perforation is high if the reaming procedure is not monitored. Repeated monitoring by the image converter in two planes during the reaming process is essential. Improper reaming can lead to heat damage to the bone and serious complications.
13. It must be possible to insert the dummy without resistance. It may be necessary to ream the medullary cavity further under monitoring by the image converter. The FITBONE Tapping Tool (60000317) is intended solely for later explantation.
14. After unpacking the intramedullary lengthening nail, check it for integrity, especially the bipolar feed line. Never use a hammer to drive the intramedullary lengthening nail FITBONE TAA into/from the medullary cavity as doing so could damage the implant.
15. Imprecisely set locking holes (oval, funnel-shaped, ragged) do not allow the implant to be fastened sufficiently. If the implant is insufficiently fastened, it may be helpful to use revision screws.
16. After unpacking the Receivers, check it for integrity, especially the bipolar feed line.
17. The connector must be inserted into the coupling completely. A marking is provided for this purpose. Do not hold the coupling or connector with surgical instruments and avoid bending the coupling, connector or supply line as this can lead to damage or dislocation.
18. The ideal distance between Receiver and Transmitter is approx. 5 mm. Power transfer is optimal at this distance. Avoid distances of more than 10 mm as such distances can negatively impact the function of the treatment system.
19. In continuous Doctor mode operation, the Transmitter can reach a maximum temperature of 47.2 °C.
20. Increased caution should be exercised for continuous energy transmission. Excessive distraction can cause neuronal damage.
21. Please advise your patients that only the attending physician is permitted to use the area under the housing cover.
22. Ultrasound and magnetic field therapies are not permitted during the active phase of distraction and while the FITBONE TAA intramedullary lengthening nail is in the bone.
23. The distraction rate for adolescents can deviate due to the accelerated rate of growth, and can lead to an earlier formation of new bone.
24. The FITBONE TAA is not a permanent implant and must be removed. Removal of the implant is recommended when the regenerated bone can support a sufficient load and the stability of the bone has been restored. In general, a period of 1 to 1 ½ years after implantation is recommended. If the explantation is delayed or not carried out, the FITBONE intramedullary lengthening nail may break.
25. Retraction must not be carried out over the initial total length of the intramedullary lengthening nail since this might cause jamming of the implant. During energy transmission, monitor the retraction with continuous X-ray monitoring. FITBONE is intended for single use (single-use implant). It must not be retracted or extended multiple times, or used for re-implantation. The FITBONE implant is not suitable for applying compression forces. If the Retraction Transmitter is positioned in an incorrect angle, an unintended distraction cannot be excluded. Stability of the bone is no longer guaranteed if the nail was retracted after a distraction and then distracted again (see table "Tbl-1").
26. If the yellow "Transmit" LED remains off, retraction might not work. In this case, please read the chapter "Malfunctions" in the Instructions for Use of the Control Set.
27. During the active phase and the early healing phase, the carrying capacity of the intramedullary lengthening nail FITBONE TAA is limited to a partial load of 20 kg (contact with the sole of the foot). Any exceeding load may cause the intramedullary lengthening nail to break.
28. The patient must receive and use walking aids.
  - The patient should be prescribed appropriate, regular physical therapy by the doctor
  - The patient must avoid unforeseen/unwanted excess weight-bearing, e.g. by falling or stumbling, as well as full weight-bearing. This may cause the intramedullary lengthening nail to break.
29. If excessive weight-bearing does occur, the patient must contact the attending physician.
  - The distraction phase should not be interrupted for longer than 2 days at most because the risk of premature bone bridging is otherwise present.
30. Non-ionising radiation is used for energy and data transmission. Electromagnetic and magnetic pulses can cause malfunctions.
  - Note that radio equipment with transmission frequencies below 500 kHz may inadvertently lengthen the intramedullary lengthening nail. Keep away from potential sources of such electromagnetic fields as, for example:
  - Industrial equipment with wireless energy transmission, including production facilities and logistics centres. Please observe any posted warnings relating to increased electromagnetic radiation.

Radio masts / radio towers used as time-signal transmitters.

Comply with the special safety precautions with regard to electromagnetic compatibility (EMC) according to the FITBONE Control Set Instructions for Use.

Only use the components (e.g. power supply units or cables) provided by the manufacturer.

31. The FITBONE TAA intramedullary lengthening nail and control electronics must not be placed directly beside or stacked with other devices.

Monitor/check the correct operation of the FITBONE TAA System if the FITBONE Control Set is nevertheless positioned in this way.

## PRECAUTIONS

The treating physician will decide about the risk and benefit of the following treatments.

Any form of electrical therapy in which current is passed through the patient's body should be avoided on the affected limb, as should therapeutic ultrasound.

- If this cannot be avoided, the function of the FITBONE TAA intramedullary lengthening nail must be carefully monitored to immediately identify any disruptions.
- If medical treatments are required in which an electric current is passed through the patient's body from an external source, turn off the control electronics and advise the patient to carefully monitor the functioning of the device during distraction over the next 4 to 5 days.

Except for diagnostic X-rays, the manufacturer has no experience with how the FITBONE TAA intramedullary lengthening nail responds to high-energy ionizing radiation. Such treatments must be avoided for the duration of the distraction phase, if possible.

- If such treatments are necessary, it is your decision to perform the MRI as the treating physician.

The leakage currents expected in the body of the patient when using defibrillators can negatively affect the FITBONE TAA intramedullary lengthening nail. Therefore, the use of defibrillators should be avoided if possible.

- If the use of a defibrillator cannot be avoided, carefully monitor the function of the implant during distraction in the following 4 to 5 days.

Please make sure to use the specific instrumentation according to the locking screw model as described in the Operative Technique.

## Only for "STANDARD LOCKING SCREW" and "REVISION LOCKING SCREW" codes:

1. To ensure proximal locking screw length, the tip of the trocar must go flat against the bone.
2. Choosing the appropriate length for the proximal locking screw is critical for the effective and safe stabilization of bone fragments. The length of the screw should be determined by using the dedicated instruments as described in the operative technique.
3. After proximal locking screw length measurements, make sure to hold the targeting handle in the same position: any movement of the targeting handle may result in a mismatch between the screw and the hole.
4. Overtightening of the screw should be avoided and the head of the screw should come just in contact with the cortex. Stop insertion when resistance is felt.

## POSSIBLE ADVERSE EVENTS

In addition to the general risks involved in surgical intervention, the following side effects may occur in some cases despite correct treatment:

- Slight tingling to severe pain in the affected limb, particularly during and after distraction
- Temporary limited mobility of the affected limb

## MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

Fitbone™ TAA nails and other components are MR unsafe. Keep away from MRI examination room.

Instrumentation have not been evaluated for safety in the MR environment. They have not been tested for heating or unwanted movement in the MR environment. The safety of instrumentations in the MR environment is unknown. Performing an MR exam on a person who has these medical devices may result in injury or device malfunction.

## SPECIFIC INFORMATION ON THE DEVICE

Detailed information on the selection of suitable implants, accessories and related devices, including locking screws, is available in the dedicated FITBONE Operative Techniques: FITBONE TAA Retrograde Femur Application, FITBONE TAA Retrograde Femur Removal, FITBONE TAA Antegrade Femur Application, FITBONE TAA Antegrade Femur Removal, FITBONE TAA Tibia Application and FITBONE TAA Tibia Removal.

## RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

### Implantable Device\*

The "SINGLE USE" implantable device\* of Orthofix is identified through symbol "⊗" reported on the product label. After the removal from the patient, the implantable device\* has to be discarded.

The re-use of an implantable device\* introduces contamination risks for users and patients.

The re-use of an implantable device\* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

### Non Implantable Device

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of a "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

## STERILE & NON-STERILE PRODUCTS

The medical devices intended for implantation are H2O2-plasma-sterilized and are supplied in a packaging system comprising a sterile barrier system and protective packaging in a sales packaging. The sterile barrier system and the protective packaging are only soft packs with a chemical indicator. Chemical indicators on both sterile bags provide proof of plasma sterilization. To open, tear the sterile bags at the respective sealed seam. The outer packaging (protective packaging) is used to forward the device to the operating room area. Check the protective packaging for damage or perforation before surgery. If it is damaged, it can no longer be considered aseptic. The implants must not be used.

Note the expiration date on the package.

For information on sterilization of Instrument Set and screws, see Chapter "Cleaning, disinfection, sterilization".

## INSTRUCTIONS FOR PROCESSING AND REPROCESSING

**For the STANDARD and REVISION LOCKING SCREWS and FITBONE reusable instruments listed in the table below please use the processing and reprocessing instructions described in the document PQRMD.**

Part #	Description	Part #	Description	Part #	Description
60001028	FITBONE CONE C13	60001854	FITBONE DUMMY TAA0960-F-220	60001024	FITBONE TUBE T12/11-XL
60001029	FITBONE CONE C13+	60001855	FITBONE DUMMY TAA0960-T-220	60001025	FITBONE TUBE T12/10-XL
60001030	FITBONE CONE C13++	60001925	FITBONE DUMMY TAA0940-F-200	60001026	FITBONE TUBE T12/09-XL
60001031	FITBONE CONE C12	60001927	FITBONE DUMMY TAA0940-T-200	60001027	FITBONE TUBE T12/08-XL
60001032	FITBONE CONE C11	60000576	FITBONE SCREWDRIVER SW3.5	60001044	FITBONE TUBE T14/13-S
60001060	FITBONE CONE C15	60000406	FITBONE SCREWDRIVER SW3.5 CANNULATED	60001045	FITBONE TUBE T13/12-S
60001061	FITBONE CONE C15+	60000218	FITBONE SPACE HOLDER TAA	60001046	FITBONE TUBE T12/11-S
60001062	FITBONE CONE C15++	60001439	FITBONE TARGETING DEVICE 45°/90°	60001047	FITBONE TUBE T12/10-S
60001063	FITBONE CONE C14	60000403	FITBONE TROCAR D4.5MM	60001048	FITBONE TUBE T12/09-S
60001888	FITBONE CONE C13 RETROPAELLAR	60000405	FITBONE TROCAR T-HANDLE D4.5MM	60001049	FITBONE TUBE T12/08-S
60000408	FITBONE DEPTH GAUGE	60001014	FITBONE TUBE T14/13-M	60001050	FITBONE TUBE T16/15-M
60000400	FITBONE DRILL SLEEVE D4.5MM	60001015	FITBONE TUBE T13-12-M	60001051	FITBONE TUBE T15/14-M
60000402	FITBONE DRILL SLEEVE D8MM	60001016	FITBONE TUBE T12/11-M	60001052	FITBONE TUBE T16/15-L
60000822	FITBONE DUMMY TAA1180-F-245	60001017	FITBONE TUBE T12/10-M	60001053	FITBONE TUBE T15/14-L
60000832	FITBONE DUMMY TAA1160-T-225	60001018	FITBONE TUBE T12/09-M	60001054	FITBONE TUBE T14/13-L
60001139	FITBONE DUMMY TAA1160-F-225	60001019	FITBONE TUBE T12/08-M	60001055	FITBONE TUBE T13/12-L
60001244	FITBONE DUMMY TAA1140-T-205	60001020	FITBONE TUBE T12/11-L	60001056	FITBONE TUBE T16/15-XL
60001248	FITBONE DUMMY TAA1140-F-205	60001021	FITBONE TUBE T12/10-L	60001057	FITBONE TUBE T15/14-XL
60001495	FITBONE DUMMY TAA1180-T-245	60001022	FITBONE TUBE T12/09-L	60001058	FITBONE TUBE T14/13-XL
60001623	FITBONE DUMMY TAA1380-F-245	60001023	FITBONE TUBE T12/08-L	60001059	FITBONE TUBE T13/12-XL

**For all other FITBONE reusable instruments and FITBONE LOCKING SCREWS follow the processing and reprocessing instructions described below and in the document PQFBR.**

These reprocessing instructions have been written in accordance with ISO 17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions provided.

## ADDITIONAL INFORMATION

### INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical followup and of the eventual removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician. When handing over the FITBONE Control Set to the patient, the HCP shall hand over also the paper copy of the "FITBONE Control Set Instructions for Use" and verbally instruct the patient on how to use the FITBONE Control Set correctly for the lengthening procedures at home. Before leaving the healthcare facility after implantation of FITBONE TAA intramedullary nail, the patient must be supplied with the Intramedullary Lengthening Nail implant card and the Subcutaneous Energy Receiver implant card. Both implant cards are supplied together with the implants. Before handing each implant card over to the patient, one of the labels available on the implant sterile packaging must be attached in the appropriate area on the implant card and the patient information must be filled-in by the HCP, following the instructions provided in the leaflet PQIPC.

### NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

### MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

Symbol	Description	
<b>MD</b>	Medical Device	
 	Consult instructions for use or consult electronic Instructions for Use	CAUTION: Consult instructions for use for important cautionary information
	Single Use. Do not re-use	Orthofix note: discard appropriately after use (treatment) on patient
	Do not resterilize	
<b>STERILE R</b>	Sterile. Sterilized using irradiation	
<b>STERILE VH202</b>	Sterilized using vaporized hydrogen peroxide	
	Non-sterile	
	Double sterile barrier system	
<b>UDI</b>	Unique Device Identifier	
<b>REF</b>	Catalogue number	Batch code
	Use-by date (year-month-day)	
<b>CE</b>  0123	CE marking in conformity to applicable European Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	The item is known to pose hazards in all MRI environments.	
<b>SN</b>	Serial number	
<b>Rx Only</b>	CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician	
	Electric voltage	
	Follow instruction for use	
	Non-ionizing electromagnetic radiation	
	<p>Disposal. In accordance with the WEEE 2012/19/UE Directive, the product and all its parts (cable, batteries, accessories, etc.) may not be treated as domestic waste.</p> <p>Please return the FITBONE Control Set to Orthofix after completion of the treatment.</p>	
	Patient name and surname	

Symbol	Description
	Surgery date
	Doctor name and surname
	Patient information website
	Temperature limit
	Humidity limitation
<b>For retraction only</b>	For retraction only

Symbol	Description
	Please connect Intramedullary Lengthening Nail and Receiver via the bipolar plug connection. Please also note the instructions for connecting the connector and coupling in the Operative Technique.
	For the Intramedullary Lengthening Nail TAA09, use 4.0 drill for the locking hole on the telescope side (freehand locking). Please also note the instructions for drilling locking holes in the Operative Technique.

**CA EN**

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online

**Important information - please read prior to use**

**See also instruction leaflet PQFBP for FITBONE Control Set and Retraction Control Set**

**See also instruction leaflet PQFBR or PQRMD for FITBONE LOCKING SCREWS and FITBONE reusable instruments.**

**This Instruction For Use (IFU) is for Canadian market**

**FITBONE™ TAA**

**Orthofix Srl**  
**Via delle Nazioni 9 - 37012 Bussolengo (VR) Italy**  
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

**MEDICAL DEVICE INFORMATION****DESCRIPTION**

The FITBONE TAA nailing system consists of an implantable intramedullary lengthening nail with a Receiver connected by a bipolar feed line. The external FITBONE Control Set consists of a control electronics station and a Transmitter. The power required for the distraction process is controlled by a hermetically enclosed motor which draws the telescope apart. The electro-magnetic field sent from the Transmitter to the Receiver is converted in the Receiver into DC-Voltage to supply the motor of the FITBONE Lengthening Nail with voltage, when actioned. There is no contact between the implanted intramedullary lengthening nail or the implanted receiver and the surface of the body.

The FITBONE TAA nails include 3 different variants in terms of diameter (TAA11, TAA09, TAA13), available in different lengths and lengthening capabilities. TAA11 and TAA09 are available for femur and for tibia. The FITBONE TAA nailing system must be fixed to the bone using FITBONE LOCKING SCREWS. FITBONE LOCKING SCREWS are available in multiple diameters, Ø4,5mm (with long and short thread) and Ø4mm, and multiple lengths. The FITBONE TAA nailing system can be implanted and removed using dedicated instruments. Further information is available in section "Specific Information on the Device", below.

**INTENDED PURPOSE AND INDICATIONS**

INTENDED PURPOSE of the FITBONE TAA System.

FITBONE TAA is an intramedullary lengthening system for limb lengthening of the femur and tibia.

The intended purpose of each individual device is provided in the table below:

Device	Device specific Intended Purpose
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Femur	The intramedullary nail is intended for lengthening of the femur.
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Tibia	The intramedullary nail is intended for lengthening of the tibia.
FITBONE SUBCUTANEOUS ENERGY RECEIVER	The receiver is intended for energy conversion and transfer to the intramedullary nail.
FITBONE CONTROL SET	The control set is an external active device intended to transmit energy to the receiver for the lengthening of the intramedullary nail.
FITBONE RETRACTION CONTROL SET	The retraction control set is an external active device intended to transmit energy to the receiver for the retraction of the intramedullary nail.

**INDICATIONS FOR USE**

FITBONE TAA is an intramedullary lengthening system for limb lengthening of the femur and tibia.

FITBONE TAA intramedullary lengthening system is indicated for adult and pediatric (greater than 12 through 21 years of age) patients.

**CONTRAINDICATIONS**

DO NOT USE the FITBONE TAA if a surgical candidate exhibits or is predisposed to any of the following contraindications

1. Patients with any open wounds / areas with poor soft-tissue coverage near the operative site
2. Patients with anatomic deformities which prevent the device from fitting
3. Patients with poor bone quality that would prevent adequate fixation of the device
4. Patients with compromised capacity for healing
5. Patients with allergies and sensitivities to metals and/or silicones
6. Patients in which the implant would cross open, healthy epiphyseal growth plates
7. Blood supply limitations, peripheral vascular disease or evidence of inadequate vascularity
8. Insufficient intramedullary space which would lead to cortical weakening or vascular damage during an implantation
9. Patients who are unable, either for physical or mental reasons, to limit weight bearing on the operated limb to a maximum of 20 kg until the end of consolidation phase
10. Differences in leg length of less than 20 mm
11. No free access for proximal insertion of the intramedullary lengthening nail (e.g. coxa valga)

12. No reliable exclusion of bone infection
13. Expected non-compliance, mentally ill patient or patient with clouded consciousness
14. Pregnancy
15. Lack of hip head roofing (hip dysplasia) with femoral extensions
16. Lack of stability in the surrounding joints

## **INTENDED PATIENTS**

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations.

FITBONE TAA intramedullary lengthening system is indicated for adult and pediatric (greater than 12 through 21 years of age) patients.

## **INTENDED USERS**

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal). One component of the system (FITBONE Control Set) is also intended to be used by the patient or caregiver. There are no specific requirements for the patient or caregiver. The HCP shall instruct the patient or caregiver on its correct use.

## **NOTES FOR USE**

### **Retraction Function**

In case of unintended over-distraction, a FITBONE Retraction Control Set can be used to activate the nail retraction function. The Retraction Control Set is only available upon request from Orthofix or your distributor. The Retraction Control Set must not be handed over to the patient. Detailed instructions on safe use of the retraction function are available in the FITBONE Retraction Control Set Quick Guide.

### **Implant Removal**

The FITBONE TAA is not a permanent implant and must be removed. Removal of the implant is recommended when, according to the treating surgeon, the regenerated bone can support a sufficient load. In general, a period of 1 to 1 ½ years after implantation is recommended. If the explantation is delayed or not carried out, the FITBONE intramedullary lengthening nail may break. Detailed instructions on implant removal are available in the dedicated FITBONE Operative Techniques: FITBONE TAA Retrograde Femur Removal, FITBONE TAA Antegrade Femur Removal and FITBONE TAA Tibia Removal.

## **DISCLAIMER**

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including intraoperative procedure post-operative care).

## **MATERIAL**

The implants are made from implant grade stainless steel, conforming to ASTM F138 and ISO-5832, and have also some parts, in contact with the body, made of other materials (NuSil MED-4870, NuSil MED-1511). All these materials meet the biocompatibility requirements for implants.

## **WARNINGS**

1. Never use a hammer to drive or remove the Intramedullary Lengthening Nail FITBONE TAA into/ from the medullary cavity since this could damage the implant.
2. If an eccentric cone has been used, do not rotate it while removing.
3. There is a risk of excessive cortical weakening and even perforation, leading to a fracture during treatment, particularly with the front cutting reamers.
4. Never use front cutting rigid reamers in the area of the final position of the tip of the FITBONE intramedullary nail.
5. Failure to adhere to instructions regarding limited weight bearing (maximum of 20 kg) on the operated limb until the end of consolidation may result in nail or bone breakage and/or other related serious complications leading to the need for immediate additional medical treatment, reoperation, revision or extended care.

## **PRECAUTIONS**

1. Check the sterile packaging for integrity and expiration date before opening. Do not use if the sterile packaging is damaged.
2. Half-pins might be used to avoid unintentional axial and torsional deviations. These must be placed in such a way that they do not interfere with the insertion of the intramedullary lengthening nail into the medullary canal.
3. Do not directly hammer the cone because that will damage the cone end and tubes will not slide on the cone as required.
4. Monitor the entire reaming process with the image intensifier in two planes to detect any reaming errors in time.
5. Never use reamers with a flexible shaft as this can lead to unnecessary weakening of the wall or cause the FITBONE TAA intramedullary lengthening nail to jam later.
6. Do not weaken the cortex as this could increase the risk of fracture during treatment.
7. After unpacking the lengthening intramedullary nail and the receiver, check them for integrity, especially the bipolar feedline.
8. Imprecisely set locking holes (oval, funnel-shaped, ragged) do not allow the implant to be fastened sufficiently. If the implant is insufficiently fastened, it may be helpful to use revision screws.
9. Make sure the white ring is in connection with the start of the coupling.
10. Do not hold the coupling or cable with surgical instruments and avoid bending the coupling or cable as this can lead to damage or unwanted disconnection.
11. Power transfer is optimal at approximately 5 mm distance from the receiver. Avoid distances of more than 10 mm under the skin as such distances can negatively affect the function of the treatment system.
12. In doctor mode continuous operation, the Transmitter can reach a maximum temperature of 47.2 °C.
13. Excessive distraction can cause neuronal damage.
14. Please advise your patients not to remove the switch cover from the control set and not touch the switches.
15. The treating physician will decide about the risk and benefit of the following treatments.
  - Electrical therapy
    - Any form of electrical therapy in which current is passed through the patient's body should be avoided on the affected limb as it could negatively affect the FITBONE TAA intramedullary lengthening nail.
    - The leakage currents expected in the body of the patient when using defibrillators can negatively affect the FITBONE TAA intramedullary lengthening nail.

- Ultrasound therapy
  - Any form of therapeutic ultrasound should be avoided on the affected limb. Effects are not tested but could potentially result in unexpected harm and malfunction of the FITBONE TAA intramedullary lengthening nail.
- Electromagnetic therapy
  - Any form of magnetic and electromagnetic field application should be avoided on the affected limb as it could negatively affect the FITBONE TAA intramedullary lengthening nail and cause unexpected harm.
- High-energy ionizing radiation therapy
  - The manufacturer has no experience with how the FITBONE TAA intramedullary lengthening nail responds to high-energy ionizing radiation. Such treatments must be avoided for the duration of the distraction phase, if possible.

If the above therapies cannot be avoided:

- The function of the FITBONETAA intramedullary lengthening nail must be carefully monitored to immediately identify any disruptions and also during distraction over the following 4 to 5 days.
- Never perform lengthening while the above therapies are applied.

X-ray and ultrasound imaging is allowed on the affected limb. The manufacturer does not know of any events on reactions of the implantable FITBONE parts due to x-rays computer tomography. However, potential artefacts might be too high for proper imaging quality.

16. The distraction rate for adolescents should be appropriate to the accelerated rate of bone growth which could lead to a faster formation of new bone.
17. Retraction must not be carried out over the initial total length of the intramedullary lengthening nail since this might cause jamming of the implant. During energy transmission, monitor the retraction with the stethoscope and continuous X-ray monitoring.
18. FITBONE must not be retracted or extended multiple times, or used for re-implantation.
19. The FITBONE implant is not suitable for applying compression forces.
20. If the Retraction Transmitter is positioned in an incorrect angle, an unintended distraction cannot be excluded.
21. During the distraction and consolidation phase, weightbearing on the operated leg should be partial and limited to 20 kg (contact with the sole of the foot). Any exceeding load may cause the FITBONE TAA intramedullary nail to break.
22. The patient must avoid unforeseen/unwanted excess weight-bearing as well as too early full weight-bearing. This may cause the FITBONE intramedullary lengthening nail to break. If excessive weight-bearing does occur, the patient must contact the treating physician.
23. The distraction phase should not be interrupted for more than 2 days at most because the risk of premature bone bridging is otherwise present.
24. Non-ionising radiation is used for energy and data transmission. Electromagnetic and magnetic pulses can cause malfunctions.
  - Note that radio equipment with transmission frequencies below 500 kHz may inadvertently lengthen the FITBONE intramedullary lengthening nail. Keep away from potential sources of such electromagnetic fields as, for example:
    - a) Industrial equipment with wireless energy transmission, including production facilities and logistics centres. Please observe any posted warnings relating to increased electromagnetic radiation.
    - b) Radio masts / radio towers used as time-signal transmitters.
  - Comply with the special safety precautions with regard to electromagnetic compatibility (EMC) according to the FITBONE Control Set Instructions for Use (PQFBP).
  - Only use the components (e.g. power supply units or cables) provided by the manufacturer.
25. FITBONE TAA devices must not be stored with other electromagnetic, magnetic, ionizing, wireless or HF devices. Monitor/check the correct operation of the FITBONE TAA System if the FITBONE Control Set is nevertheless stored in this way.
26. The drainage tube should slide without any resistance through the bony canal and the bipolar feedline should not enter and affect the adjacent joint. If not, check if the implant is inserted deeply enough (distance to notch 1cm).
27. Perform additional corrections or place blocking screws only with the dummy nail inside the bone, never while the FITBONE intramedullary nail is inserted.
28. Intraoperative functional test: prior to implanting the Fitbone nail, check the functionality of the nail, by activating it via FITBONE Control Set. The operating noise of the nail (using the supplied stethoscope) and the flashing control light will confirm the functionality.
29. The cable must enter the bony tunnel as close as possible to the point it exits from the implant. If not, it may occur that a fibre of the patella tendon is captured, which might cause severe damage to the cable.
30. If the cable is not secured, patella tendon movement may result in damage to the cable.

FITBONE devices should be used together with their corresponding FITBONE implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the FITBONE when used in conjunction with other Orthofix devices if not specifically indicated in the Operative Technique.

## **POSSIBLE ADVERSE EVENTS**

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP. Possible adverse events resulting from the usage of the FITBONE TAA intramedullary nail are:

- Superficial infection
- Deep infection
- Non-union, delayed union or malunion
- Damage to bone or surrounding tissues resulting from surgery or treatment
- Damage to local vascularization resulting in bone or surrounding tissue necrosis
- Joint contracture, dislocation, instability or loss of range of motion
- Residual deformities, persistence or recurrence of the initial condition subject to treatment
- Bone fracture during or after treatment
- Pain, discomfort or abnormal sensations due to the presence of the device
- Compartment syndrome
- Complex regional pain syndrome
- Wound healing complications
- Injury to the healthy limbs, buttocks or head due to positioning during surgery

- Injury due to postural changes during post-operative treatment
- Loosening, bending, breakage or migration of implants
- Foreign body reactions due to implantable devices
- Reoperation to replace a component or entire nail configuration
- Events caused by intrinsic risks associated with anesthesia and surgery

#### **MRI (Magnetic Resonance Imaging) SAFETY INFORMATION**

FITBONE is MR unsafe.

- Keep away from MRI examination rooms.

#### **EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE**

##### **EXPECTED CLINICAL BENEFIT**

###### **PRECISION**

- Protection of soft-tissue with minimally-invasive instrumentation
- Instruments designed not only for lengthening but for optimal limb alignment
- Safe alignment assessment and blocking screw placement with dummy (trial) nail

###### **BENEFIT TO SURGEON**

- Stability, less chance of bone on growth
- Minimally invasive, adds safety and precision when used with reamers correcting deformity

###### **BENEFIT TO PATIENT**

- Stability
- Possible reduced risk of premature consolidation
- Bone conserving
- Minimally invasive, less scarring, less chance of debris in joint/wound

##### **PERFORMANCE CHARACTERISTICS OF THE DEVICE**

###### **FIT**

- No over-reaming necessary, therefore, it is a bone-conserving procedure
- Bone is reamed to exactly fit the implant (step reamer)

###### **POWER**

- Reliable power direct to the nail
- High distraction force independent of nail size
- No soft-tissue limitations
- Power transferred to the nail with a small portable Control Set
- No magnets within the Fitbone Control Set

###### **BENEFIT TO SURGEON**

- Powerful 1200N distraction force, independent of nail size
- Reliable power direct to the nail
- The receiver is located under the skin, which allows very good energy transmission regardless of patient size
- No over-reaming necessary, facilitates exact reaming and position of implant for limb alignment
- Bone is reamed to exactly fit the implant
- Lightweight, easy to handle and transport
- No magnets within the Fitbone Control Set
- Accidental retraction not possible
- It is possible to do fast distraction or retraction
- Safe alignment assessment and blocking (poller) screw placement with dummy (trial) nail before inserting definitive implant
- 9mm nail option

###### **BENEFIT TO PATIENT**

- Easy to use/power the nail
- Lightweight, easy to handle and transport, compact – fits in briefcase
- User friendly
- No magnets within the Fitbone Control Set
- Accidental retraction not possible
- Guided procedure
- 9mm nail available and all TAA nails now indicated for adults AND adolescent patient (>12 years)

**LINK TO THE SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

This is the SSCP location after the launch of European Database on Medical Devices/Eudamed: <https://ec.europa.eu/tools/eudamed>

<b>Device</b>	<b>Basic UDI-DI</b>
FITBONE TAA System	805050436ImpActive328AJY
FITBONE SUBCUTANEOUS ENERGY RECEIVER	805050436EnRec3ACT4K
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Tibia	805050436LenTib3ACTJP
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Femur	805050436LenFem3ACTDE
FITBONE CONTROL SET	805050436CSLen3ACTNK
FITBONE RETRACTION CONTROL SET	805050436CSRet3ACTTT

**Specific Information on the Device**

Detailed information on the selection of suitable implants, accessories and related devices, including FITBONE LOCKING SCREWS, is available in the dedicated FITBONE Operative Techniques: FITBONE TAA Retrograde Femur Application, FITBONE TAA Retrograde Femur Removal, FITBONE TAA Antegrade Femur Application, FITBONE TAA Antegrade Femur Removal, FITBONE TAA Tibia Application and FITBONE TAA Tibia Removal.

**RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE****Implantable Device\***

The "SINGLE USE" implantable device of Orthofix is identified through symbol “⊗” reported on the product label or are indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

**Non Implantable Device**

The "SINGLE USE" non implantable device of Orthofix is identified through symbol reported on the label or are indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

**STERILE & NON-STERILE PRODUCTS**

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products indicators on the sterile packaging serve as proof of sterilisation. Product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

**CONTAINMENT, TRANSPORTATION AND DISPOSAL**

Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection, transportation and disposal of used devices must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility. Please note that the FITBONE Control Set must not be disposed but must be returned to Orthofix after completion of the treatment.

## INSTRUCTIONS FOR PROCESSING AND REPROCESSING

**For the FITBONE reusable instruments listed in the table below please use the processing and reprocessing instructions described in the document PQRMD.**

Part #	Description	Part #	Description	Part #	Description
60001028	FITBONE CONE C13	60001854	FITBONE DUMMY TAA0960-F-220	60001024	FITBONE TUBE T12/11-XL
60001029	FITBONE CONE C13+	60001855	FITBONE DUMMY TAA0960-T-220	60001025	FITBONE TUBE T12/10-XL
60001030	FITBONE CONE C13++	60001925	FITBONE DUMMY TAA0940-F-200	60001026	FITBONE TUBE T12/09-XL
60001031	FITBONE CONE C12	60001927	FITBONE DUMMY TAA0940-T-200	60001027	FITBONE TUBE T12/08-XL
60001032	FITBONE CONE C11	60000576	FITBONE SCREWDRIVER SW3.5	60001044	FITBONE TUBE T14/13-S
60001060	FITBONE CONE C15	60000406	FITBONE SCREWDRIVER SW3.5 CANNULATED	60001045	FITBONE TUBE T13/12-S
60001061	FITBONE CONE C15+	60000218	FITBONE SPACE HOLDER TAA	60001046	FITBONE TUBE T12/11-S
60001062	FITBONE CONE C15++	60001439	FITBONE TARGETING DEVICE 45°/90°	60001047	FITBONE TUBE T12/10-S
60001063	FITBONE CONE C14	60000403	FITBONE TROCART D4.5MM	60001048	FITBONE TUBE T12/09-S
60001888	FITBONE CONE C13 RETROPAELLAR	60000405	FITBONE TROCART HANDLE D4.5MM	60001049	FITBONE TUBE T12/08-S
60000408	FITBONE DEPTH GAUGE	60001014	FITBONE TUBE T14/13-M	60001050	FITBONE TUBE T16/15-M
60000400	FITBONE DRILL SLEEVE D4.5MM	60001015	FITBONE TUBE T13-12-M	60001051	FITBONE TUBE T15/14-M
60000402	FITBONE DRILL SLEEVE D8MM	60001016	FITBONE TUBE T12/11-M	60001052	FITBONE TUBE T16/15-L
60000822	FITBONE DUMMY TAA1180-F-245	60001017	FITBONE TUBE T12/10-M	60001053	FITBONE TUBE T15/14-L
60000832	FITBONE DUMMY TAA1160-T-225	60001018	FITBONE TUBE T12/09-M	60001054	FITBONE TUBE T14/13-L
60001139	FITBONE DUMMY TAA1160-F-225	60001019	FITBONE TUBE T12/08-M	60001055	FITBONE TUBE T13/12-L
60001244	FITBONE DUMMY TAA1140-T-205	60001020	FITBONE TUBE T12/11-L	60001056	FITBONE TUBE T16/15-XL
60001248	FITBONE DUMMY TAA1140-F-205	60001021	FITBONE TUBE T12/10-L	60001057	FITBONE TUBE T15/14-XL
60001495	FITBONE DUMMY TAA1180-T-245	60001022	FITBONE TUBE T12/09-L	60001058	FITBONE TUBE T14/13-XL
60001623	FITBONE DUMMY TAA1380-F-245	60001023	FITBONE TUBE T12/08-L	60001059	FITBONE TUBE T13/12-XL

**For all other FITBONE reusable instruments and FITBONE LOCKING SCREWS follow the processing and reprocessing instructions described below and in the document PQFBR.**

These reprocessing instructions have been written in accordance with ISO 17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions provided.

### Warnings

- Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors, including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use is the best way to determine whether a device has reached its end of serviceable life.
- Due to their constant use, the instruments are subject to natural wear and tear and usage-specific damage, which may influence the end of the product service life. In a non-exhaustive list, signs of damage and wear include: Corrosion (i.e., rust, pitting), discoloration, deep scratches, flaking, abrasions and cracks. Improperly functioning or defective and excessively worn instruments, as well as instruments with unrecognizable markings, missing or re-moved (worn away) part numbers, may not be used and must be disposed of and replaced.
- Corrosion is a form of destruction or wear which can be caused as a result of chemical reactions, including:
  - Damaged surface structures
  - The effects of human body fluids in the event of prolonged contact with the instruments
  - The excessive effects of certain solutions: Saline, iodine solutions, chloride or stronger acids, alkaline solutions and incorrectly used disinfectants
  - Insufficient water quality when cleaning, disinfecting, sterilizing using steam or rinsing instruments, e.g. through the use of corroded water pipes, the penetration of rust, metal or dirt particles in steam sterilizers etc.
  - If rust forms, this can be transferred onto other instruments. Ensure that contact is avoided, as this is very dangerous during sterilization.
  - Not observing the guidelines about concentrations and temperatures stipulated by the disinfectant manufacturer: If these concentrations and temperatures are significantly exceeded, for some materials this may lead to discoloration and/or corrosion.
  - This can also be the case if the instruments are not rinsed adequately after cleaning or disinfection.
- Equipment, users, cleaning agents and processes all contribute to the effectiveness of the treatment. The clinical institution must guarantee that the selected treatment stages are safe and effective.
- During cleaning, you must carefully check which cleaning agent may be used and with which method. Please observe the diluting and application requirements.
- Suitable protective equipment must be worn when dealing with contaminated or potentially contaminated materials, instruments and products.
- In order to guarantee the adequate transport of sterile instruments, the sterile container, the lid and the tray insert must be cleaned, disinfected and sterilized while empty.
- Locking screws (implants) should not be cleaned and disinfected together with the instruments.
- The maximum weight of a loaded tray insert during sterilization must be < 10 kg.

### POINT OF USE

The instruments should be pre-cleaned and dried within 30 minutes after use, to minimize the risk of residues drying before cleaning and thus prevent irreparable damage.

### CONTAINMENT AND TRANSPORTATION

It is recommended to cover contaminated instruments during transportation in order to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

## **PREPARATION FOR CLEANING**

Instruments with removable parts must be taken apart. The responsible person must have been appropriately trained for disassembly or assembly. For specific instructions on affected instruments refer to PQFBR.

### **Pre-cleaning**

To pre-clean the instruments thoroughly, they must be treated in an ultrasonic bath and then cleaned manually.

- Ultrasonic bath with neodisher MediZym 0.5% (v/v) and soft drinking-quality water
- Before placing an instrument in the ultrasonic bath, any bodily fluids must be removed from the instruments using absorbent, lint-free paper towels;
- Avoid bubbles on the surface and sound shadows when placing them;
- When filling the ultrasonic bath, make sure that the instruments cannot hit each other;
- Ultrasonic treatment of the instruments (10 minutes, 35 kHz, water temperature <40°C)
- Rinse the instruments under running water
- Make sure to remove all visible contaminants with
  - soft nylon brushes
  - soft round nylon brushes
- Inner diameter or holes must be cleaned with material-specific round brushes and rinsed well at least three times for at least one minute and rinse at least five times with at least one disposable 10ml syringe.

 Once the pre-cleaning step is completed, the instruments must always undergo automated cleaning and disinfection.

 The integrity of the instruments must be checked before and after cleaning and disinfection.

## **CLEANING**

### **General considerations**

It is recommended that cleaning is carried out with a correctly installed, qualified and regularly maintained washing and disinfection machine equipped with injector units, baskets and inserts. The validation in accordance with EN ISO 15883-1 / ISO 15883-1 was carried out using the Miele PG 8535 washing and disinfection machine.

- In the event of deviations from the process detailed here, equivalent cleaning agents and a procedure validated by the end user must be used. Otherwise, this may shorten the lifespan of the instruments and locking screws, and this would be the responsibility of the hospital.
- The operating instructions supplied by the manufacturer of the cleaning machines and cleaning agents must be observed.

### **Automated cleaning and disinfection**

In accordance with the operating instructions, the instruments are inserted into the washing and disinfection machine using a cleaning trolley (e.g. Miele E 450/1 69545003D). Instruments which can be taken apart must be disassembled prior to cleaning. Depth gauge for sleeves (60000408). Cannulae must be connected to the rinse connections (i.e. injector nozzle) of the washing and disinfection machine.

Load trays laying the instruments on the dedicated brackets, in accordance with the markings visible on the lower side of the tray lid.

- Avoid contact between the instruments during the washing process;
- Remove washed items immediately after the end of the wash program to avoid corrosion created by residual moisture.

The process recommended by Orthofix Srl is intended for the following cleaning and disinfection methods:

Cleaning and disinfection device: PG 8535, Miele

Cleaning program:

- a) Pre-rinsing with soft, drinking-quality water at < 10 °C for 2 minutes
- b) Cleaning with Neodisher MediClean forte 0.7% (v/v); (\*) at 55 °C for 10 minutes
- c) Rinsing I with demineralized water at < 10 °C for 1 minute
- d) Neutralizing (\*) with Neodisher Z 0.1% (v/v) at < 10 °C for 2 minutes
- e) Rinsing II with demineralized water at < 10 °C for 1 minute
- f) Thermal Disinfection at 93 °C for 5 minutes or until achieving an A° value of > 3000.
- g) Drying at 100 °C for 25 minute

(\*) Only necessary for strong alkaline cleaners, however these must also have a pH value of < 12

After the program is complete, the machine must be unloaded and the medical products checked. If the devices still show contamination residues, they must undergo another cleaning process in the machine.

 Due to structural, material and functional properties, these are medical products of the critical risk class A & B according to Robert Koch Institute (RKI) recommendations, which is why they may only be treated by the machine.

 The integrity of the instruments must be checked after the cleaning program.

## **MAINTENANCE, INSPECTION AND FUNCTION TESTING**

Before every sterilization or use, all parts of the instrument inventory and all the locking screws must be checked for cleanliness, integrity and proper functioning.

- During the visual inspection, under good light with a magnifying glass, if applicable, all instruments must be examined for contamination, corrosion and residues;
- If the visual inspection shows that the instruments have not been cleaned properly, they must be cleaned and disinfected again.
- Damaged parts may not be used and must be re-ordered for replacement from Orthofix Srl.

## PACKAGING

Disassembled instruments must be re-assembled by trained personnel. For the assembly procedure: see the "PREPARATION FOR CLEANING" section above. After the inspection and functional test, the dried instruments are placed in the tray and the locking screws are placed in the screw box. In order to prevent contamination after sterilization, Orthofix recommends using one of the following packaging systems:

- a) Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained to mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non-woven fabrics made of spun bond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA-cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b) Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every sterile barrier packaging not validated by Orthofix must be validated by individual health care facility according to instruction from the manufacturer. When using equipment and processes different from that validated by Orthofix, the health care facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10 kg.

## STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AAMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for reprocess FITBONE medical products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The gravity cycle was validated but it is suggested only when no other options are available. The gravity cycle was not validated for sterilization in rigid containers but only in wraps.

## STORAGE

After sterilization, the wrapped medical products should be stored in a location with restricted access.

- This location should be well ventilated and offer sufficient protection from dust, moisture, insects, vermin and extreme temperatures.

Before the products can be stored, they must be allowed to cool-down to ambient temperature after sterilization.

Before opening, the sterile wrap must be thoroughly examined for any damage and the expiry date must be checked.

## DISCLAIMER

The instructions provided herein have been validated by Orthofix Srl for the initial preparation of a device for first clinical use or for reprocessing of multiple use devices prior to re-use. It remains the responsibility of the individual responsible for reprocessing to ensure that the reprocessing, performed using equipment, materials and personnel in the initial preparation / reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

## ADDITIONAL INFORMATION

### INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical follow-up and of the eventual removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician. When handing over the FITBONE Control Set to the patient, the HCP shall hand over also the paper copy of the "FITBONE Control Set Instructions for Use" and verbally instruct the patient on how to use the FITBONE Control Set correctly for the lengthening procedures at home. Before leaving the healthcare facility after implantation of FITBONE TAA intramedullary nail, the patient must be supplied with the Intramedullary Lengthening Nail implant card and the Subcutaneous Energy Receiver implant card. Both implant cards are supplied together with the implants. Before handing each implant card over to the patient, one of the labels available on the implant sterile packaging must be attached in the appropriate area on the implant card and the patient information must be filled-in by the HCP, following the instructions provided in the leaflet PQIPC.

### NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

## MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability. Please contact the manufacturer if the labels on the packaging are damaged or not legible.

Symbol	Description	
<b>MD</b>	Medical Device	
 	Consult instructions for use or consult electronic Instructions for Use	CAUTION: Consult instructions for use for important cautionary information
	Single Use. Do not re-use	Orthofix note: discard appropriately after use (treatment) on patient
	Do not resterilize	
<b>STERILE</b> <b>VH202</b>	Sterilized using vaporized hydrogen peroxide	
	Non-sterile	
	Double sterile barrier system	
<b>UDI</b>	Unique Device Identifier	
<b>REF</b> <b>LOT</b>	Catalogue number	Batch code
	Use-by date (year-month-day)	
<b>CE</b> <b>CE</b> 0123	CE marking in conformity to applicable European Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	The item is known to pose hazards in all MRI environments.	
<b>SN</b>	Serial number	
<b>Rx Only</b>	CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician	
	Electric voltage	
	Follow instruction for use	
	Non-ionizing electromagnetic radiation	

Symbol	Description
	Disposal. In accordance with the WEEE 2012/19/UE Directive, the product and all its parts (cable, batteries, accessories, etc.) may not be treated as domestic waste. Please return the FITBONE Control Set to Orthofix after completion of the treatment.
	Patient name and surname
	Surgery date
	Doctor name and surname
	Patient information website
	Temperature limit
	Humidity limitation
<b>For retraction only</b>	For retraction only

Symbol	Description
<span style="float: right;"></span>	Please connect Intramedullary Lengthening Nail and Receiver via the bipolar plug connection. Please also note the instructions for connecting the connector and coupling in the Operative Technique.
<span style="float: right;"></span>	For the Intramedullary Lengthening Nail TAA09, use 4.0 drill for the locking hole on the telescope side (freehand locking). Please also note the instructions for drilling locking holes in the Operative Technique.

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**See also instruction leaflet PQFBR or PQRMD for FITBONE LOCKING SCREWS and FITBONE reusable instruments**

**This Instruction For Use (IFU) is for NOT for US and CA market**

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	<b>Device specific Intended Purpose</b>
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Femur	The intramedullary nail is intended for lengthening of the femur.
FITBONE TAA INTRAMEDULLARY LENGTHENING NAIL for Tibia	The intramedullary nail is intended for lengthening of the tibia.
FITBONE SUBCUTANEOUS ENERGY RECEIVER	The receiver is intended for energy conversion and transfer to the intramedullary nail.
FITBONE CONTROL SET	The control set is an external active device intended to transmit energy to the receiver for the lengthening of the intramedullary nail
FITBONE RETRACTION CONTROL SET	The retraction control set is an external active device intended to transmit energy to the receiver for the retraction of the intramedullary nail.

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FITBONE TAA intramedullary lengthening system is indicated for adult and pediatric (greater than 12 through 21 years of age) patients

#### CONTRAINDICATIONS

DO NOT USE the FITBONE TAA if a surgical candidate exhibits or is predisposed to any of the following contraindications:

1. Patients with any open wounds / areas with poor soft-tissue coverage near the operative site
2. Patients with anatomic deformities which prevent the device from fitting
3. Patients with poor bone quality that would prevent adequate fixation of the device
4. Patients with compromised capacity for healing
5. Patients with metal allergies and sensitivities
6. Patients in which the implant would cross open, healthy epiphyseal growth plates
7. Blood supply limitations, peripheral vascular disease or evidence of inadequate vascularity
8. Insufficient intramedullary space which would lead to cortical weakening or vascular damage during an implantation
9. Patients with a body weight of > 100 kg for TAA11/13
10. Patients with a body weight of > 50 kg for TAA09
11. Differences in leg length of less than 20 mm
12. No free access for proximal insertion of the intramedullary lengthening nail (e.g. coxa valga)

13. No reliable exclusion of bone infection
14. Expected non-compliance, mentally ill patient or patient with clouded consciousness
15. Pregnancy
16. Lack of hip head roofing (hip dysplasia) with femoral extensions
17. Lack of stability in the surrounding joints

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## **INTENDED USERS**

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal). One component of the system (the control set) is also intended to be used by the patient or caregiver. There are no specific requirements for the patient or caregiver. The HCP shall instruct the patient or caregiver on its correct use.

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## **WARNINGS**

1. Never use a hammer to drive or remove the Intramedullary Lengthening Nail FITBONE TAA into/ from the medullary cavity as doing so could damage the implant.
2. If an eccentric cone has been used, do not rotate it while removing.
3. There is a risk of excessive cortical weakening and even perforation, leading to a fracture during treatment, particularly with the front cutting reamers.
4. Never use front cutting reamers in the area of the final position of the tip of the FITBONE.

## **PRECAUTIONS**

1. Check the sterile packaging for integrity and expiration date before opening. Do not use if the sterile packaging is damaged.
2. Bone screws might be used to avoid unintentional axial and torsional deviations. These must be placed in such a way that they do not interfere with the insertion of the intramedullary lengthening nail into the medullary canal.
3. Do not directly hammer the cone because that will damage the cone end and tubes will not slide on the cone as required.
4. Monitor the entire reaming process with the image intensifier in two planes to detect any reaming errors in time.
5. Never use reamers with a flexible shaft as this can lead to unnecessary weakening of the wall or cause the FITBONE TAA intramedullary lengthening nail to jam later.
6. Do not weaken the cortex as this could increase the risk of fracture during treatment.
7. After unpacking the lengthening nail and the receiver, check them for integrity, especially the bipolar feedline.
8. Imprecisely set locking holes (oval, funnel-shaped, ragged) do not allow the implant to be fastened sufficiently. If the implant is insufficiently fastened, it may be helpful to use revision screws.
9. Make sure the white ring is in connection with the start of the coupling.
10. Do not hold the coupling or cable with surgical instruments and avoid bending the coupling or cable as this can lead to damage or unwanted disconnection.
11. Power transfer is optimal at approximately 5 mm distance from the receiver. Avoid distances of more than 10 mm under the skin as such distances can negatively affect the function of the treatment system.
12. In doctor mode continuous operation, the Transmitter can reach a maximum temperature of 47.2 °C.
13. Excessive distraction can cause neuronal damage.
14. Please advise your patients not to remove the switch cover from the control set and not touch the switches.
15. The treating physician will decide about the risk and benefit of the following treatments.
  - Electrical therapy
    - Any form of electrical therapy in which current is passed through the patient's body should be avoided on the affected limb as it could negatively affect the FITBONE TAA intramedullary lengthening nail.
    - The leakage currents expected in the body of the patient when using defibrillators can negatively affect the FITBONE TAA intramedullary lengthening nail.
  - Ultrasound therapy
    - Any form of therapeutic ultrasound should be avoided on the affected limb. Effects are not tested but could potentially result in unexpected harm and malfunction of the FITBONE TAA intramedullary lengthening nail.

- Electromagnetic therapy
  - Any form of magnetic and electromagnetic field application should be avoided on the affected limb as it could negatively affect the FITBONE TAA intramedullary lengthening nail and cause unexpected harm.
- High-energy ionizing radiation therapy
  - The manufacturer has no experience with how the FITBONE TAA intramedullary lengthening nail responds to high-energy ionizing radiation. Such treatments must be avoided for the duration of the distraction phase, if possible.

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- The function of the FITBONE TAA intramedullary lengthening nail must be carefully monitored to immediately identify any disruptions and also during distraction over the following 4 to 5 days.
- Never perform lengthening while the above therapies are applied.

X-ray and ultrasound imaging is allowed on the affected limb. The manufacturer does not know of any events or reactions of the implantable FITBONE parts due to x-rays computer tomography. However, potential artefacts might be too high for proper imaging quality.

16. The distraction rate for adolescents should be appropriate to the accelerated rate of bone growth which could lead to a faster formation of new bone.
17. Retraction must not be carried out over the initial total length of the intramedullary lengthening nail since this might cause jamming of the implant. During energy transmission, monitor the retraction with the stethoscope and continuous X-ray monitoring.
18. FITBONE must not be retracted or extended multiple times, or used for re-implantation.
19. The FITBONE implant is not suitable for applying compression forces.
20. If the Retraction Transmitter is positioned in an incorrect angle, an unintended distraction cannot be excluded.
21. During the distraction and consolidation phase, weightbearing on the operated leg should be partial and limited to 20 kg (contact with the sole of the foot). Any exceeding load may cause the FITBONE TAA to break.
22. The patient must avoid unforeseen/unwanted excess weight-bearing as well as too early full weight-bearing. This may cause the intramedullary lengthening nail to break. If excessive weight-bearing does occur, the patient must contact the treating physician.
23. The distraction phase should not be interrupted for more than 2 days at most because the risk of premature bone bridging is otherwise present.
24. Non-ionising radiation is used for energy and data transmission. Electromagnetic and magnetic pulses can cause malfunctions.
  - Note that radio equipment with transmission frequencies below 500 kHz may inadvertently lengthen the intramedullary lengthening nail. Keep away from potential sources of such electromagnetic fields as, for example:
    - a) Industrial equipment with wireless energy transmission, including production facilities and logistics centres. Please observe any posted warnings relating to increased electromagnetic radiation.
    - b) Radio masts / radio towers used as time-signal transmitters.
  - Comply with the special safety precautions with regard to electromagnetic compatibility (EMC) according to the FITBONE Control Set Instructions for Use (PQFBP).
  - Only use the components (e.g. power supply units or cables) provided by the manufacturer.
25. FITBONE TAA devices must not be stored with other electromagnetic, magnetic, ionizing, wireless or HF devices. Monitor/check the correct operation of the FITBONE TAA System if the FITBONE Control Set is nevertheless stored in this way.
26. The drainage tube should slide without any resistance through the bony canal and the bipolar feedline should not enter and affect the adjacent joint. If not, check if the implant is deep enough inserted (distance to notch 1cm).
27. Perform additional corrections or place blocking screws only with the dummy nail inside the bone, never while the Fitbone nail is inserted.
28. Intraoperative functional test: prior to implanting the Fitbone nail, check the functionality of the nail, by activating it via the control set. The operating noise of the nail (using the supplied stethoscope) and the flashing control light will confirm the functionality.
29. The cable must enter the bony tunnel as close as possible to the point it exits from the implant. If not, it may occur that a fibre of the patella tendon is captured, what might cause severe damage to the cable.
30. If the cable is not secured, patella tendon movement may result in damage to the cable.

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## **POSSIBLE ADVERSE EVENTS**

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP. Possible adverse events resulting from the usage of the FITBONE TAA are:

- Superficial infection
- Deep infection
- Non-union, delayed union or malunion
- Damage to bone or surrounding tissues resulting from surgery or treatment
- Damage to local vascularization resulting in bone or surrounding tissue necrosis
- Joint contracture, dislocation, instability or loss of range of motion
- Residual deformities, persistence or recurrence of the initial condition subject to treatment
- Bone fracture during or after treatment
- Pain, discomfort or abnormal sensations due to the presence of the device
- Compartment syndrome
- Complex regional pain syndrome
- Wound healing complications
- Injury to the healthy limbs, buttocks or head due to positioning during surgery
- Injury due to postural changes during post-operative treatment
- Loosening, bending, breakage or migration of implants

- Foreign body reactions due to implantable devices
- Reoperation to replace a component or entire nail configuration
- Events caused by intrinsic risks associated with anesthesia and surgery

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##### PERFORMANCE CHARACTERISTICS OF THE DEVICE

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###### POWER

- Reliable power direct to the nail
- High distraction force independent of nail size
- No soft-tissue limitations
- Power transferred to the nail with a small portable Control Set
- No magnets within the Fitbone Control Set.

###### BENEFIT TO SURGEON

- Powerful 1200N distraction force, independent of nail size
- Reliable power direct to the nail.
- The receiver is located under the skin, which allows very good energy transmission regardless of patient size.
- No over-reaming necessary, facilitates exact reaming and position of implant for limb alignment
- Bone is reamed to exactly fit the implant
- Lightweight, easy to handle and transport.
- No magnets within the Fitbone Control Set.
- Accidental retraction not possible
- It is possible to do fast distraction or retraction
- Safe alignment assessment and blocking (poller) screw placement with dummy (trial) nail before inserting definitive implant
- 9mm nail option

###### BENEFIT TO PATIENT

- Easy to use/power the nail
- Lightweight, easy to handle and transport, compact – fits in briefcase.
- User friendly
- No magnets within the Fitbone Control Set.
- Accidental retraction not possible
- Guided procedure
- 9mm nail available and all TAA nails now indicated for adults AND adolescent patient (>12 years)

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The "SINGLE USE" implantable device of Orthofix is identified through symbol "⊗" reported on the product label or are indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

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**Non Implantable Device**

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or are indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

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## ADDITIONAL INFORMATION

### INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical follow-up and of the eventual removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician. When handing over the FITBONE Control Set to the patient, the HCP shall hand over also the paper copy of the "FITBONE Control Set Instructions for Use" and verbally instruct the patient on how to use the FITBONE Control Set correctly for the lengthening procedures at home. Before leaving the healthcare facility after implantation of FITBONE TAA, the patient must be supplied with the Intramedullary Lengthening Nail implant card and the Subcutaneous Energy Receiver implant card. Both implant cards are supplied together with the implants. Before handing each implant card over to the patient, one of the labels available on the implant sterile packaging must be attached in the appropriate area on the implant card and the patient information must be filled-in by the HCP, following the instructions provided in the leaflet PQIPC.

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Symbol	Description
<span style="font-size: 2em; vertical-align: middle;">(1)</span> <span style="font-size: 2em; vertical-align: middle;">(2)</span> <span style="font-size: 2em; vertical-align: middle;">(3)</span> <span style="font-size: 2em; vertical-align: middle;">(4)</span> <span style="font-size: 2em; vertical-align: middle;">(5)</span>	<p>Please connect Intramedullary Lengthening Nail and Receiver via the bipolar plug connection. Please also note the instructions for connecting the connector and coupling in the Operative Technique.</p>
<span style="font-size: 2em; vertical-align: middle;">Ø 4.0</span> <span style="font-size: 2em; vertical-align: middle;">Ø 4.5</span>	<p>For the Intramedullary Lengthening Nail TAA09, use 4.0 drill for the locking hole on the telescope side (freehand locking). Please also note the instructions for drilling locking holes in the Operative Technique.</p>

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