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PQFBR B 01/23 (0432498)

INSTRUCTIONS FOR CLEANING,  
STERILIZATION AND MAINTENANCE  
OF FITBONE™ NONSTERILE  
LOCKING SCREWS AND REUSABLE  
INSTRUMENTS

## CONTENTS

<b>1. Introduction</b>	<b>3</b>
<b>2. Warnings information</b>	<b>3</b>
<b>3. Material resistance</b>	<b>4</b>
<b>4. Reprocessing Instructions</b>	<b>5</b>
4.1. Pre-treatment at the place of use	5
4.2. Containment and Transportation	5
4.3. Preparation before cleaning	5
4.4. Pre-cleaning	6
4.5. Automated cleaning and disinfection	6
4.6. Drying	7
4.7. Inspection and maintenance	7
4.8. Packaging	7
4.9. Sterilization	8
4.10. Storage	8
<b>5. Responsibilities of a hospital for hired sets</b>	<b>8</b>
<b>6. Disclaimer</b>	<b>9</b>
<b>7. Notice about serious incidents</b>	<b>9</b>
<b>8. List of instruments with material numbers</b>	<b>10</b>
<b>9. List of screws with material numbers</b>	<b>12</b>
<b>10. Interpretation of symbols</b>	<b>13</b>

## 1. Introduction

This document gives general guidance on how FITBONE™ reusable instruments manufactured by Orthofix Srl may be reprocessed to prepare them for use. It also explains how to inspect an instrument to determine if it has reached the end of its serviceable life, must not be used, and must be replaced.

FITBONE™ implants are intended for single use only and are not intended to be reprocessed after use.

The efficacy of the reprocessing also depends on the equipment, operators, cleaning agents and procedures used.

The reprocessing instructions reported in the present document have been validated and they are in accordance with the following standards and TIR AAMI:

- ISO 17664
- ISO 17665
- ISO 15883
- ISO 11607
- AAMI TIR12

Alternative processing methods which are not included in this document may be suitable for treatment, but must be validated by the end user on a case-by-case basis. The instruments may only be handled by individuals with the necessary expertise and training. The end user must be able to assess the potential risks and the corresponding consequences.

## 2. Warnings information

- Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors, including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use is the best way to determine whether a device has reached its end of serviceable life. (frase copiata senza modificarla direttamente da PQRMD)
- Due to their constant use, the instruments are subject to natural wear and tear and usage-specific damage, which may influence the end of the product service life. In a non-exhaustive list, signs of damage and wear include: Corrosion (i.e., rust, pitting), discoloration, deep scratches, flaking, abrasions and cracks. Improperly functioning or defective and excessively worn instruments, as well as instruments with unrecognizable markings, missing or re- moved (worn away) part numbers, may not be used and must be disposed of and replaced.
- Corrosion is a form of destruction or wear which can be caused as a result of chemical reactions, including:
  - Damaged surface structures
  - The effects of human body fluids in the event of prolonged contact with the instruments
  - The excessive effects of certain solutions: Saline, iodine solutions, chloride or stronger acids, alkaline solutions and incorrectly used disinfectants
  - Insufficient water quality when cleaning, disinfecting, sterilizing using steam or rinsing instruments, e.g. through the use of corroded water pipes, the penetration of rust, metal or dirt particles in steam sterilizers etc.
  - If rust forms, this can be transferred onto other instruments. Ensure that contact is avoided, as this is very dangerous during sterilization.
  - Not observing the guidelines about concentrations and temperatures stipulated by the disinfectant manufacturer: If these concentrations and temperatures are significantly exceeded, for some materials this may lead to discoloration and/or corrosion.
  - This can also be the case if the instruments are not rinsed adequately after cleaning or disinfection.
- Equipment, users, cleaning agents and processes all contribute to the effectiveness of the treatment. The clinical institution must guarantee that the selected treatment stages are safe and effective.
- During cleaning, you must carefully check which cleaning agent may be used and with which method. Please observe the diluting and application requirements.
- Suitable protective equipment must be worn when dealing with contaminated or potentially contaminated materials, instruments and products.
- In order to guarantee the adequate transport of sterile instruments, the sterile container, the lid and the tray insert must be cleaned, disinfected and sterilized while empty.
- Locking screws (implants) should not be cleaned and disinfected together with the instruments.
- The maximum weight of a loaded tray insert during sterilization must be < 10kg.

### 3. Material resistance

When selecting the cleaning agent or disinfectant, you must ensure that they do not contain the following ingredients:

- Organic, mineral and oxidizing acids (minimum permissible pH value of 7)
- Strong alkaline solutions (maximum permissible pH value of 11)
- Organic solvents (e.g. alcohol, ether, ketone, gasoline)
- Oxidizing agents (e.g. hydrogen peroxide)
- Aromatic / halogenated hydrocarbons

Recommended cleaning agents are enzymatic, neutral or mildly alkaline cleaning agents

The instruments should not be cleaned using metal brushes or steel wool.

Recommended equipment: Nylon brushes, K-wire, lint-free disposable cloth and disposable sponge

All instruments may only be exposed to temperatures of up to 142 °C (286 °F).

Passivating agents which have been validated for this purpose are exceptions to this. Passivation must take place after cleaning and intermediate rinsing. (neodisher Z)

**Instruments made from anodized aluminum should not come into contact with certain disinfectant solutions or cleaning agents; which is why you should check the solution to be used for its application and compatibility with anodized aluminum before the metal is exposed to it. Aluminum based devices are damaged by alkaline (pH>7) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.**

**The specified parameters for cleaning, disinfection and sterilization are validated according to EN ISO 17665-1 and ANSI AAMI ST79. Exceeding these parameters can shorten the service life of the instruments and is done under the sole responsibility of the end user.**

### 4. Reprocessing Instructions

Instructions for initial preparation of devices, i.e. prior to their first use, starts at Chapter 4.2. Instructions for reprocessing of instruments which have been previously used starts at Chapter 4.1.

#### 4.1. Pre-treatment at the place of use

During and/or immediately after operation, remove coarse contaminants (e.g., bodily fluids, bone meal and tissue) from the instruments.

- Wipe off contaminants using cloths;
- Rinse instruments with a lumen using sterile distilled water to prevent dry blood and bone residues in the inner lumen (saline solution is not suitable);
- Return instruments to their assigned places after use;
- Do not leave the instruments in the rinsing solution.

Prepare for transport for reprocessing.

- The instruments should be pre-cleaned and dried within 30 minutes after use, to minimize the risk of residues drying before cleaning and thus prevent irreparable damage;
- Place the instruments in the correct trays and holders.

**All instruments which were used during the operation are deemed contaminated. The instruments are to be cleaned immediately after use.**

Used instruments must be brought to be transported for cleaning in closed or covered containers to prevent any unnecessary risk of contamination and harm to personnel or the surrounding area.

#### 4.2. Containment and Transportation

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials.

Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.


**Please handle the container with care in order to maintain the sterile barrier system. The container contains sharp and pointed instruments.**

### 4.3. Preparation before cleaning

Instruments with removable parts must be taken apart.

- Affected instruments:

Material number 60000408 - Depth gauge for sleeves

Disassembly procedure:

Depth gauge for sleeves assembled

Unscrew cap with guide housing from button with slider

Remove cap with guide housing from button with slider on the left

Remove sleeve to the right

Depth gauge for sleeves disassembled

The responsible person must have been appropriately trained for disassembly or assembly.

**In general, screws, nuts, studs and other small parts should be stored together and in a well-arranged way.**

### 4.4. Pre-cleaning

**Pre-cleaning applies to all instruments listed in Chapter 7 "List of instruments". Pre-cleaning does not apply to locking screws.**

To pre-clean the instruments thoroughly, they must be treated in an ultrasonic bath and then cleaned manually.

- Ultrasonic bath with neodisher MediZym 0.5% (v/v) and soft drinking-quality water
- Before placing an instrument in the ultrasonic bath, any bodily fluids must be removed from the instruments using absorbent, lint-free paper towels;
- Avoid bubbles on the surface and sound shadows when placing them;
- When filling the ultrasonic bath, make sure that the instruments cannot hit each other;
- Ultrasonic treatment of the instruments (10 minutes, 35 kHz, water temperature <40°C)
- Rinse the instruments under running water
- Make sure to remove all visible contaminants with
  - soft nylon brushes
  - soft round nylon brushes
- Inner diameter or holes must be cleaned with material-specific round brushes and rinsed well **at least three times** for **at least one minute** and rinse **at least five times** with at least one disposable 10ml syringe.

**Once the pre-cleaning step is completed, the instruments must always undergo automated cleaning and disinfection.**

**The integrity of the instruments must be checked before and after cleaning and disinfection (see Chapter 4.4).**

### 4.5. Automated cleaning and disinfection

It is recommended that cleaning is carried out with a correctly installed, qualified and regularly maintained washing and disinfection machine equipped with injector units, baskets and inserts. The validation in accordance with EN ISO 15883-1/ISO 15883-1 was carried out using the Miele PG 8535 washing and disinfection machine.

- In the event of deviations from the process detailed here, equivalent cleaning agents and a procedure validated by the end user must be used. Otherwise, this may shorten the lifespan of the instruments and locking screws, and this would be the responsibility of the hospital.
- The operating instructions supplied by the manufacturer of the cleaning machines and cleaning agents must be observed;

In accordance with the operating instructions, the instruments are inserted into the washing and disinfection machine using a cleaning trolley (e.g. Miele E 450/1 69545003D). Instruments which can be taken apart must be disassembled prior to cleaning i.e. Depth gauge for sleeves (60000408). Cannulae must be connected to the


rinse connections (i.e. injector nozzle) of the washing and disinfection machine.

Instruments with cannulae

Load trays laying the instruments on the dedicated brackets, in accordance with the markings visible on the lower side of the tray lid. :

Instruments without cannulae

60000398	Spiral drill D4.5 L300
60000403	Trocar D4.5
60000405	Trocar D4.5, T-handle
60000384	Screw holder M2 L325
60000576	Screwdriver SW3.5
60000317	Tapping tool
60000688	Test pin D4.5
60000689	Open-end wrench SW14/17
60000408	Depth gauge – button with slider
60001184	Mounting bracket TAA
60001175	Drill guide TAA
60000175	Retaining screw M6 L9
60000003	Locking screw
60000219	Clamping nut TAA
60001439	Spot-film device 45°/90°
60001307	Wire trap 90°
60001415	Step reamer TAA1160
60001179	Step reamer TAA1180
60000832	Dummy TAA1160-T-225
60000822	Dummy TAA1180-F-245
60000513	Retaining screw M4 L400
60000410	Pliers for screws
60001039	K-wire D3 L28
60001139	Dummy TAA1160-F-225
60001185	Drill guide blocking screw TAA11
60001244	Dummy TAA1140-T-205
60001248	Dummy TAA1140-F-205
60001464	GRID plate
60001495	Dummy TAA1180-T-245
60001528	Step reamer TAA1140
60001605	Spiral drill D4.5 L300, conical

- 60001623 Dummy TAA1380-F-245
- 60001854 Dummy TAA0960-F-220
- 60001855 Dummy TAA0960-T-220
- 60001925 Dummy TAA0940-F-200
- 60001927 Dummy TAA0940-T-200
- 60001848 Spiral drill D4.0 L200
- 60001849 Step reamer TAA0960
- 60001938 Step reamer TAA0940

Screws in the screw box

- Avoid contact between the instruments during the washing process;
- Remove washed items immediately after the end of the wash program to avoid corrosion created by residual moisture.

The process recommended by Orthofix Srl is intended for the following cleaning and disinfection methods:

Cleaning and disinfection device: PG 8535, Miele Cleaning program:

	Description	Temperature	Duration
<b>Pre-rinsing</b>	Soft, drinking-quality water	< 10 °C	2 minutes
<b>Cleaning</b>	Neodisher MediClean forte 0.7% (v/v); (*)	55 °C	10 minutes
<b>Rinsing I</b>	Demineralized water	< 10 °C	1 minute
<b>Neutralizing (*)</b>	Neodisher Z 0.1% (v/v)	< 10 °C	2 minutes
<b>Rinsing II</b>	Demineralized water	< 10 °C	1 minute
<b>Thermal Disinfection</b>	-	93 °C	5 minutes or until achieving an AO value of > 3000.
<b>Drying</b>	-	100 °C	for 25 minutes

(\*) Only necessary for strong alkaline cleaners, however these must also have a pH value of < 12

After the program is complete, the machine must be unloaded and the reusable instruments checked (see 4.6). If the devices still show contamination residues, they must undergo another cleaning process in the machine.

**Due to structural, material and functional properties, these are reusable instruments of the critical risk class A & B according to Robert Koch Institute (RKI) recommendations, which is why they may only be treated by the machine.**

**The integrity of the instruments must be checked after the cleaning program (see 4.6).**

#### 4.6. Drying

Every instrument and every locking screw must be completely dry on the inside and outside to prevent the formation of rust or malfunctions. If the instruments are not dry following the drying step in the thermal disinfection process, dry the instruments. Lint-free cloths may be used.

#### 4.7. Inspection and maintenance

Before every sterilization or use, all parts of the instrument inventory and all the locking screws must be checked for cleanliness, integrity and proper functioning.

- During the visual inspection, under good light with a magnifying glass, if applicable, all instruments must be examined for contamination, corrosion and residues;
- If the visual inspection shows that the instruments have not been cleaned properly, they must be cleaned and disinfected again.
- Damaged parts may not be used and must be re-ordered for replacement from Orthofix Srl.

#### 4.8. Packaging

Disassembled instruments must be re-assembled by trained personnel.

- affected instruments:
    - 60000408 - Depth gauge for sleeves
- Disassembly procedure: see 4.2, reverse procedure

After the inspection and functional test, the dried instruments are placed in the tray and the locking screws are placed in the screw box.

In order to prevent contamination after sterilization, Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained to mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA-cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every sterile barrier packaging not validated by Orthofix must be validated by individual health care facility according to instruction from the manufacturer. When using equipment and processes different from that validated by Orthofix, the health care facility should verify that sterility can be achieved using parameters validated by Orthofix.

Do not include additional systems or instruments in the sterilization tray.

Note that sterility cannot be guaranteed if the sterilization tray is overloaded.

The total weight of a wrapped instrument tray should not exceed 10kg.

**Sharp and pointed instruments may damage the sterile barrier system.**

### 4.9. Sterilization

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended.

Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for reprocessing of FITBONE™ reusable instruments.

Use a validated, properly maintained and calibrated steam sterilizer.

The steam quality must be appropriate for the process to be effective.

Do not exceed 140°C (284°F).

Do not stack trays during sterilization.

Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	-	-	Not for use in US	WHO guidelines Not for use
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minute	4 minute	3 minute	18 minute
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated but is suggested only when no other options are available. The Gravity cycle was only validated for wraps and was not validated for sterilization in rigid containers.

### 4.10. Storage

After sterilization, the wrapped reusable instruments should be stored in a location with restricted access.

- This location should be well ventilated and offer sufficient protection from dust, moisture, insects, vermin and extreme temperatures.

**Before the products can be stored, they must be allowed to cool-down to ambient temperature after sterilization.**

**Before opening, the sterile wrap must be thoroughly examined for any damage and the expiry date must be checked.**

## 5. Responsibilities of a hospital for hired sets

Orthopedic surgical instruments generally have a long lifespan, however, incorrect handling or insufficient protection can quickly shorten this service life. Instruments or locking screws which no longer work properly due to an excessively long period of use, incorrect handling or improper maintenance should be disposed of.

Before returning hired sets to Orthofix Srl, they should be decontaminated, cleaned, disinfected, inspected and then sterilized.

Along with the returned instruments, you should also present Orthofix Srl with cleaning and sterilization evidence (PO2100-05). If instruments or locking screws from hired sets are missing or damaged, please inform the manufacturer or the distributor. This is the only way to establish that a complete and functioning set of instruments with locking screws will be available for the next operation.

The instructions contained in this handbook were validated by Orthofix Srl.

It is the responsibility of the end user to ensure that

1. the treatment is carried out with suitable equipment and materials and
2. the individuals tasked with these activities have been trained accordingly.

The equipment, as well as the processes, should be validated and monitored routinely.

To prevent any potential unwanted consequences, all deviations in processing should be properly evaluated for their effectiveness by the responsible individuals.

## 6. Disclaimer

The instructions provided herein have been validated by Orthofix Srl for the initial preparation of a device for first clinical use or for reprocessing of multiple use devices prior to re-use. It remains the responsibility of the individual responsible for reprocessing to ensure that the reprocessing, performed using equipment, materials and personnel in the initial preparation / reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

## 7. Notice about serious incidents

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

## 8. List of instruments with material numbers

Instruments with material numbers	
Code	Description
60000392	T-handle
60000398	Spiral drill D4.5 L300
60000400	Drill sleeve D4.5 black
60000402	Drill sleeve D8.0 green
60000403	Trocar D4.5
60000405	Trocar D4.5, T-handle
60000384	Screw holder M2 L325
60000406	Screwdriver SW3.5, cannulated
60000576	Screwdriver SW3.5
60000317	Tapping tool
60000688	Test pin D4.5
60000689	Open-end wrench SW14/17
60000408	Depth gauge
60000411	Reamer D8.0 L200 L480, rounded
60000412	Reamer D9.0 L100 L480, front-cutting
60000413	Reamer D9.0 L200 L480, rounded
60000414	Reamer D10.0 L100 L480, front-cutting
60000415	Reamer D10.0 L200 L480, rounded
60000416	Reamer D10.5 L200 L480, rounded
60000417	Reamer D11.0 L100 L480, front-cutting
60000418	Reamer D11.0 L200 L480, rounded
60000833	Reamer D11.5 L200 L480, rounded
60000419	Reamer D12.0 L100 L480, front-cutting
60000420	Reamer D12.0 L200 L480, rounded
60000716	Reamer D12.5 L200 L480, rounded
60000421	Reamer D13.0 L100 L480, front-cutting
60000422	Reamer D13.0 L200 L480, rounded
60000423	Reamer D13.5 L200 L480, rounded
60000834	Reamer D8.0 L200 L700, rounded
60000835	Reamer D9.0 L100 L700, front-cutting
60000836	Reamer D9.0 L200 L700, rounded
60000837	Reamer D10.0 L100 L700, front-cutting
60000838	Reamer D10.0 L200 L700, rounded
60000839	Reamer D10.5 L200 L700, rounded
60000840	Reamer D11.0 L100 L700, front-cutting
60000841	Reamer D11.0 L200 L700, rounded

Instruments with material numbers	
Code	Description
60000842	Reamer D11.5 L200 L700, rounded
60000843	Reamer D12.0 L100 L700, front-cutting
60000844	Reamer D12.0 L200 L700, rounded
60000845	Reamer D12.5 L200 L700, rounded
60000846	Reamer D13.0 L100 L700, front-cutting
60000847	Reamer D13.0 L200 L700, rounded
60000848	Reamer D13.5 L200 L700, rounded
60001184	Mounting bracket TAA
60001175	Drill guide TAA
60000175	Retaining screw M6 L9
60000003	Locking screw
60000218	Spacer TAA
60000219	Clamping nut TAA
60000310	Connection screw TAA
60001439	Spot-film device 45°/90°
60001307	Wire trap 90°
60001415	Step reamer TAA1160
60001179	Step reamer TAA1180
60000832	Dummy TAA1160-T-225
60000822	Dummy TAA1180-F-245
60000425	Reamer D14.0 L200 L480, rounded
60000426	Reamer D15.0 L200 L480, rounded
60000513	Retaining screw M4 L400
60000410	Pliers for screws
60001014	Tube T14/13-M
60001015	Tube T13/12-M
60001016	Tube T12/11-M
60001017	Tube T12/10-M
60001018	Tube T12/09-M
60001019	Tube T12/08-M
60001020	Tube T12/11-L
60001021	Tube T12/10-L
60001022	Tube T12/09-L
60001023	Tube T12/08-L
60001024	Tube T12/11-XL
60001025	Tube T12/10-XL

## Instruments with material numbers

Code	Description
60001026	Tube T12/09-XL
60001027	Tube T12/08-XL
60001028	Cone C13
60001029	Cone C13+
60001030	Cone C13++
60001031	Cone C12
60001032	Cone C11
60001033	Tube sinker TS13
60001034	Tube sinker TS12
60001035	Tube sinker TS11
60001036	Cone sinker CS15-13
60001037	Cone sinker CS12-11
60001038	Clamp
60001039	K-wire D3 L280
60001044	Tube T14/13-S
60001045	Tube T13/12-S
60001046	Tube T12/11-S
60001047	Tube T12/10-S
60001048	Tube T12/09-S
60001049	Tube T12/08-S
60001050	Tube T16/15-M
60001051	Tube T15/14-M
60001052	Tube T16/15-L
60001053	Tube T15/14-L
60001054	Tube T14/13-L
60001055	Tube T13/12-L
60001056	Tube T16/15-XL
60001057	Tube T15/14-XL
60001058	Tube T14/13-XL
60001059	Tube T13/12-XL
60001060	Cone C15
60001061	Cone C15+
60001062	Cone C15++
60001063	Cone C14
60001064	Tube sinker TS15
60001065	Tube sinker TS14
60000849	Reamer D14.0 L200 L700, rounded
60000850	Reamer D15.0 L200 L700, rounded

## Instruments with material numbers

Code	Description
60001139	Dummy TAA1160-F-225
60001185	Drill guide blocking screw TAA11
60001244	Dummy TAA1140-T-205
60001248	Dummy TAA1140-F-205
60001464	GRID plate
60001495	Dummy TAA1180-T-245
60001528	Step reamer TAA1140
60001605	Spiral drill D4.5 L300, conical
60001623	Dummy TAA1380-F-245
60001636	Impact sleeve
60001888	Cone C13 retropatellar
60001854	Dummy TAA0960-F-220
60001855	Dummy TAA0960-T-220
60001925	Dummy TAA0940-F-200
60001927	Dummy TAA0940-T-200
60001848	Spiral drill D4.0 L305
60001849	Step reamer TAA0960
60001938	Step reamer TAA0940

### 9. List of screws with material numbers

Screws with material numbers	
Code	Description
60000572	Locking screw D4.5 L80, long
60000571	Locking screw D4.5 L75, long
60000570	Locking screw D4.5 L70, long
60000569	Locking screw D4.5 L65, long
60000568	Locking screw D4.5 L60, long
60000567	Locking screw D4.5 L55, long
60000369	Locking screw D4.5 L50, long
60000368	Locking screw D4.5 L45, long
60000367	Locking screw D4.5 L40, long
60000366	Locking screw D4.5 L35, long
60000365	Locking screw D4.5 L30, long
60000364	Locking screw D4.5 L25, long
60000363	Locking screw D4.5 L20, long
60000272	Locking screw D5.8 L80
60000271	Locking screw D5.8 L75
60000270	Locking screw D5.8 L70
60000269	Locking screw D5.8 L65
60000268	Locking screw D5.8 L60
60000267	Locking screw D5.8 L55
60000266	Locking screw D5.8 L50
60000265	Locking screw D5.8 L45
60000264	Locking screw D5.8 L40
60000263	Locking screw D5.8 L35
60000262	Locking screw D5.8 L30
60000260	Locking screw D4.5 L50, short
60000259	Locking screw D4.5 L45, short
60000258	Locking screw D4.5 L40, short
60000257	Locking screw D4.5 L35, short
60000256	Locking screw D4.5 L30, short
60000255	Locking screw D4.5 L25, short
60000254	Locking screw D4.5 L20, short
60000503	Hollow screw 3.5x40mm VG
60001828	Locking screw D4.0 L20, short
60001829	Locking screw D4.0 L25, short
60001830	Locking screw D4.0 L30, short
60001831	Locking screw D4.0 L35, short

## 10. Interpretation of symbols


Symbol		Description	
		CAUTION: Consult instructions for use for important cautionary information	Consult instructions for Use
		Single use. Do not reuse	
		Non-sterile	
		Catalogue number	Lot number
		CE marking in conformity to applicable European Directives/Regulations	
		Date of manufacture	Manufacturer
		Do not use if package is opened or damaged	
		Symbol for MR conditional. The item has been demonstrated to pose no known hazards in a specified MRI environment with specified conditions of use.	
<b>Rx Only</b>		CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.	

Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9, 37012 Bussolengo  
(Verona), Italy  
Telephone +39 045 6719000,  
Fax +39 045 6719380  
[www.orthofix.com](http://www.orthofix.com)

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Šajā rokasgrāmatā ietvertās instrukcijas apstiprinājis uzņēmums Orthofix Srl.