<p>ORTHOFIX® LIMB RECONSTRUCTION SYSTEM</p> <p>2</p>	EN	<p>ORTHOFIX® SYSTEM FÖR ÅTERUPPBYGGNAD AV LEMMAR</p> <p>22</p>	SV
<p>ORTHOFIX® SISTEMA DI RICOSTRUZIONE DEGLI ARTI</p> <p>4</p>	IT	<p>ORTHOFIX® ΣΥΣΤΗΜΑ ΑΠΟΚΑΤΑΣΤΑΣΗΣ ΜΕΛΩΝ</p> <p>24</p>	EL
<p>LE SYSTÈME DE RECONSTRUCTION DE MEMBRE ORTHOFIX®</p> <p>6</p>	FR	<p>ORTHOFIX® 患肢再生システム</p> <p>26</p>	JA
<p>ORTHOFIX® SYSTEM ZUR REKONSTRUKTION DER GLIEDMASSEN</p> <p>8</p>	DE	<p>ORTHOFIX® 肢体重建外固定支架</p> <p>28</p>	ZH
<p>ORTHOFIX® SISTEMA DE RECONSTRUCCIÓN DE MIEMBROS</p> <p>10</p>	ES	<p>ORTHOFIX® SYSTÉM ZEVNÍCH FIXATÉRŮ KONČETIN</p> <p>30</p>	CS
<p>ORTHOFIX® SYSTEM TIL REKONSTRUKTION AF LEMMER</p> <p>12</p>	DA	<p>ORTHOFIX® SYSTEM DO REKONSTRUKCJI KOŃCZYN</p> <p>32</p>	PL
<p>ORTHOFIX® RAAJAN REKONSTRUKTIOJÄRJESTELMÄ</p> <p>14</p>	FI	<p>ORTHOFIX® SISTEM ZA REKONSTRUKCIJO UDOV</p> <p>34</p>	SL
<p>ORTHOFIX® SYSTEM FOR REKONSTRUKSJON AV LEMMER</p> <p>16</p>	NO	<p>ORTHOFIX® 사지 재건 시스템</p> <p>36</p>	KO
<p>ORTHOFIX® SYSTEEM VOOR LEDEMATENRECONSTRUCTIE</p> <p>18</p>	NL	 <p>ORTHOFIX®</p> <p>ORTHOFIX Srl Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy Telephone 0039-0456719000 - Fax 0039-0456719380</p>	
<p>ORTHOFIX® SISTEMA DE RECONSTRUÇÃO DE MEMBROS</p> <p>20</p>	PT	 0123	<p>PQ LRS ADV E 04/19</p>

**Keep this leaflet: you may need to read it again. Discard previous versions.
Important information - Please read prior to use**

ORTHOFIX[®]

LIMB RECONSTRUCTION SYSTEM

CE 0123

Attention, this leaflet is an addendum of instructions for use PQ EXF.

ADV Ring Hinge (53570)

1. This clamp has not been designed for lengthening procedures and should never be subjected to the forces of lengthening.
2. It may be used in gradual or acute angular corrections as long as a minimum of 3 screws are applied to the ring to guarantee rotational stability, with two screws inserted above and 1 below the ring.
3. The two outer screws should be inserted at an angle close to but no greater than 90°. This means inserting each screw into hole number 6 in a 150mm ring, or 5 in all greater ring sizes, on each side of the hole number 0 above the centre of rotation of the hinge.
4. It can also be used with a minimum of 3 tensioned Kirschner wires with a crossing angle of at least 60 degrees, but only in cases of acute correction, as long as reinforcement bars are added to the frame after the correction.

ADV T-Garches Clamp (53031)

1. It must be used always with 3 bone screws for sufficient stability.
2. This clamp must be positioned at a distance of 20mm from the bone (never greater than 30mm). Use of the spacer is necessary for correct placement of the clamp.
3. During the lengthening procedure, the patient should be allowed to weightbear no more than 30% of the body weight.
4. This clamp should not be used for correcting procurvatum/recurvatum femoral deformities.

ADV Knee Hinge (53590)


This clamp is used in association with a LRS ADV rail and either a Sheffield or TrueLok ring. It permits a similar range of knee motion to that of the Anterior and Posterior Cruciate Ligaments, from 0° to 90° in the sagittal plane. This range of motion can be limited and/or blocked when appropriate.

The knee hinge is indicated for bone stabilization in:

- Ligamentous disruption and instability
- Severe Tibial plateau injuries
- Articular Reconstruction of the knee joint with osteochondral grafting
- Temporary cross knee fixation following removal of infected knee implants

Warnings

1. The Knee Hinge must be correctly aligned with the Reference Axis of the knee joint. In order to achieve this, a K-wire should be inserted into the Reference Axis of the femur (see illustration below) and the Knee Hinge applied over it through the 2mm centering hole positioned on its proximal part.
2. Before stabilising the femur with bone screws, check under image intensification that flexion-extension of the knee is not impeded. This will confirm correct positioning of the hinge. If necessary, replace the K-wire correctly in the Reference Axis.
3. Once all bone screws have been inserted, check again the R.O.M. and if not complete, tighten the posterior locking screw to limit the hinge movement to the R.O.M.
4. In case knee movement has to be blocked, do not tighten the posterior locking screw and use a compression-distraction unit or reinforcement bars to prevent knee movement.
5. In articular fractures limit weightbearing as clinically indicated and tolerated. The surgeon should decide the amount of weightbearing according to the stability of the repair and the rigidity required to ensure congruent healing.

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