ORTHOFIX® GALAXY UNYCO DIAPHYSEAL TIBIA KIT 2				EN		ORTHOFIX® KIT PARA TÍBIA DIAFISÁRIA GALAXY UNYCO 60					
ORTHOFIX® KIT PER TIBIA DIAFISARIA GALAXY UNYCO 9		IT				ORTHOFIX® GALAXYS UNYCO-KIT FÖR DIAFYSÄRA TIBIAFRAKTURER 67					
ORTHOFIX® KIT UNYCO GALAXY POUR TIBIA DIAPHYSISAIRE 15	FR					ORTHOFIX® KIT ΔΙΑΦΥΣΗΣ ΚΝΗΜΗΣ GALAXY UNYCO 73					
ORTHOFIX® GALAXY UNYCO DIAPHYSEN-TIBIA-KIT 22	DE					オーソフィックス® GALAXY UNYCO 膝骨骨幹部用キット 80					
ORTHOFIX® KIT PARA FRACTURA DIAFISARIA DE TIBIA GALAXY UNYCO 28	ES					ORTHOFIX® GALAXY UNYCO 骨干胫骨支架 87					
ORTHOFIX® GALAXY UNYCO-SÆT TIL TIBIADIAFYSER 35	DA					SOUPLAVA PRO DIAFYZÁRNÍ ZLOMENINY HOLENNÍ KOSTI ORTHOFIX® GALAXY UNYCO 93					
ORTHOFIX® GALAXY UNYCO SÄÄRILUUN VARSIOSAN KIINNITYSPAKKAUS 41	FI	DA				ORTHOFIX® ZESTAW TRZONU KOŚCI PISZCZELOWEJ GALAXY UNYCO 99					
ORTHOFIX® GALAXY UNYCO DIAPHYSEAL TIBIA KIT 47	NO	NO				ORTHOFIX® KOMPLET GALAXY UNYCO ZA DIAFIZNI DEL GOLENICE 106					
ORTHOFIX® GALAXY UNYCO DIAFYSaire TIBIA KIT 53	NL					ORTHOFIX® GALAXY UNYCO 경골 간부 키트 112					

Orthofix Srl

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Telephone +39 045 671 9000 - Fax +39 045 671 9380

DESCRIPTION

The Orthofix Galaxy UNYCO Diaphyseal Tibia Kit consists of a series of UNYCO Screws, UNYCO Cancellous Screws, Large Multiscrew Clamp for UNYCO Screws, Rod Ø 12 mm L 350 mm and specific application tools. External fixation systems are modular, therefore different frame configurations are possible. The Orthofix components in the Galaxy UNYCO Diaphyseal Tibia Kit are not intended to replace normal healthy bone or to withstand the stresses of weight bearing. The System may be used in conjunction with the Orthofix Galaxy Fixation system for crossing the knee and the ankle during temporary fixation. All Orthofix devices are intended for professional use only. Surgeons who supervise the use of Orthofix devices must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the Orthofix modular system. To promote the proper use of its fixation system, and establish an effective promotional and training tool, Orthofix has developed several manuals or CD-ROM's containing the relevant information (i. e. general philosophy, surgical application, etc.) called "Operative Techniques". These are available in several languages as a complimentary service for surgeons who have adopted the Orthofix system. If you wish to receive a personal copy, please contact Orthofix or its local authorized representative, with a description of the medical device to be used. They are also available in the Orthofix website and Orthofix mobile apps.

INTENDED USE

The Galaxy UNYCO Diaphyseal Tibia Kit is intended to be used to provide bone stabilization in trauma, specifically lower limb fractures, that require temporary fixation prior to definitive treatment.

INDICATIONS

The indications for use include:

- Tibial fractures extending from about 8cm below the knee to about 7cm above the ankle joint, including comminuted open or closed tibial fractures and polytrauma
- Temporary stabilisation of the tibia after debridement for osteomyelitis or an infected non-union pending second stage treatment

The product is indicated for non-weight-bearing use.

NOTE: UNYCO Screws and UNYCO Cancellous Screws must always be used in association with the Large Multiscrew Clamp for UNYCO Screws in a pin cluster. If faced with a segmental fracture, the most proximal and most distal bone segments are held using the Large Multiscrew Clamp for UNYCO Screws and the intermediate segment can be held using a UNYCO Screw in a Galaxy Large Clamp attached to the same connecting rod that links the two Large Multiscrew Clamps for UNYCO Screws.

CONTRAINDICATIONS

The Orthofix Galaxy UNYCO Diaphyseal Tibia Kit is not designed or sold for any use except as indicated.

Use of the system is contraindicated in the following situations:

- Patients with mental or physiological conditions who are unwilling or incapable of following postoperative care instructions
- Patients with severe osteoporosis
- Patients with malignancy in the fracture area
- Patients with foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implant insertion

WARNINGS & PRECAUTIONS

1. Fracture stabilization has to be done following correct fracture reduction.
2. The arms of the Large Multiscrew Clamp for UNYCO Screws can be disassembled to allow the most suitable configuration (M or U) while the other components cannot be disassembled.
3. Before fixator application, make sure that the Large Multiscrew Clamps for UNYCO Screws are fully loosened.
4. UNYCO Screws must be inserted at right angles to the bone surface using a steady force to ensure optimal bone purchase.
5. Once the first screw has been inserted freehand and its correct position on the bone has been checked, apply the Galaxy UNYCO Tibia Clamp (93566) on the first screw and tighten the metal ring on the arm clockwise.
6. The system stability is guaranteed only with a minimum of 3 UNYCO screws (2 of them must be UNYCO screws 93507) coupled with the Large Multiscrew Clamp for UNYCO Screws in each segment. The use of four UNYCO screws is to be preferred. The Large Multiscrew Clamp for UNYCO Screws should never be applied across the fracture line.
7. Frame stability must be checked intra-operatively before the patient leaves the operating theatre. If one UNYCO screw is perceived as unstable, re-insert it in a new bone area, using the 10° freedom of movement to find the new insertion point.
8. Depending upon the clinical and radiological findings, the surgeon will decide on the number of UNYCO Screws (93507) and UNYCO Cancellous Screws (93508) needed to achieve the appropriate frame stability. Remember that the first screw determines the position of the remaining screws in the Large Multiscrew Clamp for UNYCO Screws with ±10° of freedom, therefore safe corridors have to be considered in advance for correct and safe screw insertion. Note: Once converging UNYCO screws have been inserted, the Large Multiscrew Clamps for UNYCO Screws can no longer slide on the screw shafts. It is therefore important to determine the final distance of the fixator from the skin before inserting the second screw in the clamp unless screws are inserted parallel to each other.
9. UNYCO Screws and UNYCO Cancellous Screws should be selected in accordance with bone type (cortical or cancellous).
10. At the end of surgery, check screw insertion both in the AP and lateral planes, using image intensification.
11. The UNYCO screw may protrude beyond the first cortex or sit close to the insertion point of a nail. In this case pay attention whenreaming or inserting a nail at a later stage.
12. Compression is never recommended in a fresh fracture.
13. Particular care should be taken that screws do not enter the joints.
14. The kit is for "single use only": therefore **IT MUST NOT BE RE-USED.**
15. For more stable fixation of a fracture with a fixator, we recommend that the nearest bone screw is applied fairly close to the fracture margin (a minimum of 1 cm is recommended) and that these distances are equal on both sides of the fracture.
16. UNYCO Screws and UNYCO Cancellous screws must be applied ONLY by drilling the screw using the Power Drill Torque Limiter (99-93568) inserted into the power drill or by hand using the Limited Torque Wrench (99-93567). When using the Limited Torque Wrench, it should be held steady so that the direction is maintained throughout the procedure. During insertion of UNYCO Screws stop drilling when the Power Drill Torque Limiter stops turning and a rattle noise is heard. During insertion of UNYCO Cancellous Screw stop drilling when its soft tissue reference line reaches the skin or, when the Power Drill Torque Limiter stops turning, whichever comes first. The power drill must be used at its lowest speed. The Power Drill Torque Limiter (99-93568) is provided in the sterile kit; the Limited Torque Wrench (99-93567) is available upon request. It is mandatory NOT to apply UNYCO screws in cancellous bone so it is advisable to use always UNYCO Cancellous screws close to the joint.
17. All equipment in the Galaxy UNYCO Diaphyseal Tibia Kit (99-93506) should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.

18. The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. It is recommended to apply the Large Multiscrew Clamp for UNYCO Screws 4 cm from the bone.
19. Additional equipment may be required for fixation application and removal.
20. Once two or more UNYCO screws have been inserted in the Large Multiscrew Clamp for UNYCO Screws, the latter should not be moved along the screw shafts. Sliding the clamp might cause loss of bone purchase.
21. In case of knee and/or ankle spanning, stabilization in the femur and foot must be performed with bicortical screws in conjunction with Galaxy Fixation System.
22. Screw and frame stability should be monitored at regular intervals. Weight-bearing or excessive loading may lead to loosening or disassembly.
23. If the system is perceived as impediment for application of the definitive treatment, remove the system parts where needed.
24. Meticulous screw site hygiene is required.
25. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
26. Removal of the device: the surgeon should make the final decision when to remove the device and replace it with definitive fixation.
27. Do not use components of the Orthofix Galaxy UNYCO Diaphyseal Tibia Kit in conjunction with products of other manufacturers, unless otherwise specified, as the combination is not covered by the necessary validation.

POSSIBLE ADVERSE EFFECTS

1. Nerve or vessel damage resulting from insertion of screws.
2. Superficial or deep bone screw tract infection, osteomyelitis, or septic arthritis, including chronic drainage of bone screw sites after device removal.
3. Oedema or swelling; possible compartment syndrome.
4. Loosening or breakage of implants.
5. Bony damage due to inappropriate implants selection or application.
6. Persistence or recurrence of the initial condition requiring treatment.
7. Foreign body reaction to implants or frame components.
8. Tissue necrosis secondary to implants insertion.
9. Pressure on the skin caused by external components when clearance is inadequate.
10. Limb length discrepancy.
11. Excessive operative bleeding.
12. Intrinsic risks associated with anesthesia.
13. Intractable pain.
14. Vascular disorders including thrombophlebitis, pulmonary embolus, wound hematomas, avascular necrosis.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Important

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require further surgical intervention to remove or replace the external fixation device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the external fixation devices are important considerations in the successful utilization of Orthofix external fixation devices by the surgeon. Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindications or is predisposed to any contraindications, DO NOT USE Orthofix Galaxy UNYCO Diaphyseal Tibia Kit.

Materials

The Orthofix Galaxy UNYCO Diaphyseal Tibia Kit is comprised of stainless steel, aluminum alloy, titanium alloy and plastic components. Those components which contact the patient are the percutaneous pins (UNYCO Screws and UNYCO Cancellous screws). These are manufactured from surgical grade stainless steel.

Sterile

Devices or kits provided STERILE are labeled as such.

Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.

IMPORTANT: The Power Drill Torque Limiter and the Limited Torque Wrench 99-93567 must be discarded if package is opened or damaged.

Non-Sterile

Unless otherwise noted, Orthofix external fixation components are provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Product integrity and performance are assured only if packaging is undamaged.

Cleaning and Sterilization

WARNINGS Aluminum based instruments are damaged by alkaline ($\text{pH} > 7$) detergents and solutions. Anodised coating is damaged by detergents with free halogen ions or sodium hydroxide. Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.

Limitations and restrictions on reprocessing PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED. Repeated reprocessing has minimal effect on reusable instruments. End of life is normally determined by wear and damage due to use.

INSTRUCTIONS FOR PROCESSING NEW DEVICES SUPPLIED "NON-STERILE" PRIOR TO THEIR FIRST USE

General

- Unless supplied sterile, all Orthofix medical devices must be sterilized prior to surgical use.
- A new product means any device taken out of its original Orthofix packaging.

Recommended decontamination process

The following sequence of processes is recommended for rendering Orthofix UNYCO devices which are supplied non-sterile safe for their first clinical use:

1. CLEANING
2. STERILIZATION

Preparation for cleaning of new products

- When applicable, the devices should be disassembled completely.
- Wherever possible, all parts of disassembled devices should be kept together in one container.

STEP 1: Cleaning of new products

- REMOVE products from their original packaging. All equipment should be carefully examined prior to use to assure proper working condition.
- THOROUGHLY CLEAN with a woven-non woven tissue soaked using a solution of 70% medical grade alcohol and 30% distilled water or with compatible detergent. Detergents with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- RINSE with sterile distilled water

STEP 2: Drying

- Carefully hand-dry using absorbent, non-shedding cloth or industrial hot air dryer, or place in a drying cabinet for 20 minutes at 110° C.

STEP 3: Preliminary inspection

- All instruments and product components should be visually inspected for cleanliness and any signs of deterioration, that may cause failure in use (such as cracks or damage to surfaces) and functions tested before being sterilized. Visually inspect for residual dirt/debris. If any dirt/debris remains repeat cleaning process until no dirt/debris can be visually detected.
- If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- When instruments form part of an assembly, check assembly with matching components.

STEP 4: Sterilization

- **PACKAGING:** Where products are to be packaged to maintain their sterility after sterilization and to prevent damage of the instrument prior to use, an appropriate medical grade packaging material should be used. The pack should be large enough to contain the instruments without stressing the seals.
- **STERILIZATION CONTAINERS:** Instruments may be loaded into a dedicated (Orthofix) instrument tray, or general-purpose sterilization tray. Cutting edges should be protected and the recommended content or maximum weight not exceeded as indicated by manufacturer.
- **STERILIZATION:** Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below. Trays must be wrapped using FDA approved sterilization wraps made of a nonwoven blend of cellulose and polyester. Do not stack trays during sterilization.

STEAM STERILIZER TYPE	GRAVITY	PRE-VACUUM
Preconditioning pulses	N/A	4 pulses
Exposure Temperature	132°C (270°F)	132°C (270°F)
Minimum Holding Time	15 minutes	4 minutes
Drying Time	30 minutes	30 minutes

STEP 5: Storage

Store the sterilized instrument in a dry, clean, and dust-free environment at room temperature.

Note: ANY DEVICE WHICH IS LABELLED "SINGLE USE ONLY" MUST NEVER BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR SAFETY AND EFFECTIVENESS FOR THE FIRST PATIENT USE OF SINGLE USE DEVICES. The institution or practitioner bears full responsibility for any subsequent use of these devices.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE IMPLANTABLE DEVICE*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol  reported on the product label.

After the removal from the patient, the implantable device* has to be dismantled.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol  reported on the label or are indicated in the "Instructions For Use" supplied with the products.

The re-use of "SINGLE USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

MRI SAFETY INFORMATION

Galaxy UNYCO Diaphyseal Tibia Kit is labeled MR CONDITIONAL according to the terminology specified in ASTM F2503 Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment.

Non-clinical testing has demonstrated that the Galaxy UNYCO Diaphyseal Tibia Kit components are MR Conditional. It can be safely scanned under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90 mT/cm)
- Maximum whole-body-averaged specific absorption rate (SAR) of $\leq 4.0 \text{ W/kg}$ (First Level Controlled Operating Mode).
- No local transmit/receive coils can must be used on the device.
- The Galaxy UNYCO Diaphyseal Tibia Kit components must be entirely outside the MR scanner bore.
- No part of the Galaxy UNYCO Diaphyseal Tibia Kit components must extend into the MR bore. Therefore MR scanning of body parts where the Galaxy UNYCO Diaphyseal Tibia Kit components are located is contraindicated.

DISPLACEMENT INFORMATION

The system will not present an additional risk or hazard to a patient in the 1.5 and 3.0 Tesla MR environment with regard to translational attraction or migration and torque.

HEATING INFORMATION

Under the scan conditions defined above, the Galaxy UNYCO Diaphyseal Tibia Kit is expected to produce a maximum temperature rise of 2°C after 15 minutes of continuous scanning.

MR PATIENT SAFETY

MRI in patients with Galaxy UNYCO Diaphyseal Tibia Kit components can only be performed under these parameters. It is not allowed to scan the Galaxy UNYCO Diaphyseal Tibia Kit components directly. Using other parameters, MRI could result in serious injury to the patient. When the Galaxy UNYCO Diaphyseal Tibia Kit components are used in conjunction with other External Fixation Systems, please be advised that this combination has not been tested in the MR environment and therefore higher heating and serious injury to the patient may occur. Because higher *in vivo* heating cannot be excluded, close patient monitoring and communication with the patient during the scan is required. Immediately abort the scan if the patient reports burning sensation or pain.

Galaxy UNYCO Diaphyseal Tibia Kit can only be guaranteed for MRI when using the following components to build a frame:

Description
2 x Large Multiscrew Clamps for UNYCO Screws
1 x Rod Ø 12 mm L 350 mm
<i>Up to 8 within these screws available in the kit:</i>
8 x UNYCO Screw QC Shaft Ø 6 mm
4 x UNYCO Cancellous Screw QC Shaft Ø 6 mm

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