# MJ-Flex

The New Metaizeau Nail™

EN	MJ-Flex The New Metaizeau Nail™ 2		MJ-Flex The New Metaizeau Nail™ 86	
IT	MJ-Flex The New Metaizeau Nail™ 9		MJ-Flex The New Metaizeau Nail™ 93	
FR	MJ-Flex The New Metaizeau Nail™ 16		MJ-Flex The New Metaizeau Nail™ 100	
DE	MJ-Flex The New Metaizeau Nail™ 23		MJ-Flex The New Metaizeau Nail™ 107	
ES	MJ-Flex The New Metaizeau Nail™ 30		MJ-Flex The New Metaizeau Nail™ 114	
DA	MJ-Flex The New Metaizeau Nail™ 37		MJ-Flex The New Metaizeau Nail™ 121	
FI	MJ-Flex The New Metaizeau Nail™ 44		MJ-Flex The New Metaizeau Nail™ 128	
NO	MJ-Flex The New Metaizeau Nail™ 51		MJ-Flex The New Metaizeau Nail™ 135	
NL	MJ-Flex The New Metaizeau Nail™ 58		„MJ-FLEX“ NAUJA „METAIZEAU NAIL™“ 142	
PT	MJ-FLEX A NOVA HASTE METAIZEAU™ 65		MJ-Flex The New Metaizeau Nail™ 149	
SV	MJ-Flex The New Metaizeau Nail™ 72		MJ-Flex The New Metaizeau Nail™ 156	
EL	MJ-Flex The New Metaizeau Nail™ 79		MJ-Flex The New Metaizeau Nail™ 163	
	JA	ZH	CS	PL
			SL	KO
			AR	BG
			LT	MS
			RO	TR

ORTHOFIX Srl  
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy  
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

0123 PQMJF D 02/22 (0424538)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online

### **Important information - please read prior to use**

#### **See also instruction leaflet PQRMD for reusable medical devices**

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

#### **DESCRIPTION**

MJ Flex is an intramedullary implant system specifically designed for Elastic Stable Intramedullary Nailing (ESIN) fracture fixation. The principle of ESIN involves balanced nailing to provide elasticity and stability at the fracture site.

#### **INTENDED USE**

The MJ-FLEX THE NEW METAIZEAU NAIL™ is intended for the treatment of diaphyseal fractures of long bones.

#### **INDICATIONS**

The MJ-FLEX THE NEW METAIZEAU NAIL™ is indicated to treat:

- upper extremity and clavicle fractures in all patients except newborns and infants;
- lower extremity fractures in pediatric patients, except newborns and infants, where the flexibility of the implant is paramount not to disrupt the growth plate;
- lower extremity fractures in small adults where the medullary canal is narrow.

#### **NOTES FOR USE**

The product is intended for professional use only. Surgeons who supervise the use of the product must have full awareness of orthopedic fixation procedures and should be familiar with the devices, instruments and surgical procedure, including the application and removal. Detailed operative technique guidance is available on request; please contact Orthofix or your local distributor.

#### **CONTRAINDICATIONS**

Use of the MJ-FLEX THE NEW METAIZEAU NAIL™ is contraindicated for patients with:

- general medical conditions not suitable for surgery;
- active or suspected latent infections in the fracture area;
- foreign body sensitivity (where material sensitivity is suspected, tests should be made prior to implant insertion);
- mental or physiological conditions, who are unwilling or incapable of following postoperative care instructions;
- compromised bone stock, potentially leading to inadequate or unstable fixation of the device;
- fractures with undiagnosed pathologies.

#### **GENERAL WARNINGS**

1. Preoperative procedures, knowledge of surgical technique and proper selection and placement of the implant are important aspects for the successful use of the device by the surgeon;
2. Do not combine MJ-FLEX THE NEW METAIZEAU NAIL™ with components from other systems;
3. Do not combine two nails of dissimilar metals, since this may cause an electrolytic reaction;
4. Never insert more than two nails in the same medullary canal;
5. Nails must not be reused. If any implant has come into contact with any body fluid it should be considered to have been used. Loosening, cracking or breakage of the implants can cause loss of stabilization;
6. This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine;
7. In skeletally immature patients, it is recommended to avoid crossing open epiphyseal plates to decrease the probability of consequent growth disturbance.

#### **GENERAL PRECAUTIONS**

1. Ensure that all components needed for the operation are available in the operating theatre;
2. Examine all components carefully prior to use. Product integrity, sterility (in the case of sterile products) and performance are assured only if the packaging is undamaged;
3. Do not use if packaging is compromised or if a component is believed to be faulty, damaged or suspect;
4. Careful handling and storage of the product is required. Scratching or damage to the component can significantly reduce the strength and fatigue resistance of the products;
5. Check the fit and functioning of nails and instruments assembly prior to implantation;
6. Ideally, two nails should be used. Where the medullary canal is very narrow, (e.g. for clavicular, metacarpal, phalangeal and metatarsal fractures), only one nail might be used taking care to fill the canal as much as possible;
7. Ideally, the two nails should have the same diameter so that the opposing forces are equal, limiting iatrogenic valgus or varus deformities;
8. The diameter of the selected nail should be no more than 50% of the width of isthmus. In case two nails with different diameters are used, the sum of the two diameters should not be more than the width of the isthmus;
9. The flat part of the bended nail must be on the concave side of the curve;

10. Avoid creating sharp bends. The shape of the bent nail must be a bow without any edge;
11. If the nail is bent more than planned, avoid bending it in reverse but consider bending a new nail;
12. When opening the cortex, insert the awl with moderate pressure. The tip of the awl should be screened to avoid inadvertently driving it further than intended. This is essential when the awl is pointing towards potentially critical locations (such as the region of the femoral artery);
13. For optimal control of the nail, make sure that one jaw of the chuck is positioned exactly over the flat part of the nail;
14. During insertion, check the progression of the nail under fluoroscopy to ensure the tip is advancing without any impediment. If the implant is not advancing properly, the surgeon should never persist in the insertion, but should review the situation and consider the following options:
  - make sure the tip of the nail has the proper orientation;
  - increase the contour angle near the tip of the nail;
  - gently hammering the nail;
  - change the nail for the next-smaller diameter;
15. Vigorous hammering should never be necessary. The nail should advance blow by blow. If the nail does not advance properly, the surgeon should never persist with hammering;
16. Make sure not to twist the tip of the second nail more than 180° to avoid the second nail spiral up around the first nail ("corkscrew effect");
17. Under fluoroscopy, verify the final position of the tip of the nail to prevent inadvertent damage to the physis. It is recommended to stop nail insertion at least 1cm before the physis to allow adequate clearance for optional final impaction of the nail;
18. Before cutting the nail, release traction and confirm bone rotation and bone alignment under fluoroscopy in two planes to avoid malalignment;
19. Improper positioning of the nails may result in loosening, cracking, or fracture of the device or bone or both;
20. Multiple re-positioning of the nail may lead to poor anchoring of the nail tip;
21. Ensure to select the appropriate cutter and cutter hole according to the selected diameter of the nail;
22. Ensure to select the appropriate impactor according to the selected diameter of the nail;
23. If callus is slow to or fails to develop, the implant may eventually break due to metal fatigue:
  - to avoid this, further measures may be required, including replacement of the implant;
24. When appropriate the implants should be removed at conclusion of treatment.

#### **POSSIBLE ADVERSE EFFECTS**

1. Early or late infection, both deep and superficial.
2. Delayed union, non-union, malunion, malalignment, increased fibrous tissue response.
3. Stiffness, loss of range of motion, limp, leg length discrepancy.
4. Intrinsic risks associated with anesthesia and surgery, including, but not limited to:
  - Soft tissue irritation and/or inflammation;
  - Nerve damage;
  - pain, heterotopic ossification at the surgery site;
  - vascular disorders such as compartment syndrome, thromboembolic events;
  - Damage to the growth plate with consequent growth disturbance.
5. Complications associated with metal sensitivity.
6. Hardware migration/displacement with or without skin or joint perforation.
7. Loss of reduction.
8. Refracture.

#### **IMPORTANT**

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the devices by the surgeon.

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE the ORTHOFIX MJ Flex Nail.

#### **MRI SAFETY INFORMATION**

ORTHOFIX MJ Flex Nail has not been evaluated for safety and compatibility in the MR environment.

It has not been tested for heating, migration, or image artifact in the MR environment. The safety of The ORTHOFIX MJ Flex Nail in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### **STERILE & NON-STERILE Product**

Orthofix provides certain devices STERILE while others are provided NON-STERILE. Please review the product label to determine the sterility of each device.

#### **Sterile**

Devices or kits provided STERILE are labeled as such. Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.

#### **Non-Sterile**

Unless otherwise noted, Orthofix external fixation components are provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Product integrity and performance are assured only if packaging is undamaged.

## INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

### Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

### Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

### POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

### CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

### PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

#### Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

### CLEANING

#### General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

## **Manual cleaning**

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

## **Manual disinfection**

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

## **Automatic cleaning and disinfection using washer-disinfector**

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
  - a. Cannulations
  - b. Long blind holes
  - c. Mating surfaces
  - d. Threaded components
  - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
  - a. Pre-cleaning for 4 minutes;
  - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
  - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
  - d. Final rinsing with deionized water for 3 minutes;
  - e. Thermal disinfection at least 90°C or 194°F (max 95 °C or 203° F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
  - f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.
- The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfecter when it completes the cycle.
13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

## MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

## PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

## STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AAMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

## STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

## DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

## **CLEANING AGENT INFORMATION**

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym  
concentration 2%
- For manual cleaning: Neodisher Mediclean  
concentration 2%
- For automated cleaning: Neodisher Mediclean  
concentration 0.5%

## **RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE**

### **IMPLANTABLE DEVICE\***

The "SINGLE USE" implantable device\* of Orthofix is identified through symbol “⊗” reported on the product label. After the removal from the patient, the implantable device\* has to be discarded.

The re-use of implantable device\* introduces contamination risks for users and patients.

The re-use of implantable device\* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

### **NON IMPLANTABLE DEVICE**

The "SINGLE USE" non implantable device of Orthofix is identified through symbol “⊗” reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

### **CAUTION**

Federal Law (USA) restricts this device to sale by or on the order of a physician.

### **Manufacturer contact**

Please contact your local Orthofix sales representative for further details.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description	
 	Consult instructions for Use	CAUTION: Consult instructions for use for important cautionary information
	Single use. Do not reuse	
	STERILE. Sterilised by irradiation	
	NON STERILE	
 	Catalogue number	Batch code
	Expiry date (year-month - day)	
 	CE marking in conformity to applicable European Directives/Regulations	
 	Date of manufacture	Manufacture
	Do not use if package is opened or damaged	
<b>Rx Only</b>	CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.	

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Orthofix Srl