JA	PREFIX™ FIXATOR 110	LV	PREFIX™ FIKSATORS 231
EL	ΣΤΑΘΕΡΟΠΟΙΗΤΗΣ PREFIX™ 101	HR	PREFIX™ FIKSATOR 222
SV	PREFIX™ FIXATOR 92	TR	PREFIX™ FIXATOR 213
PT	FIXADOR PREFIX™ 83	RO	FIXATORUL PREFIX™ 204
NL	PREFIX™ FIXATOR 74	MS	PENETAP PREFIX™ 195
NO	PREFIX™ FIKSATOR 65	LT	FIKSATORIUS „PREFIX™“ 186
FI	PREFIX™ FIXATOR 56	BG	PREFIX™ ФИКСАТОР 177
DA	PREFIX™ FIXATOR 47	AR	PREFIX™ مثبت 167
ES	FIJADOR PREFIX™ 38	KO	PREFIX™ 고정장치 157
DE	PREFIX™ FIXATEUR 29	SL	FIKSATOR PREFIX™ 148
FR	FIXATEUR PREFIX™ 20	PL	STABILIZATOR PREFIX™ 139
IT	FISSATORE PREFIX™ 11	CS	FIXÁTOR PREFIX™ 130
EN	PREFIX™ FIXATOR 2	ZH	PREFIX™ 固定器 120

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PQPFX E 05/21 (0420400)

EN

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online**Important information - please read prior to use****See also instruction leaflet PQSCR (version D or higher) for implantable devices and related instruments and PQRMD for reusable medical devices****PREFIX™ FIXATOR**

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

MEDICAL DEVICE INFORMATION**DESCRIPTION**

The PREFIX™ Fixator (PREFIX hereinafter) is a monolateral external fixator consisting of a series of components that build the external frame. The external frame is connected to the bone by means of bone screws. Application and removal of the PREFIX can be performed with Orthofix general orthopedic instrumentation.

INTENDED PURPOSE AND INDICATIONS**INTENDED PURPOSE**

The PREFIX is intended to provide bone fixation.

INDICATIONS FOR USE

The PREFIX is indicated for temporary fixation of:

- Tibia fractures, joint dislocations and soft tissue defects in adult and adolescent patients
- Ankle fractures and soft tissue defects in adult patients

CONTRAINDICATIONS

DO NOT USE the PREFIX if a surgical candidate exhibits or is predisposed to any of the following contraindications:

- Mental or physiological conditions who are unwilling or incapable of following postoperative care instructions
- Severe osteoporosis
- Malignancy in the fracture area
- Suspected or documented metal sensitivity reactions

as it could result in a treatment failure in the intended population.

INTENDED PATIENTS

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. The PREFIX is intended for adult and adolescent patients as specified in the indications of use.

INTENDED USERS

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal).

NOTES FOR USE - IMPLANT REMOVAL

Once the temporary treatment with external fixation is complete, implant must be removed and replaced with definitive treatment. HCP should consider premature removal in case of adverse events.

DISCLAIMER

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

MATERIAL

The implants are made from an implant grade material that is specified on the product label.

WARNINGS

- Compression is never recommended in a fresh fracture.
- All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. If the fixator is sited at a distance of more than 4cm from the bone the surgeon will decide on the number of rods and bone screws needed to achieve the appropriate frame stability.
- Components may not be interchangeable between all Orthofix external fixation systems. Consult individual operative technique guides for interchangeable components.

- This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
- When using the Multi-Screw Clamp, two rods must always be used to increase frame stiffness.

PRECAUTIONS

- The surgeon must evaluate the integrity of the construct at follow up visits.
- During and after insertion, ensure correct positioning of the implants under image intensification.
- Fracture stabilization has to be done following correct fracture reduction.
- Before fixator application, make sure that the clamps are fully loosened.
- Frame stability must be checked intra-operatively before the patient leaves the operating theatre.

All Orthofix devices should be used together with their corresponding Orthofix implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the PREFIX when used in conjunction with devices of other manufacturers or with other Orthofix devices if not specifically indicated in the Operative Technique.

POSSIBLE ADVERSE EVENTS

- Damage to surrounding tissues due to surgical trauma
- Superficial infection
- Deep infection
- Compartment syndrome
- Bone fracture during treatment
- Bending, breakage or migration of the device
- Pain, discomfort or abnormal sensations due to the presence of the device
- Wound healing complications
- Complex Regional Pain Syndrome
- Events caused by intrinsic risks associated with anesthesia and surgery

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP.

MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

The components of the PREFIX, listed in Table 1, have been designed as MRI conditional up to 1.5 Tesla. The MRI components have been tested according to ASTM Standards F2052, F2182, F2213 and F2119.

The components of the PREFIX marked with "MR" are MRI conditional. PREFIX can only be guaranteed for MRI when using Orthofix XCALIBER™ bone screws and the following rods to build a frame:

Table 1 – MRI conditional components of PREFIX

CLAMPS

Code	Description
92012	Blue Rod Clip Clamp
92016	Blue Screw Clip Clamp
92026	Blue Transfixing Clip Clamp
92031	Blue Multi Screw Clamp

RODS

Code	Description
92125	Rod 125mm long
92175	Rod 175mm long
92275	Rod 275mm long
92350	Rod 350mm long

XCALIBER SCREWS

Total Length (mm)	Thread Length (mm)						
	30	40	50	60	70	80	90
150	911530	911540	911550	911560	911570	911580	911590
260	912630	912640	912650	912660	912670	912680	912690

Please note that temperature changes reported below apply to the designated MRI system and characteristics used:

- MRI machine: PHILIPS Achieva 1.5T A-series
- Device Position: isocenter of MRI scanner
- Magnetic Field Strength: 1.5 Tesla
- Spatial gradient field: 9.37mT/m

If a different MRI system is used, temperature changes may vary. Contact is advised with the supplier of the MRI system and Orthofix.

It has been shown by specific MRI tests that the MRI Conditional components of PREFIX may be used for patients undergoing MRI procedures using up to 1.5 Tesla MR systems if certain specific conditions are followed. Three commonly used frames have been tested for MRI use at 1.5 Tesla. The results are as follows:

Blue Multi Screw Clamp

- A maximum temperature increase of 4.0 °C has been measured at the tip of the indicated screws (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)

Code	Description
2x92350	Rod 350mm long
4x911560	Self-Drilling XCaliber Screw, L. 150mm, thread length 60mm
4x92012	Blue Rod Clip Clamp
2x92031	Blue Multi Screw Clamp

Z Configuration


- A maximum temperature increase of 7.0 °C has been measured at the tip of the most external screw (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)

Code	Description
2x92125	Rod 125mm long
1x92275	Rod 275mm long
2x911560	Self-Drilling XCaliber Screw, L. 150mm, thread length 60mm
2x911540	Self-Drilling XCaliber Screw, L. 150mm, thread length 40mm
4x92016	Blue Screw Clip Clamp
2x92012	Blue Rod Clip Clamp

Ankle Configuration

- A maximum temperature increase of 2.1°C has been measured at the tip of the indicated screw (Maximum Whole Body averaged SAR of 2W/kg for 6 minutes of scanning)
- At the indicated clamp a temperature increase of 12.2°C (with reference to room temperature) has been measured

Code	Description
1x92125	Rod 125mm long
1x92175	Rod 175mm long
1x92275	Rod 275mm long
2x911540	Self-Drilling XCaliber Screw, L. 150mm, thread length 40mm
1x92080	Transfixing Pin, thread length 80mm, thread Ø 5mm, shaft Ø 4mm
2x92016	Blue Screw Clip Clamp
2x92012	Blue Rod Clip Clamp
2x92026	Blue Transfixing Pin Clamp

These tests have been performed in areas where the greatest temperature increase is expected with commonly used frames: at the tip and surrounding the tip of the most external screw. Due to the versatility of the system, an unlimited number of frames can be built which makes it impossible to test every construct. There are factors that can influence these results (e.g. the number of screws used). Therefore, it is recommended that each frame be evaluated by a radiologist or MR scientist before the MRI procedure is undertaken to ensure patient safety. Since different frame configurations and frame sizes might lead to higher temperature increases, Orthofix recommends that the SAR settings are minimized as much as possible. None of the components should move or migrate in the 1.5 Tesla MRI environments. Non-clinical testing has not been performed to rule out the possibility of component movement or migration at static magnetic field strengths higher than 1.5 Tesla or maximum spatial gradients higher than 9.37mT/m.

IMAGE ARTIFACTS

MR image quality may be compromised if the area of interest is in the vicinity or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. Representative components used to assemble MRI Conditional PREFIX have been evaluated in MRI chamber. Typical values of artifacts sizes are reported below.

Code	Description		
1x92275	Rod 275mm long		
1x92012	Blue Rod Clip Clamp		
1x92016	Blue Screw Clip Clamp		
1x911560	Self-Drilling XCaliber Screw, L. 150mm, thread length 60mm		
Pulse Sequence		Spin-echo	Gradient-echo
TR (ms)		500	334
TE (ms)		4,8	4,6
Slice Thickness (mm)		5	5
Matrix size		230x288	230x288
Flip Angle		90°	30°
Imaging Plane		perpendicular to static field	perpendicular to static field
Artifact Size (mm)		16	34

Code	Description		
1x92350	Rod 350mm long		
1x92012	Blue Rod Clip Clamp		
1x92031	Blue Multi Screw Clamp		
2x911560	Self-Drilling XCaliber Screw, L. 150mm, thread length 60mm		
Pulse Sequence		Spin-echo	Gradient-echo
TR (ms)		500	334
TE (ms)		4,8	4,6
Slice Thickness (mm)		5	5
Matrix size		230x288	230x289
Flip Angle		90°	30°
Imaging Plane		perpendicular to static field	perpendicular to static field
Artifact Size (mm)		20	46

The components of the PREFIX not marked with "MR" have not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment. They have not been tested for heating, migration, or image artifact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has these devices may result in patient injury.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE

EXPECTED CLINICAL BENEFITS

- Provides temporary external fixation
- Stabilizes fractures definitively by using adequate numbers of bone screws and rods to provide mechanical stability
- Rapid stabilization of complex fractures
- Reduces the fracture in order to restore alignment
- Can be applied to a single bone to stabilize a fracture, or across a joint if this is involved in the injury
- Facilitates damage control of soft tissue
- Positions screws where the condition of the bone and soft tissues permits



PERFORMANCE CHARACTERISTICS OF THE DEVICE

- Mechanical stability is assured by using an adequate numbers of bone screws and rods
- The largest diameter (14 mm) rod available for temporary external fixation
- Excellent torsional stability
- Optimal stability due to large rod and screw diameter


SPECIFIC INFORMATION ON THE DEVICE

- Dynamization and physical therapy guidelines should be followed based on each individual case and the fixation system used, and should be instituted as and when considered appropriate by the surgeon, in accordance with clinical and radiological findings.
- Additional equipment may be required for fixation application and removal such as wire cutters, mallet and power drill.
- The PREFIX does not withstand the stresses of weight-bearing.
- In order to counter-act the tightening forces, it is advisable to use the Prefix II Clamp Holder to grip the Rod Clip Clamp, the Screw Clip Clamp or the Transfixing Pin Clamp while locking them with the Universal T-Wrench or the 5mm Allen Wrench.

Correct assembly:

Incorrect assembly:

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

Implantable Device*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol "⊗" reported on the product label. After the removal from the patient, the implantable device* has to be discarded. The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

Non Implantable Device

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of a "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization release.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Aluminum based devices are damaged by alkaline ($\text{pH} > 7$) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- The contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing in order to remove soiling that accumulates in recesses.
- If a device needs particular care in pre-cleaning, a product specific IFU is available in Orthofix website, that is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on reusable fixators and instruments
- End of life is normally determined by wear and damage due to use
- Products labeled for Single Use only MUST NOT be reused regardless any reprocessing in a clinical setting

POINT OF USE

It is recommended to reprocess the reusable medical devices as soon as is reasonably practicable in order to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because this can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

It is recommended to cover contaminated instruments during transportation in order to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In case of highly contaminated reusable medical device before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described in the following paragraph) are recommended.

Manual Pre-cleaning

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution in order to displace trapped air.
5. Scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using a clean, lint free cloth.

CLEANING

General considerations

In these instructions Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to cleaning agents used. Staff shall follow the safety precautions to comply with the procedure of the healthcare facility using protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent's manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft bristled nylon brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in ultrasonic device with degassed cleaning solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends on the basis of the validation performed to use an ultrasound frequency of 35kHz, power 300Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the used and the concentration shall be in compliance with the detergent's manufacturer technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present it is possible to use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device and had to be removed with the brush, the cleaning step must be repeated as described above.
14. Carefully hand-dry using a clean, lint-free cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three time with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using a clean, lint-free cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a Washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices in order to locate the cannulations in vertical position and blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 min.
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of the detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 min at 55°C.
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 min.
 - d. Final rinsing with deionised water for 3 min.
 - e. Thermal disinfection at least 90°C or 194°F (max 95°C or 203°F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110°C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.
- The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent's manufacturer technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion on the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfecter when it complete the cycle.
13. If necessary, drain off excess water and dry by using clean, linf.free cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use. All functional checks and inspections described below also cover the interfaces with other instruments or components. The failure modes below may be caused by end of life of the product, improper use or improper maintenance. Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use. Specific instructions may be available for some product codes. These instructions are linked to the product code and are available in a dedicated Orthofix website. Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained to mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS).The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instruction from the manufacturer. When equipment and processes differ from those validated by Orthofix the health care facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always to use a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated but it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers but only in wraps.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations. These cleaning agents are not listed in preference to other available cleaning agents which may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0,5%

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

ADDITIONAL INFORMATION

INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation: the patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need of periodical medical follow-up and of the removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performances to his/her physician.

NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

Symbol	Description	
MD	Medical Device	
 	Consult instructions for use or consult electronic instructions for use	Caution: Consult instructions for use for important cautionary information
	Single Use. Do not re-use	Orthofix note: discard appropriately after the use (treatment) on patient
		Non-sterile
REF 	LOT 	Catalogue number
	Use-by date (year-month-day)	
	CE marking in conformity to applicable European Medical Device Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	Symbol for MR conditional. The item has been demonstrated to pose no known hazards in a specified MRI environment with specified conditions of use	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	
UDI	Unique Device Identifier	

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