Contours PHP

Proximal Humeral Plate™

EN	ORTHOFIX® CONTOURS PHP PROXIMAL HUMERAL PLATE SYSTEM 2	JA	ORTHOFIX® CONTOURS PHP 上腕骨近位部プレートシステム 74
IT	ORTHOFIX® CONTOURS PHP SISTEMA DI PLACCHE PER Omero PROSSIMALI 8	ZH	ORTHOFIX® CONTOURS PHP 肱骨近端接骨板系统 81
FR	ORTHOFIX® CONTOURS PHP SYSTÈME DE PLAQUE HUMÉRALE PROXIMALE 14	CS	ORTHOFIX® CONTOURS PHP SYSTÉM PROXIMÁLNÍCH HUMERÁLNÍCH DLAH 87
DE	ORTHOFIX® CONTOURS PHP PROXIMALES HUMERUSPLATTENSYSYSTEM 20	PL	ORTHOFIX® CONTOURS PHP SYSTEM PŁYTEK DO ZESPALANIA ZŁAMAŃ BLIŻSZEGO KOŃCA KOŚCI RAMIENNEJ 93
ES	ORTHOFIX® CONTOURS PHP SISTEMA DE PLACA HUMERAL PROXIMAL 26	SL	ORTHOFIX® CONTOURS PHP SISTEM PROKSIMALNE HUMERALNE PLOŠČICE 99
DA	ORTHOFIX® CONTOURS PHP PROKSIMALT HUMERUS PLADESYSYSTEM 32	KO	ORTHOFIX® CONTOURS PHP 근단 상완판 시스템 105
FI	ORTHOFIX® CONTOURS PHP PROXIMAL HUMERAL PLATE SYSTEM 38		
NO	ORTHOFIX® CONTOURS PHP PROKSIMALT HUMERUSPLATESYSYSTEM 44		
NL	ORTHOFIX® CONTOURS PHP PROXIMAAL HUMERUSPLAATSYSYSTEM 50		
PT	SISTEMA ORTHOFIX® CONTOURS PHP PROXIMAL HUMERAL PLATE 56		
SV	ORTHOFIX® CONTOURS PHP PROXIMALT HUMERUSPLATTYSYSTEM 62		
EL	ORTHOFIX® CONTOURS PHP ΣΥΣΤΗΜΑ ΕΓΓΥΣ ΒΡΑΧΙΟΝΙΑΣ ΠΛΑΚΑΣ 68		

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

0123 PQPHP C 02/22 (0424730)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online

Important information - please read prior to use

See also instruction leaflet PQRMD for reusable medical devices

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

Symbol	 Attention, see instructions for use PQ PHP
--------	--

DESCRIPTION

The Orthofix Contours PHP consists of a left and a right plate that has to be used in combination with at least two Diaphyseal Screws, a main locking screw and at least two titanium hexagonal Fine Threaded Screws. In addition titanium/polymeric cerclage wires can be used to achieve further stabilization of the bone fragments.

These implants are not intended to replace normal healthy bone. Load-bearing and activity levels will dictate the longevity of the implant. The patient should understand that stress on an implant can be produced without actual weight-bearing. In the absence of solid bony union, the weight of the limb alone, muscular forces associated with moving a limb, or repeated stresses of apparent relatively small magnitude, can result in failure of the implant. Therefore, the patient should follow the post-operative instructions given by the surgeon.

The product is intended for professional use only. Surgeons who supervise the use of the product must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the product. Prior to surgery, surgeons should be familiar with the devices, instruments and surgical procedure, including the application and removal. Detailed operative technique guidance is available on request; please contact Orthofix or your local distributor.

INDICATIONS

The Orthofix Contours PHP Proximal Humeral Plate is intended for fractures, osteotomies and non unions of the proximal humerus, particularly in osteopenic bone.

CONTRAINDICATIONS

The Orthofix Contours PHP is not designed or sold for any use except as indicated.

Use of the Orthofix Contours PHP is contraindicated in the following situations:

1. Active infection at the site of surgery
2. General medical conditions including: impaired blood supply, pulmonary insufficiency (i.e. ARDS, fat embolism) or latent infection.
3. Patients who are unwilling or incapable of following post-operative care instructions.
4. Obese patients.
5. Suspected or documented metal allergy or intolerance.

WARNINGS AND PRECAUTIONS

1. All equipment should be carefully examined prior to use. Product integrity and performance are assured only if the packaging is undamaged. DO NOT USE if packaging is compromised or if a component or instrument is believed to be faulty, damaged or suspect.
2. The Orthofix Contours PHP should not be used if the fracture is open or if the overlying skin is damaged.
3. Plates, screws, and hexagonal Fine Threaded Screws as well as certain parts of the instrumentation (where indicated on the label) are for single-use only and must never be reused.
4. Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been shown to affect bone healing and could potentially have an adverse effect on bone repair during fracture healing.
5. Care should be used in the handling and storage of the components. The implants should not be scratched, notched or otherwise damaged as these may reduce the functional strength of the component.
6. Fracture distraction for any time should be avoided during the operation, and fractures should never be locked in distraction, since this may cause delayed union.
7. Bending of the hexagonal Fine Threaded Screws during insertion must be avoided, since this may lead to breakage of the implant. When the chamfer is close to the plate, the speed of insertion must be reduced.
8. The selection of the main locking screw and hexagonal Fine Threaded Screws lengths are very important, they MUST NOT penetrate the articular surface.
9. All implants and wires used in the humeral head, such as the central screws, the hexagonal fine threaded screws and the Kirschner wires have to be positioned 5 - 10 mm from the second cortex.
10. If the thread of the hexagonal head of the Fine Threaded Screws is not in contact with the plate, the insertion of the hexagonal Fine Threaded Screws should be completed manually or the hexagonal Fine Threaded Screws removed.
11. The Fine Threaded Screws should have a maximum angulation of +/-15° with respect to the plane of the plate in order to ensure good purchase in their threaded hole.
12. It is recommended to apply the main locking screw, AT LEAST two hexagonal Fine Threaded Screws and AT LEAST two Diaphyseal Screws for stabilization of the plate.
13. All implants, Kirschner wires and cannulated drill bits are all single-use devices and should NEVER be reused.

14. The Orthofix Contours PHP is available in a RIGHT and LEFT version, which are anatomically shaped. It is therefore important to select the correct plate. DO NOT use a plate labelled "RIGHT" for a left humerus and vice versa.
15. As plates and screws are made of titanium, DO NOT use them in association with stainless steel implants, since they may cause an electrolytic reaction. Wires for cerclage of the humeral head may be made of titanium or of a polymeric material.
16. Do not use components of the Orthofix Contours PHP in conjunction with products of other manufacturers, unless otherwise specified.
17. Wherever a Kirschner-wire is used to guide a screw or instrument into position, the following precautions should be taken:
 - A) The wire should always be NEW and should NEVER be reused.
 - B) The wire should be checked before insertion to exclude any scratches or bends in the wire; if any are found it should be discarded.
 - C) During the introduction of any instrument or implant over a wire, the surgeon should screen the wire tip as continuously as possible to exclude inadvertently driving the wire further than intended.
 - D) During each pass of the instrument or implant, the surgeon should check that there is no bony or other debris built up on the wire or inside the instrument or implant which might cause it to bind on the wire and push it forward.
18. Careful monitoring of the progress of healing must be undertaken in all patients. If callus is slow to develop, other measures may be required to promote its formation, such as a bone graft, or exchanging the implant.
19. Any fixation device may break if subjected to the increased loading caused by delayed or non-union.
20. During treatment with the Orthofix Contours PHP the patient should not use crutches or other walking aids that could apply load on the plate.
21. It is very important that the operating surgeon has received instruction on the insertion technique for the Orthofix Contours PHP, and that the full set of instruments and implants are available.
22. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
23. Removal of the device: the surgeon should make the final decision whether an implant can be removed. Implant removal should be followed by adequate post-operative management to avoid refracture. If the patient is older with a low activity level, the surgeon may decide not to remove the implant, thus eliminating the risks involved in a second surgical procedure.
24. Additional equipment may be required for application and removal such as soft tissue retractors, etc.

Warning: MRI (Magnetic Resonance Imaging) image quality may be affected by the presence of the Orthofix Contours PHP. The Orthofix Contours PHP has not been evaluated for safety and compatibility in the MR environment. The Orthofix Contours PHP has not been tested for heating or migration in the MR environment.

POSSIBLE ADVERSE EFFECTS

1. Loosening, bending or breakage of implanted components.
2. Loss of anatomic positioning with malunion.
3. Scar formation possibly causing pain and/or neurological compromises around nerves.
4. Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke excessive bleeding, phlebitis, wound or bone necrosis, wound infection or damage to blood vessels or nerves.
5. Non-union or delayed union, which may lead to implant breakage.
6. Metal sensitivity, or allergic reaction to a foreign body.
7. Pain, discomfort, or abnormal sensations due to the presence of the device.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Important

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require surgical re-intervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the devices by the surgeon. Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE the Orthofix Contours PHP.

Materials

The Orthofix Contours PHP system is comprised of titanium alloy and stainless steel components. Implantable components are made of titanium alloy, while those components which come into contact with patients are made of surgical grade stainless steel.

STERILE and NON-STERILE Products

Orthofix provides titanium humeral plates, Main Locking- and Diaphyseal Screws, Kirschner wires and hexagonal Fine Threaded Screws STERILE.

The instrumentation is provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures.

RISKS DUE TO THE RE-USE OF "SINGLE-USE" DEVICE

Implantable Device*

The "SINGLE-USE" implantable device* of Orthofix is identified through symbol (⊗) reported on the product label. After the removal from the patient, the implantable device* has to be dismantled. The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device. Any device intended to be totally / partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days.

Non Implantable Device

The "SINGLE-USE" non implantable device of Orthofix is identified through symbol (⊗) reported on the label or are indicated in the "Instructions

For Use" supplied with the products. The re-use of "SINGLE-USE" non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.

3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
 - d. Final rinsing with deionized water for 3 minutes;
 - e. Thermal disinfection at least 90°C or 194° F (max 95 °C or 203° F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.

The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfector when it completes the cycle.
13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132° C (270° F)	132° C (270° F)	134° C (273° F)	134° C (273° F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description	
 	Consult instructions for use or consult electronic instructions for use	CAUTION: Consult instructions for use for important cautionary information
	Single Use. Do not re-use	Orthofix note: discard appropriately after use (treatment) on patient
	Sterile. Sterilized using irradiation	
	NON STERILE	
 	Catalogue number	Batch code
	Use-by date (year-month-day)	
	CE marking in conformity to applicable European Medical Device Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	

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Orthofix Srl