# LRS Advanced™

Limb Reconstruction System

|    |                                  |    |                                   |
|----|----------------------------------|----|-----------------------------------|
| EN | RECONSTRUCTION SYSTEMS<br>2      | ZH | 重建系统<br>111                       |
| IT | SISTEMI DI RICOSTRUZIONE<br>10   |    | REKONSTRUKČNÍ SYSTÉMY<br>119      |
| FR | SYSTÈMES DE RECONSTRUCTION<br>18 |    | SYSTEMY DO REKONSTRUKCJI<br>127   |
| DE | REKONSTRUKTIONSSYSTEME<br>26     |    | SISTEMI ZA REKONSTRUKCIJU<br>135  |
| ES | SISTEMAS DE RECONSTRUCCIÓN<br>35 |    | 재건 시스템<br>143                     |
| DA | REKONSTRUKTIONSSYSTEMER<br>44    |    | أنظمة إعادة البناء<br>152         |
| FI | KORJAUSJÄRJESTELMÄT<br>52        |    | СИСТЕМЫ ЗА РЕКОНСТРУКЦИЯ<br>160   |
| NO | REKONSTRUKSJONSSYSTEMER<br>60    |    | REKONSTRUKCIJOS SISTEMOS<br>169   |
| NL | RECONSTRUCTIESYSTEMEN<br>68      |    | SISTEM PEMBENTUKAN SEMULA<br>177  |
| PT | SISTEMAS DE RECONSTRUÇÃO<br>76   |    | SISTEME DE RECONSTRUCTIE<br>185   |
| SV | REKONSTRUKTIONSSYSTEM<br>84      |    | REKONSTRÜKSİYON SİSTEMLERİ<br>193 |
| EL | ΣΥΣΤΗΜΑΤΑ ΑΝΑΔΟΜΗΣΗΣ<br>92       |    | SUSTAVI ZA REKONSTRUKCIJU<br>201  |
| JA | 骨再建システム<br>101                   | LV | REKONSTRUKCIJAS SISTĒMAS<br>209   |
|    |                                  | TR |                                   |
|    |                                  | RO |                                   |
|    |                                  | MS |                                   |
|    |                                  | LT |                                   |
|    |                                  | BG |                                   |
|    |                                  | AR |                                   |
|    |                                  | KO |                                   |
|    |                                  | SL |                                   |
|    |                                  | CS |                                   |
|    |                                  | ZH |                                   |

ORTHOFIX Srl  
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy  
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PQ REC B 06/22 (0425684)

# LRS Advanced™

## Limb Reconstruction System

CE

**Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online**

**Important information - please read prior to use**

**See also instruction leaflet PQSCR (version D or higher) for implantable devices and related instruments and PQRMD for reusable medical devices**

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### MEDICAL DEVICE INFORMATION

#### DESCRIPTION

This leaflet refers to the following Reconstruction Systems: Adult Limb Reconstruction System (Adult LRS hereinafter), Paediatric Limb Reconstruction System (Paediatric LRS hereinafter), and Advanced Limb Reconstruction System (Advanced LRS hereinafter). The Reconstruction Systems are monolateral external fixators consisting of a series of components that build the external frame. The external frame is connected to the bone by means of bone screws. The Reconstruction Systems can be used in conjunction with Orthofix circular external fixators and Kirschner wires. Application and removal of the Reconstruction Systems can be performed with Orthofix general orthopedic instrumentation.

#### INTENDED PURPOSE AND INDICATIONS

##### INTENDED PURPOSE

The Reconstruction Systems are intended to provide bone fixation.

##### INDICATIONS FOR USE

The Reconstruction Systems are indicated for fractures, joint distraction, bone transport, lengthening and angular corrections in long bones.

#### CONTRAINDICATIONS

DO NOT USE the Reconstruction Systems if a surgical candidate exhibits or is predisposed to any of the following contraindications:

- mental or physiological conditions who are unwilling or incapable of following postoperative care instructions
  - inflammatory arthropathies and age over 45 years for arthrodiatasis of the hip
  - severe osteoporosis
  - HIV positive
  - severe, poorly controlled diabetes mellitus
  - suspected or documented metal sensitivity reactions
- as it could result in a treatment failure in the intended population.

#### INTENDED PATIENTS

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. The Reconstruction Systems are intended for adult and pediatric patients with the exception of newborns.

#### INTENDED USERS

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal).

#### NOTES FOR USE - IMPLANT REMOVAL

Once the treatment with external fixation is complete, implant must be removed. HCP should consider premature removal in case of adverse events.

#### DISCLAIMER

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

#### MATERIAL

The implants are made from an implant grade material that is specified on the product label.

## **WARNINGS**

1. All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
2. Compression is never recommended in a fresh fracture.
3. Axial displacement may occur if the body of the fixator is not in line with and parallel to the bone.
4. Medial or lateral translation may occur if the body of the fixator is not placed parallel to the diaphysis.
5. The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. If the fixator is sited at a distance of more than 4cm from the bone, the use of 3 screws per clamp is advisable.
6. Components may not be interchangeable between all Orthofix external fixation systems. Consult individual operative technique guides for interchangeable components.
7. This device is not approved for fixation or attachment by means of screws to posterior elements (pedicles) of the cervical, thoracic or lumbar spinal column.

### ADULT LRS ADV HINGE CONNECTION (SS53537) and PAEDIATRIC LRS HINGE CONNECTION (SS55537)

1. Position the Hinge Pin in the Center of rotation of the knee. Wrong positioning could limit the R.O.M. of the knee.

### ADV RADIOLUCENT LRS RAIL 500MM (SS53640R) and ADV RADIOLUCENT LRS RAIL 600MM (SS53660R)

1. Additional stabilization of the limb is required in all cases. At least two clamps should be used in each bone segment to guarantee stability during the entire bone healing process. The rail has a central area that precludes free sliding of the clamp. The clamp can nevertheless pass through this area by removing its clamp screw temporarily.

### ADV RING HINGE (53570)

1. This clamp has not been designed for lengthening procedures and should never be subjected to the forces of lengthening.
2. A minimum of three independent single screws are applied to the ring to guarantee rotational stability.

### ADV T-GARCHES CLAMP (53031)

1. It must be used always with 3 bone screws for sufficient stability.
2. Never position this clamp at a distance greater than 30mm from the bone.
3. During the lengthening procedure, the patient should be allowed to weightbear no more than 30% of the body weight.
4. This clamp should not be used for correcting procurvatum/recurvatum femoral deformities.

### KNEE HINGE (53590)

1. In case knee movement has to be blocked, do not tighten the posterior locking screw and use a compression-distraction unit or reinforcement bars to prevent knee movement.

## **PRECAUTIONS**


1. The surgeon must evaluate the integrity of the construct at follow-up visits.
2. The fracture or bone gap must be checked periodically during treatment, making any necessary adjustments to the fixation. An excessive or persistent gap can delay consolidation.
3. In patients undergoing callus distraction, the regenerated bone must be checked regularly and monitored radiologically.
4. Insert the screw in the Advanced Inclination Clamp only when the bone-fixator distance is correct (once the screws are inserted, the fixator distance from the bone can no longer be adjusted due to the offset of the screws).
5. When performing angular deformity correction with the ADV Micrometric Translation-Angulation Clamp, the weight of the limb should be supported to prevent excessive strain on the bone screws.
6. During and after insertion of the implants, ensure their correct positioning under image intensification.

### ADULT LRS ADV HINGE CONNECTION (SS53537) and PAEDIATRIC LRS HINGE CONNECTION (SS55537)

1. The Hinge Pin must always be inserted into one of the Dyna Block holes closest to the distal diaphyseal bone screws.
2. Once all bone screws have been inserted, check that flexion-extension of the knee is not impeded. This will confirm correct positioning of the hinge. If necessary, replace the K-wire and Hinge Pin.
3. The SS53537 and SS55537 have not been tested in traumatic and severely unstable knee conditions. These could require an increased stiffness to ensure congruent healing and a more anatomical rotation to avoid worsening of joint conditions.
4. An additional stabilizing rod should be applied anteriorly, between the anterior face of the Dyna Block and the Ring, to span the hinge and protect from a knee flexion contracture when knee movement is not controlled (e.g. patient is asleep). (Use only the central threaded hole present in the Dyna Block to fix the stabilizing rod).
5. Lengthening of the femur along mechanical axis could generate anatomical axis mal-alignment. Adequate preoperative planning of corrective procedures should be performed to avoid this condition.

### KNEE HINGE (53590)

1. The Knee Hinge must be correctly aligned with the Reference Axis of the knee joint. In order to achieve this, a K-wire should be inserted into the Reference Axis of the femur (see illustration below) and the Knee Hinge applied over it through the 2mm centering hole positioned on its proximal part.

2. Before stabilizing the femur with bone screws, check under image intensification that flexion-extension of the knee is not impeded. This will confirm correct positioning of the hinge. If necessary, replace the K-wire correctly in the Reference Axis.
3. Once all bone screws have been inserted, check again the R.O.M. and if not complete, tighten the posterior locking screw to limit the hinge movement to the R.O.M.

All Orthofix devices should be used together with their corresponding Orthofix implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the Reconstruction Systems when used in conjunction with devices of other manufacturers or with other Orthofix devices if not specifically indicated in the Operative Technique.

#### **POSSIBLE ADVERSE EVENTS**

- Non-union, delayed union or malunion
- Superficial infection
- Deep infection
- Loss of fixation
- Bending, breakage or migration of the device
- Reoperation to replace a component or entire frame configuration
- Bone fracture during or after treatment
- Damage to bone or surrounding tissues resulting from surgery or treatment
- Tension affecting the soft tissues and/or the frame during callus manipulation (e.g. corrections of bony deformity and/or bone lengthening)
- Joint contracture, dislocation, instability or loss of range of motion
- Proximal migration of the fibula, due to non-anchorage to the adjacent tibia
- Pain, discomfort or abnormal sensations due to the presence of the device
- Residual deformities, persistence or recurrence of the initial condition requiring treatment
- Premature bone callus consolidation during distraction
- Wound healing complications
- Stiffness at the surgery site
- Compartment syndrome
- Events caused by intrinsic risks associated with anesthesia and surgery

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP.

#### **MRI (Magnetic Resonance Imaging) SAFETY INFORMATION**

The Orthofix Reconstruction Systems have not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Orthofix Reconstruction Systems in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### **EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE**

##### **EXPECTED CLINICAL BENEFITS**

- Stability and safety in corrections
- Easy to apply. Versatile fixator, suitable for the treatment of various fracture patterns
- Reduces surgery time
- For the patient: comfort and ease with a monolateral fixator
- Small sizes specifically developed for patient comfort
- Accurate deformity correction using a monolateral fixator
- Facilitates ease of application
- For the surgeon: short learning curve
- Stable fixation during angular correction. Versatility of use in angular corrections
- No need for a second surgery for removal (compared to internal fixation)

##### **PERFORMANCE CHARACTERISTICS OF THE DEVICE**

- Possibility of positioning screws in different planes matching various bone curvatures
  - Availability of radiolucent components
  - Compression distraction unit with simple turning mechanism
- Advanced Dyna-Ring (53536)
- To be locked to the rail with its silicone cushion just in contact with the clamp which has been unlocked to enable dynamization

#### **SPECIFIC INFORMATION ON THE DEVICE**

- Additional equipment may be required for fixation application and removal such as wire cutters, mallet and power drill
- Dynamization and physical therapy guidelines should be followed based on each individual case and the fixation system used, and should be instituted as and when considered appropriate by the surgeon, in accordance with clinical and radiological findings
- In case Advanced LRS is used in conjunction with TrueLok or TL-HEX ring refer to instructions available in leaflet PQTLK and PQWTN

#### **ADULT LRS ADV HINGE CONNECTION (SS53537) and PAEDIATRIC LRS HINGE CONNECTION (SS55537)**

- The Hinge Connections are available in two models: adult (SS53537) and pediatric (SS55537). SS53537 is used in association with a LRS ADV rail, while the SS55537 with a Pediatric LRS rail. Hinge Connections are used either with a Sheffield or a TrueLok Ring.

#### **ADV RING HINGE (53570)**

- This clamp may be used in gradual or acute angular corrections as long as a minimum of 3 screws are applied to the ring to guarantee rotational stability, with two screws inserted above and 1 below the ring.

- The two outer screws should be inserted at an angle close to but no greater than 90°. This means inserting each screw into hole number 6 in a 150mm ring, or 5 in all greater ring sizes, on each side of the hole number 0 above the centre of rotation of the hinge.
- It can also be used with a minimum of 3 tensioned Kirschner wires with a crossing angle of at least 60 degrees, but only in cases of acute correction, as long as reinforcement bars are added to the frame after the correction.

#### ADV T-GARCHE CLAMP (53031)

- This clamp must be positioned at a distance of 20mm from the bone (never greater than 30mm). Use of the spacer is necessary for correct placement of the clamp.

#### KNEE HINGE (53590)

- This clamp is used in association with a LRS ADV rail and either a Sheffield or TrueLok ring. It permits a similar range of knee motion to that of the Anterior and Posterior Cruciate Ligaments, from 0° to 90° in the sagittal plane. This range of motion can be limited and/or blocked when appropriate.

### RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

#### **Implantable Device\***

The "SINGLE USE" implantable device\* of Orthofix is identified through symbol "⊗" reported on the product label. After the removal from the patient, the implantable device\* has to be discarded. The re-use of implantable device\* introduces contamination risks for users and patients. The re-use of implantable device\* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\* ) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

#### **Non Implantable Device**

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of a "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

### STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

### INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions provided.

#### **Warnings**

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization release.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Aluminum based devices are damaged by alkaline (pH>7) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- The contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing in order to remove soiling that accumulates in recesses.
- If a device needs particular care in pre-cleaning, a product specific IFU is available in Orthofix website, that is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

#### **Limitations on reprocessing**

- Repeated reprocessing has minimal effect on reusable fixators and instruments.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless any reprocessing in a clinical setting.

#### **POINT OF USE**

It is recommended to reprocess the reusable medical devices as soon as is reasonably practicable in order to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because this can cause the fixation of residue.

#### **CONTAINMENT AND TRANSPORTATION**

It is recommended to cover contaminated instruments during transportation in order to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

## **PREPARATION FOR CLEANING**

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In case of highly contaminated reusable medical device before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described in the following paragraph) are recommended.

### **Manual Pre-cleaning**

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution in order to displace trapped air.
5. Scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using a clean, lint free cloth.

## **CLEANING**

### **General considerations**

In these instructions Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices andto cleaning agents used.

Staff shall follow the safety precautions to comply with the procedure of the healthcare facility using protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent's manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

### **Manual cleaning**

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft bristled nylon brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in ultrasonic device with degassed cleaning solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends on the basis of the validation performed to use an ultrasound frequency of 35kHz, power 300Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the used and the concentration shall be in compliance with the detergent's manufacturer technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present it is possible to use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device and had to be removed with the brush, the cleaning step must be repeated as described above.
14. Carefully hand-dry using a clean, lint-free cloth.

### **Manual disinfection**

1. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three time with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using a clean, lint-free cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

## **Automatic cleaning and disinfection using washer-disinfector**

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
    - a. Cannulations
    - b. Long blind holes
    - c. Mating surfaces
    - d. Threaded components
    - e. Rough surfaces
  2. Use a Washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
  3. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
  4. Ensure that the washer-disinfector and all services are operational.
  5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
  6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
  7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
  8. Arrange medical devices in order to locate the cannulations in vertical position and blind holes incline downwards to promote the leakage of any material.
  9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
    - a. Pre-cleaning for 4 min.
    - b. Cleaning with the appropriate solution. Orthofix recommends the use of the detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 min at 55°C.
    - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 min.
    - d. Final rinsing with deionised water for 3 min.
    - e. Thermal disinfection at least 90°C or 194°F (max 95°C or 203°F) for 5 minutes or until A0=30000 is reached. The water used for thermal disinfection must be purified.
    - f. Drying at 110°C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.
- The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent's manufacturer technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
  11. On completion on the cycle, ensure that all stages and parameters have been achieved.
  12. Wearing protective equipment unload the washer disinfecter when it complete the cycle.
  13. If necessary, drain off excess water and dry by using clean, linf.free cloth.
  14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

## **MAINTENANCE, INSPECTION AND FUNCTION TESTING**

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use. All functional checks and inspections described below also cover the interfaces with other instruments or components. The failure modes below may be caused by end of life of the product, improper use or improper maintenance. Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use. Specific instructions may be available for some product codes. These instructions are linked to the product code and are available in a dedicated Orthofix website. Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

Cams and bushes of Reconstruction Systems are for SINGLE USE ONLY. They MUST be discarded and replaced every time a fixator is cleaned after use and prior to sterilization.

## **PACKAGING**

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained to mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instruction from the manufacturer. When equipment and processes differ from those validated by Orthofix the health care facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

## STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

| Steam sterilizer type        | Gravity           | Pre-vacuum    | Pre-vacuum        | Pre-vacuum     |
|------------------------------|-------------------|---------------|-------------------|----------------|
| Notes                        | Not for use in EU | -             | Not for use in US | WHO guidelines |
| Minimum exposure Temperature | 132°C (270°F)     | 132°C (270°F) | 134°C (273°F)     | 134°C (273°F)  |
| Minimum exposure Time        | 15 minutes        | 4 minutes     | 3 minutes         | 18 minutes     |
| Drying Time                  | 30 minutes        | 30 minutes    | 30 minutes        | 30 minutes     |
| Number of pulses             | N/A               | 4             | 4                 | 4              |

Orthofix recommends always to use a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated but it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers but only in wraps.

## CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations. These cleaning agents are not listed in preference to other available cleaning agents which may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0,5%

## STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

## DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

## ADDITIONAL INFORMATION

### INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation: the patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need of periodical medical follow-up and of the removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performances to his/her physician.

The HCP shall instruct the patient to:

- Make the adjustments or get help in making the adjustments as needed
- Report if adjustment schedule cannot be met
- Report any adverse or unexpected effects

### NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

### MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

| Symbol  | Description  |  |
|---|--|--|
| <b>MD</b>   | Medical Device   |  |
|   | Consult instructions for use or consult electronic instructions for use                    | Caution: Consult instructions for use for important cautionary information |
|    | Single Use. Do not re-use  | Orthofix note: discard appropriately after the use (treatment) on patient  |
|    | Non-sterile  |  |
| <b>REF</b>  | <b>LOT</b>   | Catalogue number      Batch code   |
|    | Use-by date (year-month-day)   |  |
|    | CE marking in conformity to applicable European Medical Device Directives/Regulations      |  |
|   | Date of manufacture  | Manufacturer   |
|    | Do not use if package is damaged and consult instructions for use                          |  |
| <b>Rx Only</b>  | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician |  |

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