VeroNail™

System

EN	ORTHOFIX® TITANIUM TROCHANTERIC NAILING SYSTEM 2	JA	ORTHOFIX® チタニウム・転子用・ネイリング・システム 81
IT	SISTEMA DI CHIODI TROCANTERICI IN TITANIO ORTHOFIX® 8	ZH	ORTHOFIX® 钛质转子骨钉系 88
FR	TROCHANTERIQUE EN TITANE ORTHOFIX® RX UNIQUEMENT 15	CS	TITANOVÁ SOUPRAVA TROCHANTERICKÝCH FIXAČNÍCH HŘEBŮ ORTHOFIX® 94
DE	ORTHOFIX® TITAN-TROCHANter-NAGELUNGSSYSTEM 22	PL	SYSTEM KREȚARZOWYCH GWOŹDZI TYTANOWYCH ORTHOFIX® 100
ES	SISTEMA DE ENCLAVAMIENTO TROCANTÉRICO ORTHOFIX® TITANIUM 29	SL	SISTEM TITANOVIH ŽEBLJEV ORTHOFIX® 107
DA	ORTHOFIX® TROCHANTÆRE MARVSØMSYSTEM AF TITANIUM 36	KO	ORTHOFIX® 티타늄 전자부 골수정 기구 113
FI	REISILUUN PÄÄN ORTHOFIX®-TITAANIKIINNITYSJÄRJESTELMÄ 42	AR	جهاز التسمير المدورى من التيتانيوم من ORTHOFIX® 120
NO	ORTHOFIX® TROKANTÄRT TITANSTIFTESYSTEM 48	BG	ORTHOFIX® ТИТАНОВА СИСТЕМА ЗА ЗАКОВАВАНЕ ЗА ТРОХАНТЕР 128
NL	ORTHOFIX® TITANIUM TROCHANTER-NAGELSYSTEEM 54	LT	„ORTHOFIX®“ TITANO GÜBRIÖ VINIŲ SISTEMA 135
PT	SISTEMA ORTHOFIX® DE HASTES TROCANTÉRICAS DE TITÂNIO 61	MS	ORTHOFIX® SISTEM PEMAKUAN TROKANTERIK TITANIUM 141
SV	ORTHOFIX® TROKANTÄRT TITANPIKSYSTEM 68	RO	SISTEMUL DE OSTEOSINTEZĂ TROHANTERIANĂ CU TIJĂ DE TITAN ORTHOFIX® 148
EL	TPOXANTHPA TITANIOY ORTHOFIX® 74	TR	ORTHOFIX® TİTANYUM TROKANTERİK ÇİVİLEME SİSTEMİ 154

 **ORTHOFIX®**

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

PQTTN H 02/22 (0424725)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online

Important information - please read prior to use

See also instruction leaflet PQRMD for reusable medical devices

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

INDICATIONS

The Orthofix Titanium Trochanteric Nailing System is a set of intramedullary nails used for treatment of intertrochanteric and subtrochanteric fractures of the femur alone or when these fractures occur in combination with diaphyseal fractures. The use of the Titanium Trochanteric Nailing System is indicated for the alignment, stabilization and fixation of various types of fractures or deformities caused by trauma or disease (provided that the epiphyses are closed). These include: traumatic fractures, re-fractures, non-union, reconstruction, malunion, malalignment, pathological fractures and impending pathological fractures.

DESCRIPTION

The Orthofix Titanium Trochanteric Nailing System consists of Intramedullary Nailing Implants for the femur with respective end-caps, locking screws. These implants are not intended to replace normal healthy bone. Loads produced by weight bearing and activity levels will dictate the longevity of the implant. The patient should understand that stress on an implant can be produced without actual weight bearing. In the absence of solid bony union, the weight of the limb alone, muscular forces associated with moving a limb, or repeated stresses of apparent relatively small magnitude, can result in failure of the implant. Therefore, the patient should follow the post-operative instructions given by the surgeon.

The product is intended for professional use only. Surgeons who supervise the use of the product must have full awareness of orthopedic fixation procedures and should have received adequate training in the use of the product. Prior to surgery, surgeons should be familiar with the devices, instruments and surgical procedure, including the application and removal. Detailed operative technique guidance is available on request; please contact Orthofix or your local distributor.

CONTRAINDICATIONS

The Orthofix Titanium Trochanteric Nailing System is not designed or sold for any use except as indicated. Use of the Intramedullary Nailing Implants is contraindicated in the following situations:

1. Active or latent infection in the affected area.
2. General medical conditions including: impaired blood supply, pulmonary insufficiency (i.e. ARDS, fat embolism).
3. Patients who are unwilling or incapable of following post-operative care instructions.
4. Obese patients.
5. Suspected or documented metal allergy or intolerance.

GENERAL WARNINGS

- 1) Ensure that all components needed for the operation are available in the operating theatre.
- 2) Implants, nail end caps, locking screws as well as certain parts of the instrumentation (where indicated on the label), are for single use only and must never be reused. If any implant has come into contact with any body fluid it should be considered to have been used.
- 3) Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been shown to affect bone healing and could potentially have an adverse effect on bone repair during fracture healing.
- 4) This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

GENERAL PRECAUTIONS

The nail has been developed for use in pectrochanteric, intertrochanteric and subtrochanteric fractures of the femur alone or in combination with diaphyseal fractures. There are two versions: **Short**, 200 mm long, and **Long**, in left and right version, 280-440 mm long, in 20 mm increments. Both versions are green and are cannulated, and have a diameter of 15 mm proximally and 10 mm distally. Proximally they are locked with specific cephalic screws, distally with 4.8 mm standard locking screws. It might be physically possible to insert other types of screws, but Orthofix strongly advises against their use.

- 1) Cephalic screws are to be inserted into the neck of the femur and are available as sliding screws for parallel insertion and convergent screws for convergent insertion. Cephalic sliding screws are yellow, have an integrated outer tube and should be used in type 31.A1 and 31.A2 fractures. When treating this fracture type in combination with a diaphyseal fracture of the femur, a long nail should be used. Green cephalic convergent screws are for type 31.A3 fractures. When treating this fracture type in combination with a diaphyseal fracture of the femur, a long nail should be used. In case of applying a long nail in 32A, 32B and 32C type fracture, without per- or inter-trochanteric fractures (31.A1 and 31.A2), convergent cephalic screws should be used. Each should be used for the appropriate configuration and they should NEVER be used together.
- 2) Both the cephalic screws must be used for each version of the nail (short or long) and for each configuration (sliding or convergent).
- 3) When measuring the length of the cephalic screws, the length of both the proximal and distal locking screw is read on the scale of the ruler and note that the proximal sliding cephalic screw should be 10 mm shorter than the distal sliding cephalic screw. When using convergent cephalic screws, the proximal screw must be 10 mm shorter than the distal screw.
- 4) Cephalic screws are inserted using the long screw guide.
- 5) Distal locking in the long nail is performed using the free hand technique.
- 6) Implants are made of titanium alloy. Do not use with implants of dissimilar metals, since this may cause an electrolytic reaction.

- 7) Do not use components of the Orthofix nailing system in conjunction with products of other manufacturers, unless otherwise specified.
- 8) Care should be used in the handling and storage of the components. The implants should not be scratched, notched or otherwise damaged as these may reduce the functional strength of the component. Contours or bending of an implant should be avoided.
- 9) All components should be carefully examined prior to use. Product integrity, sterility (in the case of sterile products) and performance are assured only if the packaging is undamaged. DO NOT USE if packaging is compromised or if a component is believed to be faulty, damaged or suspect.
- 10) It is recommended to assemble the system prior to implantation.
- 11) Correct selection of the implant model and size is extremely important. An unstable or severely comminuted fracture, or a fracture in the metaphyseal area, may need a longer nail for adequate stabilization. The appropriate type and size should be selected for the patient. Improper positioning of nail and screws may result in loosening, bending, cracking, or fracture of the device or bone or both.
- 12) Care should be taken to attach all components correctly, ensuring that any location marks correspond, and locking them securely. It is important to check that the required joint between instrument remain firmly fixed during the procedure.
- 13) The length of locking screws should be determined by using a proper K-wire in a corresponding K-wire Guide, or a drill bit and a drill stop in a corresponding drill guide, and measuring, with the aid of the correct scale provided with the system, the amount of K-wire protruding from the K-wire guide or the drill bit protruding beyond the drill stop from the drill guide.
- 14) When locking screws are inserted, care should be taken not to penetrate the articular surface.
- 15) It is important that the locking screws are not inserted close to the fracture line, as effective screw fixation may be compromised, causing implant or fixation failure. In very distal fractures, it is important to ensure that the more proximal of the distal locking screws is at least 1 cm distal to the fracture line. Exceptions are the sliding parallel screws that are designed to cross a fracture line.
- 16) Care should be taken to achieve and maintain reduction before nail insertion, and to choose the correct locking screw length.
- 17) Fracture distraction for any time should be avoided during the operation, and fractures should never be locked in distraction, since this may cause delayed union and/or compartment syndrome.
- 18) Nails may be inserted with or without reaming, depending on patient, fracture type and bone diameter and quality. This decision is up to the operating surgeon.
- 19) Hammering of the nails may only be performed through the Insertion Knob. Do not hammer any of the aluminium instruments, as they will be damaged. Hammering should always be gentle; vigorous hammering should never be necessary. The surgeon should never persist with hammering if the implant is not advancing, but should review the situation, and consider further reaming.
- 20) Cannulated drill bits should be inspected prior to use to confirm that the lumen is free from obstruction. The correctly sized K-wire should be passed through it, to check that it slides easily, before and after each use.
- 21) Prior to any use, wires that are used to guide cannulated instruments or implants should be checked and should be discarded if found to be scratched or bent.
- 22) During the introduction of any instrument or implant over a wire, the wire tip should be screened, using fluoroscopy, as continuously as possible to exclude inadvertently driving the wire further than intended. This is particularly important when the wires are pointing towards potentially dangerous locations, such as the pelvis. The surgeon should check that there is no bony or other debris built up on the wire or inside the instrument or implant which might cause it to bind on the wire and push it forward.
- 23) After the procedure, check that fracture reduction and the position of all implants are correct using the image intensifier.
- 24) General advice on weightbearing (unless otherwise specified): full weightbearing can be commenced when there is radiological evidence of bridging callus. In stable fractures (e.g. types 31A1, AO classification), where there is good contact between the two intact segments of bone, so that load sharing can be anticipated, early weightbearing as tolerated should be encouraged. In unstable fractures (e.g. types 31A2 and 31A3, AO classification), weightbearing is allowed when good bony contact has been achieved. When the bone is comminuted or there is bone loss, so that load sharing will not be possible until callus is formed, initial weightbearing should be only partial. The exact amount of load to be carried depends on the size of implant inserted and the stature of the patient. Hip and knee mobility should always be encouraged within the pain limits. However, the best clinical results are obtained by encouraging early mobility and full weightbearing as tolerated as early as possible and according to the patient's condition. A patient may have a diaphyseal fracture in combination with any of the above types of fracture. There will be a wide variation in fracture stability and the potential for early load sharing in these situations, and each individual patient must be advised as appropriate. Generally, early weightbearing as tolerated should be encouraged. If the fracture is very comminuted and unstable, weightbearing should be limited and protected by the use of crutches or a walking frame until early signs of healing are visible on the X-rays.
- 25) Careful monitoring of the progress of healing must be undertaken in all patients. If callus is slow to develop, other measures may be required to promote its formation, such as dynamisation of the implant, a bone graft, or replacing the implant.
- 26) Additional equipment may be required for fixation application and removal such as soft tissue retractors, flexible reaming set, cannulated drills, etc.
- 27) Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
- 28) The surgeon should make the final decision whether an implant can be removed. Generally intramedullary devices are removed if necessary about one year after bone healing. Implant removal should be followed by adequate post-operative management to avoid refracture.

POSSIBLE ADVERSE EFFECTS

1. Loosening, bending or breakage of implanted components.
2. Loss of anatomic positioning with malunion.
3. Scar formation possibly causing pain and/or neurological compromises around nerves.
4. Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke excessive bleeding, phlebitis, wound or bone necrosis, wound infection or damage to blood vessels or nerves.
5. Non-union or delayed union, which may lead to implant breakage.
6. Metal sensitivity, or allergic reaction to a foreign body.
7. Pain, discomfort, or abnormal sensations due to the presence of the device.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the devices by the surgeon.

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE the Orthofix Titanium Trochanteric Nailing System.

MRI SAFETY INFORMATION

The Orthofix Veronail Titanium Trochanteric Nail has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Orthofix Veronail Titanium Trochanteric Nail in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

MATERIALS

The Orthofix Titanium Trochanteric Nailing System is comprised of implantable titanium parts while the instrumentation is made of stainless steel or aluminum alloy. Those components of the instrumentation that come into contact with patients are made of surgical grade stainless steel.

STERILITY

All implantable components of the Orthofix Titanium Trochanteric Nailing System are supplied STERILE. The instrumentation is supplied NON-STERILE and requires cleaning and sterilization prior to use, following the respective recommended procedures. Please review the product label to determine the sterility of each device.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
 - d. Final rinsing with deionized water for 3 minutes;
 - e. Thermal disinfection at least 90°C or 194°F (max 95 °C or 203°F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.

f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.

The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.

10. Select and start a cycle according to the recommendations of the washer manufacturer.

11. On completion of the cycle, ensure that all stages and parameters have been achieved.

12. Wearing protective equipment unload the washer disinfecter when it completes the cycle.

13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.

14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym
concentration 2%
- For manual cleaning: Neodisher Mediclean
concentration 2%
- For automated cleaning: Neodisher Mediclean
concentration 0.5%

INFORMATION

For further information please contact your local representative or call 800.266.3349 (within US) or +39 0456719000 (Customer Service International, outside US).

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description	
 	Consult instructions for Use	CAUTION: Consult instructions for use for important cautionary information
	Single use. Do not reuse	
 R	STERILE. Sterilised by irradiation	
	NON STERILE	
 	Catalogue number	Batch code
	Expiry date (year-month-day)	
 	CE marking in conformity to applicable European Directives/Regulations	
 	Date of manufacture	Manufacturer
	Do not use if package is opened or damaged	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	

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Orthofix Srl