EN	GALAXY UNYCO™	ZH	GALAXY UNYCO™
IT	GALAXY UNYCO™	CS	GALAXY UNYCO™
FR	GALAXY UNYCO™	PL	GALAXY UNYCO™
DE	GALAXY UNYCO™	SL	GALAXY UNYCO™
ES	GALAXY UNYCO™	KO	GALAXY UNYCO™
DA	GALAXY UNYCO™	AR	GALAXY UNYCO™ جهاز التثبيت
FI	GALAXY UNYCO™	BG	GALAXY UNYCO™
NO	GALAXY – UNYCO™	LT	GALAXY UNYCO™
NL	GALAXY UNYCO™	MS	GALAXY UNYCO™
PT	GALAXY UNYCO™	RO	GALAXY UNYCO™
SV	GALAXY UNYCO™	TR	GALAXY UNYCO™
EL	GALAXY UNYCO™	HR	GALAXY UNYCO™
JA	GALAXY UNYCO™	LV	GALAXY UNYCO™

EN

Instructions For Use (IFU) are subject to change; the most current version of each IFU is always available online

Important information - please read prior to use

See also instruction leaflet PQSCR (version D or higher) for implantable devices and related instruments and PQRMD for reusable medical devices

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MEDICAL DEVICE INFORMATION

DESCRIPTION

The GALAXY UNYCO™ Fixator (GALAXY UNYCO hereinafter) is a monolateral external fixator consisting of a series of components that build the external frame. The external frame is connected to the bone by means of monocortical bone screws. The GALAXY UNYCO is provided to the user in sterile sets including clamps, rods, screws and specific application and removal tools. The GALAXY UNYCO is compatible with the GALAXY FIXATION™ System (GALAXY FIXATION hereinafter).

INTENDED PURPOSE AND INDICATIONS

INTENDED PURPOSE

The GALAXY UNYCO is intended to provide bone fixation.

INDICATIONS FOR USE

The GALAXY UNYCO is indicated for temporary fixation of tibia and ankle fractures, soft tissue defects and dislocations.

CONTRAINDICATIONS

DO NOT USE the GALAXY UNYCO if a surgical candidate exhibits or is predisposed to any of the following contraindications:

- Mental or physiological conditions who are unwilling or incapable of following postoperative care instructions
- Severe osteoporosis
- Malignancy in the fracture area
- Suspected or documented metal sensitivity reactions

as it could result in a treatment failure in the intended population.

INTENDED PATIENTS

Proper patient selection and the patient's ability to comply with physician instructions and follow the prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. The GALAXY UNYCO is intended for adult patients.

INTENDED USERS

The product is intended for use by Healthcare Professionals (HCP) only and such HCP must have full awareness of the appropriate orthopedic procedures and must be familiar with the devices, instruments and surgical procedures (including application and removal).

NOTES FOR USE - IMPLANT REMOVAL

Once the temporary treatment with external fixation is complete, implant must be removed and replaced with definitive treatment. HCP should consider premature removal in case of adverse events. If the GALAXY UNYCO is perceived as impediment for application of the definitive treatment, remove the system parts where needed.

DISCLAIMER

The HCP is fully responsible for the selection of the appropriate treatment and of the relevant device for the patient (including post-operative care).

MATERIAL

The implants are made from an implant grade material that is specified on the product label.

WARNINGS

- The GALAXY UNYCO does not withstand the stresses of weight-bearing. Excessive loading may lead to loosening or disassembly.
- If faced with a segmental fracture, the intermediate segment can be held using additional uni-cortical screws in different planes in GALAXY FIXATION Large Clamps attached to the connecting rod.

- This device is not approved for fixation or attachment by means of screws to posterior elements (pedicles) of the cervical, thoracic or lumbar spinal column.
- All equipment in the GALAXY UNYCO should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.

PRECAUTIONS

- Fracture stabilization has to be done following correct fracture reduction.
- The arms of the clamps can be disassembled to allow the most suitable configuration while the other components cannot be disassembled.
- Before fixator application, make sure that the clamps are fully loosened.
- Once the first screw has been inserted freehand and its correct position on the bone has been checked, apply the clamp on the first screw and tighten the metal ring on the arm clockwise.
- A cluster of uni-cortical screws should never sit across the fracture line. Consider the four screws inserted in the clamps as a cluster as well as the three screws connected to the Radiolucent Foot Unit.
- The system stability is guaranteed only with 4 UNYCO screws coupled with the Large Multiscrew Clamp in each segment. The Large Multiscrew Clamp for UNYCO Screws should never be applied across the fracture line.
- Frame stability must be checked intra-operatively before the patient leaves the operating theatre.
- If one UNYCO screw is perceived as unstable or if its repositioning is necessary, re-insert it in a new bone area avoiding any overlap with existing hole and using the 10° freedom of movement to find the new insertion point.
- Remember that the first screw determines the position of the remaining screws in the same bone segment, therefore safe corridors have to be considered in advance for correct and safe screw insertion.
- Once converging uni-cortical screws have been inserted, the distance of the clamp from the soft tissue surface cannot be altered. It is therefore important to determine the final distance of the fixator from the skin before inserting the second screw in the clamp.
- Stability of the clamp relies on correct spacing between the screws at the bone surface. To guarantee the system stability, screws must be accurately inserted aiming at the center of the bone.
- The Radiolucent Foot Unit must be used with three uni-cortical screws.
- During and after insertion, ensure correct positioning of the implants under image intensification. Check screw insertion both in the AP and lateral planes.
- Compression is never recommended in a fresh fracture.
- Once two or more UNYCO screws have been inserted in the Large Multiscrew Clamp for UNYCO Screws, the latter should not be moved along the screw shafts. Sliding the clamp might cause loss of bone purchase.
- The UNYCO screw may protrude beyond the first cortex or sit close to the insertion point of a nail. In this case pay attention when reaming or inserting a nail at a later stage.
- Ensure there is a sufficient length of Rod on either side of both proximal and distal clamps so as to enable reduction maneuvers which may require distraction (lengthening) between the two clamps.
- The surgeon must evaluate the integrity of the construct at follow-up visits.

All Orthofix devices should be used together with their corresponding Orthofix implants, components, accessories and instrumentation following the Operative Technique recommended by the manufacturer. Orthofix does not guarantee the safety and effectiveness of the GALAXY UNYCO when used in conjunction with devices of other manufacturers or with other Orthofix devices if not specifically indicated in the Operative Technique.

POSSIBLE ADVERSE EVENTS

- Damage to surrounding tissues due to surgical trauma
- Superficial infection
- Deep infection
- Compartment syndrome
- Bone fracture during or after treatment
- Bending, breakage or migration of the device
- Loss of fixation
- Pain, discomfort or abnormal sensations due to the presence of the device
- Wound healing complications
- Complex Regional Pain Syndrome
- Events caused by intrinsic risks associated with anesthesia and surgery

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the device are important considerations in the successful utilization of the device by the HCP.

MRI (Magnetic Resonance Imaging) SAFETY INFORMATION

GALAXY UNYCO is labeled MR CONDITIONAL  according to the terminology specified in ASTM F2503 "Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment".

All components of GALAXY UNYCO must be identified as MR Conditional prior to being placed in or near an MR environment.

Non-clinical testing has demonstrated that the GALAXY UNYCO components are MR Conditional. GALAXY UNYCO can be safely scanned under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla.
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90 mT/cm).
- Maximum whole-body-averaged Specific Absorption Rate (SAR) of ≤ 4.0 W/kg (First Level Controlled Operating Mode).
- No local transmit/receive coils must be used on the device.
- The GALAXY UNYCO components must be entirely outside the MR scanner bore.
- No part of the GALAXY UNYCO components must extend into the MR bore. Therefore MR scanning of body parts where the GALAXY UNYCO components are located is contraindicated.

Displacement Information

The GALAXY UNYCO will not present an additional risk or hazard to a patient in the 1.5 and 3.0 Tesla MR environment with regard to translational attraction or migration and torque.

Heating Information

Under the scan conditions defined above, the GALAXY UNYCO is expected to produce a maximum temperature rise of 2°C after 15 minutes of continuous scanning.

PATIENT SAFETY

MRI in patients with GALAXY UNYCO components can only be performed under these parameters. It is not allowed to scan the GALAXY UNYCO components directly. Using other parameters, MRI could result in serious injury to the patient. When the GALAXY UNYCO components are used in conjunction with the GALAXY FIXATION, please be advised that this combination has not been tested in the MR environment and therefore higher heating and serious injury to the patient may occur. Because higher in vivo heating cannot be excluded, close patient monitoring and communication with the patient during the scan is required. Immediately abort the scan if the patient reports burning sensation or pain.

GALAXY UNYCO can only be guaranteed for MRI when using the following components to build the frames:

Tibia frame

2 x Large Multiscrew Clamps for UNYCO Screws

1 x Rod Ø 12mm L 350mm

Up to 8 within these screws available in the kit:

8 x UNYCO Screw QC Shaft Ø 6mm

4 x UNYCO Cancellous Screw QC Shaft Ø 6mm

Ankle spanning – Delta frame

1 x Large Multiscrew Clamps for Unyco Screws

1 x Transfixing Pin Thread L80mm Shaft Ø 4mm, Thread Ø 5mm

2 x Rod Ø 12mm L 350mm

4 x UNYCO Cancellous Screw QC Shaft Ø 6mm

1 x Galaxy Large Clamp

2 x Galaxy Large-Medium Transition Clamp

Ankle spanning – Unilateral frame

1 x Large Multiscrew Clamps for Unyco Screws


1 x Radiolucent Foot Unit

2 x Rod Ø 12mm L 350mm

7 x UNYCO Cancellous Screw QC Shaft Ø 6mm

6 x Galaxy Large Clamp

Tested clamp configurations: M+Z, U+Z, M+S, and U+S as illustrated below

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS OF THE DEVICE

EXPECTED CLINICAL BENEFITS

- Temporary stabilization when immediate definitive treatment is not preferable
- Minimally invasive
- Ideal for DCO (Damage Control Orthopedics) in ready to use sterile kit
- Reduced X-ray exposure during application
- Skeletal stabilization in hemodynamically or systematically compromised polytrauma

PERFORMANCE CHARACTERISTICS OF THE DEVICE

- Designed to avoid contamination of the medullary canal
- Designed to facilitate conversion from temporary to definitive fixation
- Designed for a minimally invasive approach
- Fewer steps in the operative technique
- Versatility
- Compatible with GALAXY FIXATION system
- Fast and easy
- Stability
- Sterile Kits ready to be used

SPECIFIC INFORMATION ON THE DEVICE

- In case of knee and/or ankle spanning, stabilization in the femur and foot must be performed with bicortical screws in conjunction with GALAXY FIXATION.
- For more stable fixation of a fracture with a fixator, we recommend that the nearest bone screw is applied fairly close to the fracture margin (a minimum distance of 1cm is recommended) and that these distances are equal on both sides of the fracture.
- The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. It is recommended to apply the Large Multiscrew Clamp for UNYCO Screws 4cm from the bone.

- Make sure the Torque Limiter is coaxially inserted within the chuck holder of the power drill. During insertion of Uni-cortical screws stop drilling when the Torque Limiter stops turning and a rattle noise is heard. In cancellous bone use reference line on the screw to ensure that the screw is inserted maximum 15mm. Stop drilling when this insertion depth is reached or when the torque limiter stops turning, whichever comes first. The power drill must be used at its lowest speed.
- Depending upon the clinical and radiological findings, the surgeon will decide on the number of UNYCO Screws (93507) and UNYCO Cancellous Screws (93508) needed to achieve the appropriate frame stability.
- Additional equipment may be required for fixation application and removal.
- Tibia and Ankle Kits can also be used with bicortical screws.

RISKS DUE TO THE RE-USE OF “SINGLE USE” DEVICE

Implantable Device*

The “SINGLE USE” implantable device* of Orthofix is identified through symbol “⊗” reported on the product label. After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*) Implantable device: any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

Non Implantable Device

The “SINGLE USE” non implantable device of Orthofix is identified through symbol “⊗” reported on the label or indicated in the “Instructions For Use” supplied with the products. The re-use of a “SINGLE USE” non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such. In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged. The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled “FOR SINGLE USE ONLY” can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization release.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Aluminum based devices are damaged by alkaline (pH>7) detergents and solutions. Refer to PQALU for the list of Orthofix aluminum based devices.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- The contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing in order to remove soiling that accumulates in recesses.
- If a device needs particular care in pre-cleaning, a product specific IFU is available in Orthofix website that is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on reusable fixators and instruments
- End of life is normally determined by wear and damage due to use
- Products labeled for Single Use only MUST NOT be reused regardless any reprocessing in a clinical setting

POINT OF USE

It is recommended to reprocess the reusable medical devices as soon as is reasonably practicable in order to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because this can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

It is recommended to cover contaminated instruments during transportation in order to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated.

Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In case of highly contaminated reusable medical device before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described in the following paragraph) are recommended.

Manual Pre-cleaning

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution in order to displace trapped air.
5. Scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using a clean, lint free cloth.

CLEANING

General considerations

In these instructions Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to cleaning agents used. Staff shall follow the safety precautions to comply with the procedure of the healthcare facility using protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent's manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment following safety precautions to comply with the procedure of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft brush until all visible soiling is removed. Use a soft bristled nylon brush to remove residuals from lumens, rough or complex surfaces using a twisting motion.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in ultrasonic device with degassed cleaning solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends on the basis of the validation performed to use an ultrasound frequency of 35kHz, power 300W_{eff}, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent's manufacturer technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present it is possible to use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device and had to be removed with the brush, the cleaning step must be repeated as described above.
14. Carefully hand-dry using a clean, lint-free cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution in order to displace trapped air; it is important to ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using a clean, lint-free cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces

2. Use a Washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry, no visible foreign material can be present. Ensure that the washer-disinfector and all services are operational.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices in order to locate the cannulations in vertical position and blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 min.
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of the detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 min at 55°C.
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 min.
 - d. Final rinsing with deionized water for 3 min.
 - e. Thermal disinfection at least 90°C or 194°F (max 95°C or 203°F) for 5 minutes or until A0=3000 is reached, The water used for thermal disinfection must be purified.
 - f. Drying at 110°C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.

The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent's manufacturer technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion on the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfector when it complete the cycle.
13. If necessary, drain off excess water and dry by using clean, linf.free cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use. All functional checks and inspections described below also cover the interfaces with other instruments or components. The failure modes below may be caused by end of life of the product, improper use or improper maintenance. Orthofix does not typically specify the maximum number of uses for re-useable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use. Specific instructions may be available for some product codes. These instructions are linked to the product code and are available in a dedicated Orthofix website. Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained to mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non-woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS). The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instruction from the manufacturer. When equipment and processes differ from those validated by Orthofix the health care facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products. Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always to use a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated but it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers but only in wraps.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations. These cleaning agents are not listed in preference to other available cleaning agents which may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0,5%

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

ADDITIONAL INFORMATION

INFORMATION FOR THE PATIENT

The HCP shall inform the patient that the medical device does not replicate a normal healthy bone and counsel him/her about the correct behavior after implantation. The patient must pay attention to premature weight bearing, load carrying and excessive activity levels. The HCP shall inform the patient about any known or possible restrictions regarding the exposure to reasonably foreseeable external influences or environmental condition and in undergoing specific diagnostic investigations, evaluation or therapeutic treatment after implantation. The HCP shall inform the patient about the need for periodic medical follow-up and of the eventual removal of the medical device in the future. The HCP shall warn the patient about the surgical and residual risks and make him/her aware of possible adverse events. A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure that require surgical re-intervention to remove or replace the medical device. The HCP shall instruct the patient to report any unusual changes to the operative site or in device performance to his/her physician.

NOTICE ABOUT SERIOUS INCIDENTS

Report any serious incident involving a device to Orthofix Srl and the appropriate governing body in which the user and/or patient is established.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

MANUFACTURER CONTACT

Please contact your local Orthofix sales representative for further details and ordering.

Symbols presented below may apply or not to a specific product: refer to its label for applicability.

Symbol	Description	
	Medical Device	
	Consult instructions for use or consult electronic instructions for use	Caution: Consult instructions for use for important cautionary information
	Single Use. Do not re-use	Orthofix note: discard appropriately after the use (treatment) on patient
	Do not re-sterilize	
	Sterilized using irradiation	
	Double sterile barrier system	
	Catalogue number	Batch code
	Use-by date (year-month-day)	
	CE marking in conformity to applicable European Medical Device Directives/Regulations	
	Date of manufacture	Manufacturer
	Do not use if package is damaged and consult instructions for use	
	Symbol for MR conditional. The item has been demonstrated to pose no known hazards in a specified MRI environment with specified conditions of use.	
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	
	Unique Device Identifier	

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ORTHOFIX Srl