

Agile Nail™

System



**Very Small,
Absolutely Strong**





Not an actual patient

For femoral fractures and deformity treatments of adolescents and small adults

Diameter from 7mm up to 10mm

Double locking option for the proximal and the distal screws

(static or dynamic)



Procurvation design for easy insertion

Jig with patented locking mechanism





Benefits to Surgeons


- Diameter from 7mm up to 10mm
- Titanium alloy implants
- Procurvation design for easy insertion
- Optimized and lean instrumentation
- Jig with patented locking mechanism



Benefits to Patients

- Minimally invasive (with a small diameter for the nails and the screws)
- Early weight-bearing as tolerated by the patient and under surgeon discretion
- Titanium alloy implants to avoid allergic reaction to nickel
- 10° proximal bending to facilitate the lateral insertion of the nail and to reduce the impact on the growth plate
- Small proximal diameter of the nail for reduced invasiveness

**Anodization to increase
the fatigue resistance**



All implants available in sterile packaging

Agile Nail™ - Left, Sterile				
Nail Length (mm)	ø7	ø8	ø9	ø10
180	99-T837180L			
200	99-T837200L			
220	99-T837220L			
240	99-T837240L	99-T838240L	99-T839240L	99-T830240L
260	99-T837260L	99-T838260L	99-T839260L	99-T830260L
280	99-T837280L	99-T838280L	99-T839280L	99-T830280L
300	99-T837300L	99-T838300L	99-T839300L	99-T830300L
320	99-T837320L	99-T838320L	99-T839320L	99-T830320L
340	99-T837340L	99-T838340L	99-T839340L	99-T830340L
360	99-T837360L	99-T838360L	99-T839360L	99-T830360L
380	99-T837380L	99-T838380L	99-T839380L	99-T830380L
400		99-T838400L	99-T839400L	99-T830400L

Agile Nail™ - Right, Sterile				
Nail Length (mm)	ø7	ø8	ø9	ø10
180	99-T837180R			
200	99-T837200R			
220	99-T837220R			
240	99-T837240R	99-T838240R	99-T839240R	99-T830240R
260	99-T837260R	99-T838260R	99-T839260R	99-T830260R
280	99-T837280R	99-T838280R	99-T839280R	99-T830280R
300	99-T837300R	99-T838300R	99-T839300R	99-T830300R
320	99-T837320R	99-T838320R	99-T839320R	99-T830320R
340	99-T837340R	99-T838340R	99-T839340R	99-T830340R
360	99-T837360R	99-T838360R	99-T839360R	99-T830360R
380	99-T837380R	99-T838380R	99-T839380R	99-T830380R
400		99-T838400R	99-T839400R	99-T830400R



4.0mm Locking Screws, Sterile	
Length (mm)	Description
18	99-T834018
20	99-T834020
22	99-T834022
24	99-T834024
26	99-T834026
28	99-T834028
30	99-T834030
32	99-T834032
34	99-T834034
36	99-T834036
38	99-T834038
40	99-T834040
42	99-T834042
44	99-T834044
46	99-T834046
48	99-T834048
50	99-T834050
55	99-T834055
60	99-T834060
65	99-T834065
70	99-T834070
75	99-T834075
80	99-T834080
85	99-T834085

3.2mm Peg Screws, Sterile (For 7mm Solid Nail Distal Locking only)	
Length (mm)	Description
18	99-T833018
20	99-T833020
22	99-T833022
24	99-T833024
26	99-T833026
28	99-T833028
30	99-T833030
32	99-T833032
34	99-T833034
36	99-T833036
38	99-T833038
40	99-T833040
42	99-T833042
44	99-T833044
46	99-T833046
48	99-T833048
50	99-T833050
55	99-T833055
60	99-T833060

End Caps	
Length (mm)	Description
0	99-T830000
5	99-T830005
10	99-T830010
15	99-T830015
20	99-T830020

183990 Agile Nail™ Instrument Box, Empty

Part #	Description
183990	Agile Nail Sterilization Box, Empty
Can accommodate:	
183260	1 Cannulated Awl
183230	2 Entry Reamer Guide
193325	3 6mm Hex Screwdriver
193271	4 6mm Hex Extension
183100	5 Targeting Handle
183212	6 Trocar
183320	7 Short Screwdriver
183321	8 Screwdriver
183277	9 Distal Screw Sizing Gauge
183213	10 Distal Tissue Guide
183285	11 Distal Drill Bit D 3.2mm
183270	12 Entry Reamer
183275	13 Sizing Guide
177380	14 Slotted Mallet
193115	15 Impaction Rod
183110	16 Locking Bolt
183302	17 Tissue Guide
183274	18 Screw Sizing Gauge
183264	19 Reduction Tool
173276	20 Ruler Support
183276	21 Ruler
80	99-T834080
85	99-T834085




Sterile packaged instruments


Part #	Description
99-183288	Guide Wire L 300mm D 3.2mm
99-183281	Guide Wire with olive L 980mm D 2.8mm

183991 Agile Nail™ Extraction Box, Empty

Part #	Description
183991	Agile Nail Sterilization Box, Extraction, Empty
Can accommodate:	
193325	1 6mm Hex Screwdriver
193271	2 6mm Hex Extension
183320	3 Short Screwdriver
183336	4 M7 Conical Mallet Adaptor
183978	5 Plier Cutter D 15mm
183337	6 Screw Extractor Size 3.5-9mm
177395	7 Needle Nose Vise Grips
SMN173370	8 Slap Hammer
17955	9 Universal Chuck with T-Handle
91017	10 Universal Allen Wrench
177380	11 Slotted Mallet



The Complete Orthofix Pediatric Solution Portfolio



TL-HEX TrueLok Hexapod System™



Galaxy Fixation™ System



Fitbone™ Intramedullary Lengthening System



Agile Nail™




JuniOrtho™ Plating System



Eight-Plate Guided Growth System +™




TrueLok™ EVO Ring Fixation System



MJ-Flex
The New Metaizeau Nail™



LRS Pediatric



Small Daf



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo
(Verona), Italy
Tel.: +39 045 6719000
Fax: +39 045 6719380

Rx Only



Orthofix products or services referenced herein are trademarks or registered trademarks of Orthofix SRL and its group of companies. All rights reserved.

TrueLok is a trademark or registered trademark exclusively licensed to Orthofix Medical Inc., and its affiliate companies. All rights reserved.

Electronic Instructions for use available at the website

<http://ifu.orthofix.it>

Electronic Instructions for use - Minimum requirements for consultation:

- Internet connection (56 kbps)
- Device capable to visualize PDF (ISO/IEC 32000-1) files
- Disk space: 50Mbytes

Free paper copy can be requested to customer service (delivery within 7 days):

tel.: +39 045 6719301, fax: +39 045 6719370

e-mail: customerservice@orthofix.it

Distributed by:



Proper surgical procedure is the responsibility of the medical professional. This Manual is furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the Agile Nail System Instructions for Use (PQ ANS) supplied with the products for specific information on indications for use, contraindications, warnings, precautions, adverse effects and sterilization.

Orthofix.com
OrthofixKids.com



AG-1701-PL-E0 E 04/23

ORTHOFIX®
ORTHOPEDICS