


Galaxy Fixation[®]

System

In a world of trauma
a Galaxy of solutions



Components

The Galaxy Fixation components are MR conditional when marked with "MR".

Rods Diameter 12mm



Code	Description
932100	Rod 100mm long
932150	Rod 150mm long
932200	Rod 200mm long
932250	Rod 250mm long
932300	Rod 300mm long
932350	Rod 350mm long
932400	Rod 400mm long
99-932450	Rod 450mm long, sterile*
99-932500	Rod 500mm long, sterile*
99-932550	Rod 550mm long, sterile*
99-932600	Rod 600mm long, sterile*
99-932650	Rod 650mm long, sterile*

* Special order only

Rods Diameter 9mm



Code	Description
939100	Rod 100mm long
939150	Rod 150mm long
939200	Rod 200mm long
939250	Rod 250mm long
939300	Rod 300mm long

Rods Diameter 6mm



Code	Description
936060	Rod 60mm long
936080	Rod 80mm long
936100	Rod 100mm long
936120	Rod 120mm long
936140	Rod 140mm long
936160	Rod 160mm long
936180	Rod 180mm long
936200	Rod 200mm long

L Rod Diameter 6mm



Code	Description
936010	6mm L Rod

Semi-Circular Rods Diam. 12mm (Aluminium)



Code	Description
932010	Semi-Circular Rod small 180mm long
932020	Semi-Circular Rod medium 215mm long
932030	Semi-Circular Rod large 250mm long

Semi-Circular Rods Diam. 9mm

939010	Semi-Circular Rod small 115mm long
939020	Semi-Circular Rod medium 140mm long
939030	Semi-circular Rod large 165mm long

Note: All components of Galaxy Fixation System frames must be identified as MR Conditional prior to being placed in or near an MR Environment.

The Threaded Wires (93100), the Wire Locking Clamps (93620), and the L Rod (936010) and Semi-Circular Rods (939010, 939020, 939030) are not MR Conditional. Any construct/frame that is using Threaded Wires, the Wire Locking Clamps, the L-Rod and Semi-Circular Rods must therefore be considered as MR Unsafe.

THE SYSTEM IS MRI CONDITIONAL at 1.5 and 3 Tesla
More details are available on relevant Instructions
For Use - PQGAL



Orthofix ABS – A touch of Orthofix Knowledge

We are now covering all the mobile devices:

iOS version (App Store) for iPad


Android version (Google Play Store) for Android tablet

Website: <https://abs.orthofix.it/> for smartphone and PC


Components

The Galaxy Fixation components are MR conditional when marked with "MR".

Clamps



Multiscrew Clamps



Modules






HUMERUS

N.	Code	Description
3	99-93100	Threaded Wire Sterile (pack of 2)
3	93620	Wire Locking Clamp
4	93310	Small Clamp
1	936080	Rod Ø 6mm L 80mm
1	936140	Rod Ø 6mm L 140mm



N.	Code
6	93110
2	939100
1	939150
4	911530




WRIST

N.	Code	Description
5	93310	Small Clamp
1	936140	Rod Ø 6mm L 140mm
1	936080	Rod Ø 6mm L 80mm
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC




N.	Code
6	93310
1	936100
1	936080
1	936060
4	947320



PELVIS

N.	Code	Description
8	93010	Large Clamp
2	932200	Rod Ø 12mm L 200mm
1	932300	Rod Ø 12mm L 300mm
1	932400	Rod Ø 12mm L 400mm
4	912650	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 260/50



FEMUR

N.	Code
6	93010
1	932100
2	932200
4	912640

Simple components,

Description

Medium Clamp
 Rod Ø 9mm L 100mm
 Rod Ø 9mm L 150mm
 Self-drilling XCaliber Screws
 Shaft Ø 6mm
 Thread Ø 6.0-5.6mm
 150/30



ELBOW

N.	Code	Description
1	93020	Multiscrew Clamp
1	932200	Rod Ø 12mm L 200mm
1	93410	Elbow Hinge
1	939150	Rod Ø 9mm L 150mm
2	93110	Medium Clamp
2	911530	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/30
2	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20

Description

Small Clamp
 Rod Ø 6mm L 100mm
 Rod Ø 6mm L 80mm
 Rod Ø 6mm L 60mm
 Self Drilling XCaliber
 Cylindrical Screw
 Shaft Ø 4mm
 Thread 3mm L 100/20 QC



N.	Code	Description
1	93350	Wrist Module
2	93330	Small Multiscrew Clamp-Short
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC


Description

Large Clamp
 Rod Ø 12mm L 100mm
 Rod Ø 12mm L 200mm
 XCaliber Screws
 260/40 6.0-5.6mm



KNEE

N.	Code	Description
8	93010	Large Clamp
2	932350	Rod Ø 12mm L 350mm
1	932250	Rod Ø 12mm L 250mm
4	911540	XCaliber Screw 150/40 6-5.6mm
2	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50



FOOT & ANKLE

N.	Code	Description
8	93010	Large Clamp
1	932200	Rod Ø 12mm L 200mm
1	932100	Rod Ø 12mm L 100mm
1	932250	Rod Ø 12mm L 250mm
1	932400	Rod Ø 12mm L 400mm
1	911540	XCaliber Screw 150/40 6-5.6mm
2	911530	XCaliber Screw 150/30 6-5.6mm

Infinite solutions



FOREARM

N.	Code	Description
5	93110	Medium Clamp
1	939100	Rod Ø 9mm L 100mm
1	939150	Rod Ø 9mm L 150mm
4	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20



N.	Code	Description
1	936010	6mm L Rod
1	93330	Small Multiscrew Clamp-Short
1	93320	Small Multiscrew Clamp-Long
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC




PELVIS



TIBIA

N.	Code	Description
6	93010	Large Clamp
1	932100	Rod Ø 12mm L 100mm
1	932150	Rod Ø 12mm L 150mm
1	932200	Rod Ø 12mm L 200mm
4	911540	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/40




N.	Code	Description
4	93010	Large Clamp
2	99-93030	Large-Medium Transition Clamp (Sterile)
1	932350	Rod Ø 12mm L 350mm
1	932250	Rod Ø 12mm L 250mm
1	932200	Rod Ø 12mm L 200mm
1	92080	Transfix Pin 80mm, Ø 4mm
2	911540	XCaliber Screw 150/40 6-5.6mm

Possibility to use six Large Clamps
and Ø 6mm Transfix Pin (93080)



N.	Code	Description
6	93010	Large Clamp
2	932250	Rod Ø 12mm L 250mm
1	932350	Rod Ø 12mm L 350mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50



N.	Code	Description
9	93010	Large Clamp
3	932300	Rod Ø 12mm L 300mm
2	932250	Rod Ø 12mm L 250mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50




N.	Code	Description
1	99-93794	Galaxy UNYCO Diaphyseal Tibia Box
		or
1	99-93574	Galaxy UNYCO Mini Kit Tibia Sterile
2	99-932350	Rod D12mm L 350mm Sterile
1	99-93509	Galaxy UNYCO Mini Kit Instruments Sterile


Get it all sterile!

Besides the sterile pre-packaged kits, Galaxy Fixation System offers all **clamps and rods single-packed and sterile**, except for the 6mm L-rod cod. 936010. They can be ordered using the item code number preceded by 99- (e.g. 99-939300).

99-93505 Galaxy Shoulder Sterile Kit



99-93504 Galaxy Elbow Sterile Kit




99-93501 Pelvis Sterile Kit



99-93498 Galaxy Small Z-Configuration Sterile Kit




99-93502 Lower Limb Diaphyseal Sterile Kit




99-93499 Ankle Sterile kit - Transfix Pin Ø4mm




99-93503 Ankle Sterile Kit




99-93520 Galaxy Medium - Paediatric Sterile Kit - Ø5mm screw thread



99-93601 Galaxy Wrist Sterile Kit



99-93521 Galaxy Medium - Paediatric Sterile Kit - Ø4mm screw thread




Manufactured by:
 ORTHOFIX Srl
 Via Delle Nazioni 9, 37012
 Bussolengo (Verona), Italy
 Telephone +39 045 6719000,
 Fax +39 045 6719380

orthofix.com

GF-1202-PL-E0 P 2/25

Distributed by:



 **ORTHOFIX**[®]
 ORTHOPEDICS