


Galaxy Fixation™

Wrist System

In a world of trauma
a Galaxy of solutions

- Modular components compatible with the Galaxy Fixation System
- Early mobilisation option for quick wrist function recovery
- Allows for application of combined internal synthesis (i.e. K-wires, FFS, plates and screws)
- Radiolucent

Bridging



Non Bridging



Motion range

selector ($\pm 20^\circ / \pm 40^\circ$) for controlled, limited flexion-extension

Joint mobilization

Option (up to $\pm 40^\circ$) for improved functional recovery and patient comfort

Compression and distraction option



Radiolucency

for fracture site visibility


MRI Conditional at 1.5 and 3 Tesla

More details are available on relevant Instructions For Use - PQGAL

Comprehensive

stand-alone instrumentation tray

93999 Galaxy Wrist Tray, Empty




99-93601 Galaxy Wrist Sterile Kit



For details about Tray content please refer to *Operative Technique GF-1101-OPT* or *Orthofix General Catalogue PX GEN*

Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo
(Verona), Italy
Tel.: +39 045 6719000
Fax: +39 045 6719380

Rx Only



Orthofix products or services referenced herein are trademarks or registered trademarks of Orthofix SRL and its group of companies. All rights reserved.

www.orthofix.com



GF-1309-PL-E0 D 03/22

ORTHOFIX[®]
ORTHOPEDICS