


In a world of trauma a Galaxy of solutions



- Modular components compatible with the Galaxy Fixation System
- Early mobilisation option for quick wrist function recovery
- Allows for application of combined internal synthesis (i.e. K-wires, FFS, plates and screws)
- Radiolucent

Bridging



Non Bridging



Joint

Motion range

selector ($\pm 20^\circ / \pm 40^\circ$) for controlled, limited flexion-extension

mobilisation

Option (up to $\pm 40^\circ$) for improved functional recovery and patient comfort

Compression and distraction option



MRI Conditional at 1.5 and 3 Tesla
More details are available on relevant Instructions For Use - PQGAL

Radiolucency

for fracture site visibility

Comprehensive

stand-alone instrumentation tray


93999C Galaxy Wrist Tray



99-93601 Galaxy Wrist Sterile Kit



Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

www.orthofix.com

GF-1309-PL-E0 C 12/16


ORTHOFIX®