


PRODUCT SHEET

Galaxy Fixation™

Shoulder System



**Minimally
invasive and
stable. How?**


Simple, standardized and reproducible

frame application for most common proximal humeral fracture patterns*

Minimally invasive

six wires externally stabilized with Galaxy Shoulder components

Subchondral and lateral bone purchase




Stable

7cm threaded paired bicortical wires

ORDERING INFORMATION

Description	Code
Sterilisation Tray	93998C
Sterile Kit	99-93505
Threaded wire sterile (pack of 2)	99-93100

Sterilisation Tray




Sterile Kit



Case Study

Two-part fracture (AO 11-A2)



Post operative



4 months FU with good functional recovery



Case study provided by Dr. D.Blonna MD, Mauriziano "Umberto I" University Hospital, Turin, Italy

References

- * See IFU PQ GAL

Prodotto da:
ORTHOFIX Srl
 Via Delle Nazioni 9, 37012
 Bussolengo (Verona) Italia
 Telefono +39 045 6719000
 Fax +39 045 6719380
 www.orthofix.com



www.galaxyfixation.com

orthofix.com
 GF-1801-PL-EO C 07/22

Distributed by:



Instructions for Use: See actual package insert for Instructions for Use.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for specific information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.

 **ORTHOFIX**[®]
 ORTHOPEDICS