



**UNYCO**™ **GALAXY**


One Giant Leap  
for Damage Control Orthopaedics

Designed to  
**facilitate the conversion  
from temporary  
to definitive fixation**



Conversion into IM Nail with Galaxy UNYCO in situ


Designed to  
**avoid contamination  
of the medullary canal**



Uni-cortical screws

Designed for a  
**minimally invasive  
approach**

**Compatible with  
Galaxy Fixation System**



## Stability

By clustering minimum 3  
Uni-cortical screws together



The System is **MRI Conditional** at 1.5 and 3 Tesla  
(refer to the relevant IFU for MRI safety information and Indications For Use)

## Fast Locking

Bar snap in and pre-closure by hand; definitive locking in one step

## Quick Connection




## Flexibility

$\pm 10^\circ$  polyaxial movement on each screw



## Speed and Ease of Use

Power torque limited insertion for correct screw application



## For the Surgeon

Ideal for emergency situations  
in ready-to-use sterile kits

## For the Hospital

Sterile packs provide cost savings  
and a simple, comprehensive stock  
management while helping with  
accurate budget planning



## For the Sterilization Department

Minimum amount of items  
to be reprocessed

## For the Patient

Fewer steps in the operative  
technique and designed for a  
minimally invasive approach

All sterile components can be ordered separately or  
grouped by procedure through the **Galaxy UNYCO Boxes**



99-93574  
Galaxy UNYCO  
Mini Kit Tibia Sterile

- 1x Large Multiscrew Clamp
- 4x UNYCO Cancellous Screw



99-93509  
Galaxy UNYCO  
Mini Kit Instruments Sterile

- 1x Allen Wrench
- 1x Power Drill Torque Limiter



99-93477  
Galaxy UNYCO  
Mini Kit Ankle Unilateral Sterile

- 1x Radiolucent Foot Unit
- 3x UNYCO Cancellous Screw



99-93570  
Galaxy UNYCO  
Cancellous Screw Kit Sterile

- 2x UNYCO Cancellous Screw

99-93794  
Galaxy UNYCO  
Diaphyseal Tibia Box

- 2x Mini Kit Tibia Sterile
- 1x Mini Kit Instruments Sterile
- 1x Galaxy Rod Ø 12mm  
L 350mm sterile

99-93793  
Galaxy UNYCO  
Ankle Bridging Box -  
Delta Frame

- 1x Mini Kit Tibia Sterile
- 1x Mini Kit Instruments Sterile
- 2x Transitional Clamp Sterile
- 1x Large Single Clamp Sterile
- 2x Galaxy Rod Ø 12mm  
L 350mm sterile
- 1x Transfixing Pin Ø 4mm Sterile

99-93792  
Galaxy UNYCO  
Ankle Bridging Box -  
Unilateral Frame

- 1x Mini Kit Tibia Sterile
- 1x Mini Kit Ankle  
Unilateral Sterile
- 1x Mini Kit Instruments Sterile
- 6x Large Single Clamp Sterile
- 2x Galaxy Rod Ø 12mm  
L 350mm sterile

Distributed by:



Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380



@galaxyunyco  
[www.galaxyunyco.com](http://www.galaxyunyco.com)

[www.orthofix.com](http://www.orthofix.com)

MC-1510-PL-E0 A 10/15

Instructions for Use: See actual package insert for Instructions for Use.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for specific information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.

