

ORTHOFIX TRAUMA STERILE KITS



HUMERUS | ELBOW | FOREARM



WRIST | HAND



PELVIS




FEMUR | KNEE | TIBIA



FOOT & ANKLE



PEDIATRICS



Summary


● Humerus	03
● Elbow	03
● Forearm	03
● Wrist & Hand	04
● Pelvis	05
● Femur	05
● Knee	06
● Tibia	07
● Foot & Ankle	08
● Pediatrics	10

EXTERNAL FIXATION INDICATIONS IN TRAUMA

The aim of an external fixator for **damage control** is to provide skeletal stability and for this to be achieved efficiently; that is, **the objective is accomplished quickly and does not hinder future definitive care.**

The Galaxy Fixation™ and Galaxy UNYCO Systems have been designed with this characteristic in mind. While Galaxy UNYCO has been specifically developed for damage control and emergency - therefore it is purely for temporary stabilisation - Galaxy Fixation System can also be used for definitive treatment of long bone fractures. The external fixator configuration will depend on the amount of bone contact at the fracture site, the fracture pattern and the segment or segments of bone involved.

Humerus



99-93505 Galaxy Shoulder Sterile Kit


Code	Description
4x93310	SMALL CLAMP
3x93620	WIRE LOCKING CLAMP
1x936080	ROD D 6MM L 80MM
1x936100	ROD D 6MM L 100MM
1x936120	ROD D 6MM L 120MM
1x936140	ROD D 6MM L 140MM
1x30017	ALLEN WRENCH 5MM
1x19970	WIRE GUIDE
1x91150	UNIVERSAL "T" WRENCH
1x81031	OPEN END WRENCH
6x93100	300MM THREADED WIRE



99-91600UE XCaliber Metadiaphyseal Complete Sterile Kit

Code	Description
1x91150	UNIVERSAL T WRENCH
1x91017	UNIVERSAL ALLEN WRENCH
1x11106	DRILL GUIDE 3.2X40MM
1x11104	DRILL GUIDE 4.8X40MM
1x91000	XCALIBER META-DIAPHYSEAL FIXATOR
3x11102	SCREW GUIDE L 60MM
1x91015	DYNAMIC COMPRESSION/DISTRACTION UNIT
1x99-611560	XCALIBER OSTEOTITE SCREW 150/60
1x99-611550	XCALIBER OSTEOTITE SCREW 150/50
1x99-611540	XCALIBER OSTEOTITE SCREW 150/40
3x99-611530	XCALIBER OSTEOTITE SCREWS 150/30
1x1101101	CANNULATED DRILL BIT 3.2X200MM
2x11014	X-WIRES WITHOUT OLIVE 1.5X250MM
1x1100101	DRILL BIT 4.8X180MM

Elbow



99-93504 Galaxy Elbow Sterile Kit

Code	Description
1x93020	MULTISCREW CLAMP
2x93110	MEDIUM CLAMP
1x93410	ELBOW HINGE
1x932200	ROD Ø 12MM L 200MM
1x939150	ROD Ø 9MM L 150MM
1x30017	ALLEN WRENCH 5MM
1x91150	UNIVERSAL "T" WRENCH
1x1-1100201	DRILL BIT Ø 4.8X240MM
1x11138	DRILL GUIDE Ø 4.8MM
2x11137	SCREW GUIDE 80MM
1x19950	DRILL GUIDE Ø 3.2MM L 100MM
1x1-1100301	DRILL BIT Ø 3,2X200MM QC
2x911530	XCALIBER SCREWS 150/30
2x10137	CORTICAL SCREWS 120/20 Ø 4,5/3,5 SHAFT 6MM
1x11146	X-WIRE WITHOUT OLIVE Ø 2MM LENGTH 150MM

Forearm



99-93498 Galaxy Small Z-Configuration Sterile Kit

Code	Description
6x93310	SMALL CLAMP
1x936120	ROD Ø 6MM L 120MM
1x936100	ROD Ø 6MM L 100MM
1x936080	ROD Ø 6MM L 80MM
2x19995	SLEEVE 4.5MM
2x13530	DRILL GUIDE D 2.7MM
1x91017	ALLEN WRENCH
4x947320	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 4MM THREAD 3MM L 100/20 QC
2x1-1355001	DRILL BIT D 2,7MM
1x93160	"T" WRENCH AO QC

Wrist & Hand



Galaxy Fixation™
System

99-93498 Galaxy Small Z-Configuration Sterile Kit

Code	Description
6x93310	SMALL CLAMP
1x936120	ROD Ø 6MM L 120MM
1x936100	ROD Ø 6MM L 100MM
1x936080	ROD Ø 6MM L 80MM
2x19995	SLEEVE 4.5MM
2x13530	DRILL GUIDE D 2.7MM
1x91017	ALLEN WRENCH
4x947320	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 4MM THREAD 3MM L 100/20 QC
2x1-1355001	DRILL BIT D 2,7MM
1x93160	"T" WRENCH AO QC



Galaxy Fixation™
System

99-93601 Galaxy Wrist Sterile Kit


Code	Description
2x93330	SMALL MULTISCREW CLAMP-SHORT
1x93350	WRIST MODULE
2x19995	SCREW GUIDE
2x13530	DRILL GUIDE Ø 2.7MM
1x91017	ALLEN WRENCH
4x947320	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 4MM THREAD 3MM L 100/20 QC
2x1-1355001	DRILL BIT Ø 27MM
1x93160	SCREW T WRENCH QC
3x13715	KWIRE Ø 1,5X150MM

Wrist & Hand



99-36501 RadioLucent Wrist Fixator Sterile Kit

Code	Description
1x	RADIO-TRANSPARENT WRIST FIXATOR
1x	COMPRESSION-DISTRACTION UNIT
1x	ALLEN WRENCH 4MM
2x	DRILL GUIDES
1x	DRILL BIT Ø 2.7MM
1x	SCREW TEMPLATE
1x	T-WRENCH
2x	SELF-DRILLING SCREWS, SHAFT Ø 4MM, THREAD Ø 3.3-3MM, TOTAL LENGTH 70MM, THREAD LENGTH 20MM
2x	SELF-DRILLING SCREWS, SHAFT Ø 4MM, THREAD Ø 3.3-3MM, TOTAL LENGTH 80MM, THREAD LENGTH 25MM
4x	BONE SCREW COVERS





Pennig
Minifixator

99-MIL10 Mini-Pennig Fixator Std Kit Complete Sterile

Qty	Description
1x	ALLEN WRENCH, 3 MM
1x	THREADED, WIRE EXTRACTOR
6x	PMF THREADED WIRES D 1.6 MM L 70/15 MM
1x	PENNIG MINIFIXATOR, STD, BODY
1x	CLAMPS, STANDARD, (ONE PAIR) FOR PENNIG MINIFIXATOR
1x	NUT, COMPRESSION/DISTRACTION
2x	PENNIG MINIFIXATOR SPACER
1x	LENGTHENING BAR, STD, 100 MM FOR PENNIG MINIFIXATOR
6x	PMF THREADED WIRES D 2.0 MM L 100/15 MM
1x	REDUCTION FORCEPS (2 PAIRS) FOR PENNIG MINIFIXATOR



Pelvis



Galaxy Fixation®
System

99-93501 Pelvis Sterile Kit

Code	Description
8x93010	LARGE CLAMP
1x932350	ROD Ø 12MM L 350MM
1x932300	ROD Ø 12MM L 300MM
2x932200	ROD Ø 12MM L 200MM
4x912640	SELF-DRILLING XCALIBER SCREWS, L 260MM, THREAD LENGTH 40MM
1x11138	DRILL GUIDE Ø 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT Ø 4.8X180MM COATED - QUICK CONNECT
1x91150	UNIVERSAL "T" WRENCH




Galaxy Fixation
Gemini® System

99-94710 Galaxy Fixation Gemini Pelvis Set Complete Sterile

Code	Description
8x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
2x932200	ROD D 12MM L 200MM
1x932300	ROD D 12MM L 300MM
1x932350	ROD D 12MM L 350MM
4x942540	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/40 QC
1x11137	SCREW GUIDE, L 80MM
1x93160	T-WRENCH QC


Femur



Galaxy Fixation®
System

99-93502 Lower Limb Diaphyseal Sterile Kit

Code	Description
6x93010	LARGE CLAMP
1x932300	ROD Ø 12MM L 300MM
2x932150	ROD Ø 12MM L 150MM
2x912650	SELF-DRILLING XCALIBER SCREWS, L 260MM, THREAD LENGTH 50MM
2x911550	SELF-DRILLING XCALIBER SCREWS, L 150MM, THREAD LENGTH 50MM
1x11138	DRILL GUIDE Ø 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT Ø 4.8X180MM COATED - QUICK CONNECT
1x91150	UNIVERSAL "T" WRENCH




Galaxy Fixation
Gemini® System

99-94720 Galaxy Fixation Gemini Z Frame Set Complete Sterile

Code	Description
6x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
2x932150	ROD D 12MM L 150MM
1x932300	ROD D 12MM L 300MM
2x 941550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 180/50 QC
2x942550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/50 QC
1x11137	SCREW GUIDE, L 80MM
1x93160	T-WRENCH QC

Femur



Galaxy Fixation
Gemini® System

99-94760 Galaxy Fixation Gemini Quadrilateral Set Complete Sterile*

Code	Description
2x94040	GALAXY FIXATION GEMINI DOUBLE MULTIPIN CLAMP
2x932400	ROD D 12MM L 400MM
6x942540	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/40 QC
2x11137	SCREW GUIDE, LENGTH 80MM
1x93160	T-WRENCH QC

* Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.




TrueLok™ EVO
Ring Fixation System

99-8800X TL-EVO Double Ring Frame Set Sterile*

Code	Description
2x99-882XXX	RX 5/8 MODULAR RING DXXXMM
1x99-881XXX	RX 3/8 MODULAR RING DXXXMM
3x99-88600X	TRUELOK EVO RX STRUT (LONG or MEDIUM) WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
2x99-885003	TLEVO BOLT L 16.5MM
2x99-885000	TLEVO NUT WITH WASHER

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-EVO operative technique TL-2113-OPT

Femur




TL-HEX™ Trauma
TrueLok Hexapod System

99-880XX TL-HEX Double Ring Frame Set Sterile*

Code	Description
2x99-56-214X0	TL-HEX 5/8 MODULAR RING DxxxMM
1x99-56-23XX0	TL-HEX 3/8 MODULAR RING DxxxMM STERILE
3x99-88600X	TRUELOK EVO RX STRUT (LONG or MEDIUM) WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
2x99-885003	TLEVO BOLT L 16.5MM
2x99-885000	TLEVO NUT WITH WASHER

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-HEX operative technique TH-1905-OPT

Knee




Galaxy Fixation®
System

99-93502 Lower Limb Diaphyseal Sterile Kit

Code	Description
6x93010	LARGE CLAMP
1x932300	ROD Ø 12MM L 300MM
2x932150	ROD Ø 12MM L 150MM
2x912650	SELF-DRILLING XCALIBER SCREWS, L 260MM, THREAD LENGTH 50MM
2x911550	SELF-DRILLING XCALIBER SCREWS, L 150MM, THREAD LENGTH 50MM
1x11138	DRILL GUIDE Ø 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT Ø 4.8X180MM COATED - QUICK CONNECT
1x91150	UNIVERSAL "T" WRENCH

Knee




Galaxy Fixation

Gemini™ System

99-94720 Galaxy Fixation Gemini Z Frame Set Complete Sterile

Code	Description
6x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
2x932150	ROD D 12MM L 150MM
1x932300	ROD D 12MM L 300MM
2x 941550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 180/50 QC
2x942550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/50 QC
1x11137	SCREW GUIDE, L 80MM
1x93160	T-WRENCH QC

Tibia




XCaliber™

Meta-Diaphyseal Fixator

99-91600UE XCaliber Metadiaphyseal Complete Sterile Kit

Code	Description
1x91150	UNIVERSAL T WRENCH
1x91017	UNIVERSAL ALLEN WRENCH
1x11106	DRILL GUIDE 3.2X40MM
1x11104	DRILL GUIDE 4.8X40MM
1x91000	XCALIBER META-DIAPHYSEAL FIXATOR
3x11102	SCREW GUIDE L 60MM
1x91015	DYNAMIC COMPRESSION/DISTRACTION UNIT
1x99-611560	XCALIBER OSTEOTITE SCREW 150/60
1x99-611550	XCALIBER OSTEOTITE SCREW 150/50
1x99-611540	XCALIBER OSTEOTITE SCREW 150/40
3x99-611530	XCALIBER OSTEOTITE SCREWS 150/30
1x1101101	CANNULATED DRILL BIT 3.2X200MM
2x11014	X-WIRES WITHOUT OLIVE 1.5X250MM
1x1100101	DRILL BIT 4.8X180MM

Tibia




Galaxy Fixation®

System

99-93502 Lower Limb Diaphyseal Sterile Kit

Code	Description
6x93010	LARGE CLAMP
1x932300	ROD Ø 12MM L 300MM
2x932150	ROD Ø 12MM L 150MM
2x912650	SELF-DRILLING XCALIBER SCREWS, L 260MM, THREAD LENGTH 50MM
2x911550	SELF-DRILLING XCALIBER SCREWS, L 150MM, THREAD LENGTH 50MM
1x11138	DRILL GUIDE Ø 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT Ø 4.8X180MM COATED - QUICK CONNECT
1x91150	UNIVERSAL "T" WRENCH




Galaxy Unyco™

System

99-93794 Galaxy UNYCO Diaphyseal Tibia Box

Code	Description
2x99-93574	GALAXY UNYCO MINI KIT TIBIA STERILE
1x99-939250	ROD Ø 12MM L 350MM STERILE
1x99-93509	GALAXY UNYCO MINI KIT INSTRUMENTS STERILE

Tibia




Galaxy Fixation

Gemini® System

99-94720 Galaxy Fixation Gemini Z Frame Set Complete Sterile

Code	Description
6x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
2x932150	ROD D 12MM L 150MM
1x932300	ROD D 12MM L 300MM
2x 941550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 180/50 QC
2x942550	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/50 QC
1x11137	SCREW GUIDE, L 80MM
1x93160	T-WRENCH QC



Galaxy Fixation

Gemini® System

99-94760 Galaxy Fixation Gemini Quadrilateral Set Complete Sterile*

Code	Description
2x94040	GALAXY FIXATION GEMINI DOUBLE MULTIPIN CLAMP
2x932400	ROD D 12MM L 400MM
6x942540	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 260/40 QC
2x11137	SCREW GUIDE, LENGTH 80MM
1x93160	T-WRENCH QC

* Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.

Tibia



TrueLok® EVO

Ring Fixation System

99-880XX TL-EVO Double Ring Frame Set Sterile*

Code	Description
2x99-882XXX	RX 5/8 MODULAR RING DXXXMM
1x99-881XXX	RX 3/8 MODULAR RING DXXXMM
3x99-88600X	TRUELOK EVO RX STRUT (LONG or MEDIUM) WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
2x99-885003	TLEVO BOLT L 16.5MM
2x99-885000	TLEVO NUT WITH WASHER

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-EVO operative technique TL-2113-OPT



TL-HEX® Trauma



TrueLok Hexapod System

99-880XX TL-HEX Double Ring Frame Set Sterile*

Code	Description
2x99-56-214X0	TL-HEX 5/8 MODULAR RING DXXXMM
1x99-56-23XX0	TL-HEX 3/8 MODULAR RING DXXXMM STERILE
3x99-88600X	TRUELOK EVO RX STRUT (LONG OR MEDIUM) WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
2x99-885003	TLEVO BOLT L 16.5MM
2x99-885000	TLEVO BOLT L 16.5MM

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-HEX operative technique TH-1905-OPT


Foot & Ankle



Galaxy Fixation®
System

99-93503 Ankle Sterile Kit - Transfix Pin Ø6mm

Code	Description
6x93010	LARGE CLAMP
1x932300	ROD D 12MM L 300MM
1x932200	ROD D 12MM L 200MM
1x932150	ROD D 12MM L 150MM
2x911540	SELF-DRILLING XCALIBER SCREWS, L 150MM, THREAD L 40MM
1x11138	DRILL GUIDE D 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT D 4.8X180MM COATED - QUICK CONNECT
1x91150	UNIVERSAL "T" WRENCH
1x93080	TRANSFIX PIN 80MM - SHAFT D 6MM/THREAD D 7MM




Galaxy Fixation®
System

99-93499 Ankle Sterile Kit - Transfix Pin Ø4mm

Code	Description
4x93010	LARGE CLAMP
2x93030	TRANSITIONAL CLAMP
1x932300	ROD Ø 12MM L 300MM
1x932200	ROD Ø 12MM L 200MM
1x932150	ROD Ø 12MM L 150MM
2x911540	XCALIBER SCREWS
1x11138	DRILL GUIDE D 4.8MM
1x11137	SCREW GUIDE 80MM
1x1-1100101	DRILL BIT Ø 4.8X180MM
1x91150	UNIVERSAL "T" WRENCH
1x92080	TRANSFIX PIN THREAD LENGHT 60MM, THREAD Ø 5MM, SHAFT 4MM


Foot & Ankle



Galaxy Fixation
Gemini® System

99-94730 Galaxy Fixation Gemini Ankle 6mm Set Complete Sterile

Code	Description
6x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
1x932150	ROD D 12MM L 150MM
1x932200	ROD D 12MM L 200MM
1x932300	ROD D 12MM L 300MM
2x941540	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 180/40 QC
1x1-93080	TRANSFIX PIN 80MM QC D6MM
1x11102	SCREW GUIDE, L 60MM
1x93160	T-WRENCH QC



Galaxy Fixation
Gemini® System

99-94740 Galaxy Fixation Gemini Ankle 4mm Set Complete Sterile

Code	Description
4x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
2x94030	GALAXY FIXATION GEMINI TRANSITIONAL SINGLE CLAMP
1x932150	ROD D 12MM L 150MM
1x932200	ROD D 12MM L 200MM
1x932300	ROD D 12MM L 300MM
2x941540	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT D 6MM THREAD 5MM L 180/40 QC
1x1-92080	TRANSFIX PIN 80MM QC D4MM
1x93160	T-WRENCH QC
1x11102	SCREW GUIDE, L 60MM

Foot & Ankle



Galaxy Fixation


Gemini™ System

99-94750 Galaxy Fixation Gemini Metatarsal Set Complete Sterile*,**

Code	Description
2x94010	GALAXY FIXATION GEMINI SINGLE CLAMP
1x932150	ROD D 12MM L 150MM
1x945425	SELFDRILLING XCALIBER CYLINDRICAL SCREW SS L150/25MM D 6/4MM QC
1x1-1300301	DRILL BIT D 3.2MM L 140MM TIN COATED - QC
1x11106	DRILL GUIDE D 3.2MM L 40MM

* Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.

** To be used in combination with ankle kit (99-94730 or 99-94740).



Galaxy Fixation

Gemini™ System

99-94001 Galaxy Fixation Gemini Ankle 6mm Set Multiscrew Clamp Complete Sterile

Code	Description
1x94040	GALAXY FIXATION GEMINI DOUBLE MULTIPIN CLAMP STERILE
4x94010	GALAXY FIXATION GEMINI SINGLE CLAMP STERILE
2x932400	ROD D 12 MM L 400 MM
1x932200	ROD D 12 MM L 200 MM
3x941540	SELF DRILLING XCALIBER, CYLINDRICAL SCREW, SHAFT D.6MM THREAD 5MM, L.180/40 QC
1x1-93050	TRANSFIX PIN 50 MM
2x11137	SCREW GUIDE, LENGTH 80 MM
1x93160	T WRENCH QC

Foot & Ankle




Galaxy Fixation

Gemini™ System

99-94004 Galaxy Fixation Gemini Foot Support Set Complete Sterile

Code	Description
1x932030	SEMI CIRCULAR ROD D12MM LARGE
2x94010	GALAXY FIXATION GEMINI SINGLE CLAMP STERILE



TrueLok™ EVO


Ring Fixation System

99-8800X TL-EVO Footplate Frame Set Sterile*

Code	Description
1x99-882XXX	RX 5/8 MODULAR RING DXXXMM
1x99-883XXX	RX FOOTPLATE DXXXMM
1x99-884XXX	RX FOOT ARCH DXXXMM
1x99-881XXX	RX 3/8 MODULAR RING DXXXMM
3x99-886005	TRUELOK EVO RX STRUT LONG WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
2x99-885003	TLEVO BOLT L 16.5MM
4x99-885000	TLEVO NUT WITH WASHER

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-EVO operative technique TL-2113-OPT

Foot & Ankle



TL-HEX[®] Trauma
TrueLok Hexapod System

99-8801X TL-HEX Footplate Frame Set Sterile*

Code	Description
1x99-56-214X0	TL-HEX 5/8 MODULAR RING DXXXMM
1x99-56-220X0	TL-HEX FOOTPLATE DXXXMM
1x99-56-116X0	TL HALF RING DXXXMM
1x99-56-230X0	TL-HEX 3/8 MODULAR RING DXXXMM STERILE
3x99-886005	TRUELOK EVO RX STRUT LONG WITH DYNA
3x99-885008	TLEVO SPEED NUT
3x99-885009	TLEVO SPEED BOLT
4x99-885003	TLEVO BOLT L16.5MM
2x99-885000	TLEVO NUT WITH WASHER

* For all the specific reference codes with related ring and strut sizes, please refer to the TL-HEX operative technique TH-1905-OPT




XCaliber[™]
Articulated Ankle Fixator

99-91647UE XCaliber Ankle Complete Sterile Kit

Code	Description
1x91150	UNIVERSAL T WRENCH
1x91017	UNIVERSAL ALLEN WRENCH
1x11106	DRILL GUIDE 3.2X40MM
1x11104	DRILL GUIDE 4.8X40MM
1x91047	XCALIBER ANKLE FIXATOR
1x11947	RADIOLUCENT ANKLE PIN GUIDE
3x11102	SCREW GUIDE L 60MM
1x91015	DYNAMIC COMPRESSION/DISTRACTION UNIT
2x99-611540	XCALIBER OSTEOTITE SCREW 150/40
3x99-611530	XCALIBER OSTEOTITE SCREWS 150/30
1x1101101	CANNULATED DRILL BIT 3.2X200MM
2x11014	X-WIRES WITHOUT OLIVE 1.5X250MM
1x1100101	DRILL BIT 4.8X180MM


Foot & Ankle



CalcFix Plus
Calcaneal Fixator System

99-M1460 CalcFix Plus Calcaneal Fixator Set Sterile

Code	Description
2x10012	ALLEN WRENCH, 3MM
1xM210	T-WRENCH FOR BONE SCREWS
1xM211	DRIVER FOR 3MM SHAFT SCREWS
1xM146	CALCFIX PLUS CALCANEAL FIXATOR
2xM315	SELF-DRILLING, CORTICAL SCREW, LENGTH 70/25MM D 3.0-2.5MM
6xM317	SELF-DRILLING, CORTICAL SCREW, LENGTH 100/30MM D 3.0-2.5MM
3x1371501	KIRSCHNER WIRE LENGTH 150MM DIAMETER 1.5MM




Pennig
Minifixator

99-MIL10 Mini-Pennig Fixator Std Kit Complete Sterile


Qty	Description
1x	ALLEN WRENCH, 3 MM
1x	THREADED, WIRE EXTRACTOR
6x	PMF THREADED WIRES D 1.6 MM L 70/15 MM
1x	PENNIG MINIFIXATOR, STD, BODY
1x	CLAMPS, STANDARD, (ONE PAIR) FOR PENNIG MINIFIXATOR
1x	NUT, COMPRESSION/DISTRACTION
2x	PENNIG MINIFIXATOR SPACER
1x	LENGTHENING BAR, STD, 100 MM FOR PENNIG MINIFIXATOR
6x	PMF THREADED WIRES D 2.0 MM L 100/15 MM
1x	REDUCTION FORCEPS (2 PAIRS) FOR PENNIG MINIFIXATOR

Pediatrics



99-93520 Galaxy Medium - Pediatric Sterile Kit - Ø5mm Screw Thread


Code	Description
2x93140	DOUBLE MULTIPIN CLAMP MEDIUM
2x939250	ROD Ø 9MM L250MM
4x944540	SELFDILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 6MM THREAD 5MM L150/40QC
2x11137	SCREW GUIDE 80MM
1x30017	ALLEN WRENCH 5MM
1x93160	QC WRENCH



99-93521 Galaxy Medium - Pediatric Sterile Kit - Ø4mm Screw Thread


Code	Description
2x93140	DOUBLE MULTIPIN CLAMP MEDIUM
2x939250	ROD Ø 9MM L250MM
4x944530	SELFDILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 6MM THREAD 4MM L150/30QC
2x11137	SCREW GUIDE 80MM
1x30017	ALLEN WRENCH 5MM
1x93160	QC WRENCH

Pediatrics



99-93520 Galaxy Medium - Pediatric Sterile Kit - Ø5mm Screw Thread


Code	Description
2x93140	DOUBLE MULTIPIN CLAMP MEDIUM
2x939250	ROD Ø 9MM L250MM
4x944540	SELFDILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 6MM THREAD 5MM L150/40QC
2x11137	SCREW GUIDE 80MM
1x30017	ALLEN WRENCH 5MM
1x93160	QC WRENCH



99-93521 Galaxy Medium - Pediatric Sterile Kit - Ø4mm Screw Thread

Code	Description
2x93140	DOUBLE MULTIPIN CLAMP MEDIUM
2x939250	ROD Ø 9MM L250MM
4x944530	SELFDILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 6MM THREAD 4MM L150/30QC
2x11137	SCREW GUIDE 80MM
1x30017	ALLEN WRENCH 5MM
1x93160	QC WRENCH

Pediatrics



Galaxy Fixation®
System

99-93498 Galaxy Small Z-Configuration Sterile Kit

Code	Description
6x93310	SMALL CLAMP
1x936120	ROD Ø 6MM L 120MM
1x936100	ROD Ø 6MM L 100MM
1x936080	ROD Ø 6MM L 80MM
2x19995	SLEEVE 4.5MM
2x13530	DRILL GUIDE D 2.7MM
1x91017	ALLEN WRENCH
4x947320	SELF DRILLING XCALIBER CYLINDRICAL SCREW SHAFT Ø 4MM THREAD 3MM L 100/20 QC
2x1-1355001	DRILL BIT D 2,7MM
1x93160	"T" WRENCH AO QC



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

Rx Only

Distributed by:

Orthofix.com

OF-1602-PL-E0 H 09/23



 **ORTHOFIX[®]**
ORTHOPEDICS