

# Oscar Pro™

System

## GENERATOR

- 7-inch touch screen
- Multi-language
- Audible feedback
- Visual notifications
- Integrated Wi-Fi
- Data log

## ULTRASONIC ARTHROPLASTY REVISION INSTRUMENT

Prominence

Proactivity

Professionalism

Progress

# UNIVERSAL HANDSET



Magnetic button with LED light


Universal handset:

- Cement removal
- Cementless prosthesis removal
- Osteotomies




Accommodate connection of irrigation nozzle


Aluminum



Optional Handset handle for enhanced maneuverability and grip



Handset holder to stabilise universal handset and protect probe while not in use




## FOCUSED INSTRUMENTATION


- New tray design allowing an easy probe assembly inside the sterile field
- Labelled instrumentation for easy identification of corresponding parts

## SINGLE USE PROBES

- Extended range of probes
- No extension bar needed
- Single use
- M6 connection (reduces the probability of breaking)
- Easy probe coupling



Distributed by:



Instructions for Use: See actual package insert for Instructions for Use.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for specific information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.

Manufactured by:  
**ORTHOFIX Srl**  
 Via Delle Nazioni 9, 37012 Bussolengo  
 (Verona), Italy  
 Telephone +39 045 6719000,  
 Fax +39 045 6719380  
 www.orthofix.com



orthofix.com  
 www.oscar-ultrasonic.com

05-2004-PL-E0 A 07/22

