

Oscar Pro™

System

Ultrasonic Arthroplasty Revision Instrument

Prominence

Proactivity

Professionalism

Progress



GENERATOR

- 7-inch touch screen
- Multi-language
- Audible feedback
- Visual notifications
- Integrated Wi-Fi
- Data log

Universal Handset



Universal handset:


- Cement removal
- Cementless prosthesis removal
- Osteotomies

Accommodate connection of irrigation nozzle



Magnetic button with LED light

Aluminum



Optional Handset handle for enhanced maneuverability and grip




Handset holder to stabilise universal handset and protect probe while not in use



Focused Instrumentation

- New tray design allowing an easy probe assembly inside the sterile field
- Labelled instrumentation for easy identification of corresponding parts



Single Use Probes

- Extended range of probes
- No extension bar needed
- Single use
- M6 connection (reduces the probability of breaking)
- Easy probe coupling

O4SL09D09

O4FL08B06

O4PL09D06

O4EL09B06

O4SCL08B04

O4AL09DNA

O4PSFB8


O4EXP

O4OL09B10

O4HL18D06

O4CL09B10

Distributed by:



Instructions for Use: See actual package insert for Instructions for Use.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for specific information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9,
37012 Bussolengo (Verona), Italy
Tel.: +39 045 6719000
Fax: +39 045 6719380

CE 0123



orthofix.com
www.oscar-ultrasonic.com

