




PENNIG MINIFIXATOR



# Pennig™ Minifixator



ORTHOFIX®

# For fractures, joint fusions, joint distractions, lengthening and angular corrections



**Secure fracture fixation**  
with adjacent joints free to move

**Minimally invasive**  
percutaneous wire insertion

**Straightforward technique**

- double ball-joint makes fracture reduction possible in all planes
- reduction forceps facilitate fracture alignment with the fixator in place

**Versatile**

- clamp design allows for choice of wire placement within same clamp (convergent or parallel)
- straight bars designed for lengthening
- the system allows precise, controlled compression / distraction




**MINI-PENNIG FIXATOR STD KIT COMPLETE STERILE 99-MIL10**

Description	Qty
ALLEN WRENCH, 3 MM	1
THREADED, WIRE EXTRACTOR	1
PMF THREADED WIRES D 1.6 MM L.70/15 MM	6
PENNIG MINIFIXATOR, STD, BODY	1
CLAMPS, STANDARD, (ONE PAIR) FOR PENNIG MINIFIXATOR	1
NUT, COMPRESSION/DISTRACTION	1
PENNIG MINIFIXATOR SPACER	2
LENGTHENING BAR, STD, 100 MM FOR PENNIG MINIFIXATOR	1
PMF THREADED WIRES D 2.0 MM L.100/15 MM	6
REDUCTION FORCEPS (2 PAIRS) FOR PENNIG MINIFIXATOR	1

For further ordering information see section "Pennig Minifixator" in the Orthofix General Catalogue. See Quick Reference Guide PG 090

References:

- 1) Pennig D. et al.: "The Use of Minimally Invasive Fixation in Fractures of The Hand - The Minifixator Concept". Injury, 31 (2000) Suppl. 1: 102-112.
- 2) Pennig D. et Gausepohl T.: "Metacarpal Fractures, Phalangeal Fractures and Reconstructive Procedures: the Pennig MiniFixator in the Hand" Orthofix External Fixation in Trauma and Orthopaedics, G. De Bastiani et al. Springer 2000, Ch. 19: 195-218.

  
Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380

 0123

[www.orthofix.com](http://www.orthofix.com)

PB MIN EO G 05/22

Distributed by:

