


Prefix²

Fixator

The new beat of a
multifunctional fixator

Temporary, spanning fixation




Fracture stability assured

Optimal stability due to large rod and screw diameter
Excellent torsional stability

Screw Clip Clamp

- > Free pin placement for rapid assembly of desired frame construction
- > Easy, fast connection of 6 mm screws
- > Transfixing pin clamp also available for use in ankle-spanning configurations



Rod Clip Clamp

- > Quick, stable connection of two rods in all planes
- > Strong radiolucent rod in four different lengths



MRI
conditional
up to 1.5 Tesla*

* See instruction leaflet (PQ PFX) and Prefix II Operative Technique.

Prefix²
Fixator

Multi-Screw Clamp

- > Parallel screw positioning in T- or straight clamp configuration for later exchange with Orthofix External Fixators
- > Four possible connections for best placement of Rod Clip Clamp


Fast and easy stabilisation
Facilitates damage control of soft tissues

Rapid,
versatile
assembly

99-92501 PREFIX² DIAPHYSEAL COMPLETE KIT STERILE SHORT SCREWS

Consisting of:


Part #	Description
2X92125	ROD 125 MM LONG
1X92275	ROD 275 MM LONG
2X92012	BLUE ROD CLIP CLAMP
4X92016	BLUE SCREW CLIP CLAMP
1X11137	SCREW GUIDE, 80 MM
1X11138	DRILL BIT GUIDE, Ø 4.8 MM, LENGTH 60 MM
1X1100101	DRILL BIT 4.8 MM
1X91150	UNIVERSAL T-WRENCH
2X911560	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 60 MM
2X911540	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 40 MM



99-92502 PREFIX² DIAPHYSEAL COMPLETE KIT STERILE NO SCREWS

Consisting of:


Part #	Description
2X92125	ROD 125 MM LONG
1X92275	ROD 275 MM LONG
2X92012	BLUE ROD CLIP CLAMP
4X92016	BLUE SCREW CLIP CLAMP
1X11137	SCREW GUIDE, 80 MM
1X11138	DRILL BIT GUIDE, Ø 4.8 MM, LENGTH 60 MM
1X1100101	DRILL BIT 4.8 MM
1X91150	UNIVERSAL T-WRENCH
2X911560	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 60 MM
2X911540	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 40 MM



99-92503 PREFIX² DIAPHYSEAL COMPLETE KIT STERILE LONG SCREWS

Consisting of:

Part #	Description
2X92125	ROD 125 MM LONG
1X92275	ROD 275 MM LONG
2X92012	BLUE ROD CLIP CLAMP
4X92016	BLUE SCREW CLIP CLAMP
1X11137	SCREW GUIDE, 80 MM
1X11138	DRILL BIT GUIDE, Ø 4.8 MM, LENGTH 60 MM
1X1100101	DRILL BIT 4.8 MM
1X91150	UNIVERSAL T-WRENCH
2X911550	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 50 MM
2X912650	SELF-DRILLING XCALIBER SCREW, L. 260 MM, THREAD LENGTH 50 MM



Tibia



Femur



Humerus




Bibliography

- Giotakis N., Narayan B. Stability with unilateral external fixation in the tibia. Strat Traum Limb Recon (2007) 2:13–20.
- Behrens F, Johnson W. Unilateral External Fixation Methods to Increase and Reduce Frame Stiffness. Clin Orthop Relat Res 1989;241:48–56.

99-92504 PREFIX² ANKLE COMPLETE KIT STERILE

Consisting of:


Part #	Description
1x92125	ROD 125 MM LONG
1x92275	ROD 275 MM LONG
1x92175	ROD 175 MM LONG
2x92012	BLUE ROD CLIP CLAMP
2x92016	BLUE SCREW CLIP CLAMP
2x92026	BLUE TRANSFIXING PIN CLAMP
1x11137	SCREW GUIDE, 80 MM
1x11138	DRILL BIT GUIDE, Ø 4.8 MM, LENGTH 60 MM
1x1100101	DRILL BIT 4.8 MM
1x91150	UNIVERSAL T-WRENCH
2x911540	SELF-DRILLING XCALIBER SCREW, L. 150 MM, THREAD LENGTH 40 MM
1x92080	TRANSFIXING PIN, THREAD LENGTH 80 MM, THREAD Ø 5 MM, SHAFT Ø 4 MM




92995 PREFIX² STERILIZATION BOX, EMPTY

Can accommodate:


Part #	Description
2x92125	ROD 125 MM LONG
2x92175	ROD 175 MM LONG
2x92275	ROD 275 MM LONG
2x92350	ROD 350 MM LONG
2x19930	MULTI SCREW CLAMP GUIDE
2x11138	DRILL BIT GUIDE, Ø 4.8 MM, LENGTH 60 MM
1X30017	ALLEN WRENCH 5 MM
1x91150	UNIVERSAL T-WRENCH
4x92012	BLUE ROD CLIP CLAMP
8x92016	BLUE SCREW CLIP CLAMP
2x92026	BLUE TRANSFIXING PIN CLAMP
2x1100201	DRILL BIT, Ø 4.8 MM, LENGTH 240 MM
2x11137	SCREW GUIDE, 80 MM
1x92080	TRANSFIXING PIN, THREAD LENGTH 80 MM, THREAD Ø 5 MM, SHAFT Ø 4 MM
1x92050	TRANSFIXING PIN, THREAD LENGTH 50 MM, THREAD Ø 5 MM, SHAFT Ø 4 MM
2x92031	BLUE MULTI SCREW CLAMP
1x92017	PREFIX ² CLAMP HOLDER



Ankle



Cross Knee Application



Cross Ankle Application



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo
(Verona), Italy
Telephone +39 045 6719000,
Fax +39 045 6719380

CE 0123

Rx Only

Orthofix products or services referenced herein are trademarks or registered trademarks of Orthofix SRL and its group of companies. All rights reserved.

www.orthofix.com



PF-0901-PL-E0 C 01/22

Distributed by:



 **ORTHOFIX**[®]
ORTHOPEDICS