

A. OSCAR PRO™



1 OS4000 - OSCAR PRO Generator

See Operative Technique OS-2001-OPT
for application instructions.



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



1. Handset



1
O4HAND
OSCAR PRO UNIVERSAL HANDSET



2. SHOULDER - Cement Removal



1



2



3



4

SINGLE USE PROBES

1

99-04SL18D06

OSCAR PRO TITANIUM CEMENT SCRAPER
L180MM D6MM STERILE

2

99-04PL18D06

OSCAR PRO TITANIUM CEMENT PIERCER
L180MM D6MM STERILE

3


99-04HL18D06

OSCAR PRO TITANIUM CEMENT HOE L180MM
D6MM STERILE

4

99-04HL18D08

OSCAR PRO TITANIUM CEMENT HOE L180MM
D8MM STERILE



OSCAR PRO

3. ELBOW - Cement Removal



1



2



3



4

SINGLE USE PROBES

1

99-04SL18D06

OSCAR PRO TITANIUM CEMENT SCRAPER
L180MM D6MM STERILE

2

99-04PL18D06

OSCAR PRO TITANIUM CEMENT PIERCER
L180MM D6MM STERILE

3


99-04HL18D06

OSCAR PRO TITANIUM CEMENT HOE L180MM
D6MM STERILE

4


99-04HL18D08

OSCAR PRO TITANIUM CEMENT HOE L180MM
D8MM STERILE



OSCAR PRO

4. HIP



SINGLE USE PROBES

Cement Removal - SCRAPER

- 1**
99-04SL09D06
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L90MM D6MM STERILE
- 99-04SL18D06**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L180MM D6MM STERILE
- 99-04SL09D08**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L90MM D8MM STERILE
- 99-04SL18D08**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L180MM D8MM STERILE
- 99-04SL27D08**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L270MM D8MM STERILE
- 99-04SL09D10**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L90MM D10MM STERILE
- 99-04SL18D10**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L180MM D10MM STERILE
- 99-04SL27D10**
 OSCAR PRO TITANIUM CEMENT SCRAPER
 L270MM D10MM STERILE

Cement Removal - PIERCER

- 2**
99-04PL09D06
 OSCAR PRO TITANIUM CEMENT PIERCER
 L90MM D6MM STERILE

99-04PL18D06

OSCAR PRO TITANIUM CEMENT PIERCER
 L180MM D6MM STERILE

99-04PL09D08

OSCAR PRO TITANIUM CEMENT PIERCER
 L90MM D8MM STERILE

99-04PL18D08

OSCAR PRO TITANIUM CEMENT PIERCER
 L180MM D8MM STERILE

99-04PL27D08

OSCAR PRO TITANIUM CEMENT PIERCER
 L270MM D8MM STERILE

99-04PL09D10

OSCAR PRO TITANIUM CEMENT PIERCER
 L90MM D10MM STERILE

99-04PL18D10

OSCAR PRO TITANIUM CEMENT PIERCER
 L180MM D10MM STERILE

99-04PL27D10

OSCAR PRO TITANIUM CEMENT PIERCER
 L270MM D10MM STERILE

11

99-04EXPL27

OSCAR PRO TITANIUM CEMENT PLUG
 EXTRACTOR L270MM STERILE

Cement Removal - HOE

3

99-04HL18D06

OSCAR PRO TITANIUM CEMENT HOE L180MM
 D6MM STERILE

99-04HL18D08

OSCAR PRO TITANIUM CEMENT HOE L180MM
 D8MM STERILE

Cementless Prosthesis Removal

4

99-04AL09DNA

OSCAR PRO TITANIUM CURVED CEMENT
 PIERCER L90MM STERILE

5

99-04PSFB8

OSCAR PRO TITANIUM FLAT CEMENT PIERCER
 L90MM STERILE

6

99-04OL09B10

OSCAR PRO TITANIUM FLAT BLADE SMOOTH
 L90MM W10MM STERILE

7

99-04CL09B10

OSCAR PRO TITANIUM CURVED BLADE
 SMOOTH L90MM W10MM STERILE

8

99-04FL08B06

OSCAR PRO TITANIUM FLAT BLADE SMOOTH
 L80MM W6MM STERILE

9

99-04SCL08B04


OSCAR PRO TITANIUM CURVED BLADE
 SERRATED L80MM W4MM STERILE

Osteotomies

10

99-04EL09B06

OSCAR PRO TITANIUM OSTEOTOME WEDGE TIP
 L90MM W6MM STERILE



5. KNEE



1



3



4



6



8



2



5



7

SINGLE USE PROBES

Cement Removal - SCRAPER

1

99-04SL09D06

OSCAR PRO TITANIUM CEMENT SCRAPER
L90MM D6MM STERILE

99-04SL18D06

OSCAR PRO TITANIUM CEMENT SCRAPER
L180MM D6MM STERILE

99-04SL09D08

OSCAR PRO TITANIUM CEMENT SCRAPER
L90MM D8MM STERILE

99-04SL18D08

OSCAR PRO TITANIUM CEMENT SCRAPER
L180MM D8MM STERILE

99-04SL27D08

OSCAR PRO TITANIUM CEMENT SCRAPER
L270MM D8MM STERILE

99-04SL09D10

OSCAR PRO TITANIUM CEMENT SCRAPER
L90MM D10MM STERILE

99-04SL18D10

OSCAR PRO TITANIUM CEMENT SCRAPER
L180MM D10MM STERILE

99-04SL27D10

OSCAR PRO TITANIUM CEMENT SCRAPER
L270MM D10MM STERILE

Cement Removal - PIERCER

2

99-04PL09D06

OSCAR PRO TITANIUM CEMENT PIERCER
L90MM D6MM STERILE

99-04PL18D06

OSCAR PRO TITANIUM CEMENT PIERCER
L180MM D6MM STERILE

99-04PL09D08

OSCAR PRO TITANIUM CEMENT PIERCER
L90MM D8MM STERILE

99-04PL18D08

OSCAR PRO TITANIUM CEMENT PIERCER
L180MM D8MM STERILE

99-04PL27D08

OSCAR PRO TITANIUM CEMENT PIERCER
L270MM D8MM STERILE

99-04PL09D10

OSCAR PRO TITANIUM CEMENT PIERCER
L90MM D10MM STERILE

99-04PL18D10

OSCAR PRO TITANIUM CEMENT PIERCER
L180MM D10MM STERILE

99-04PL27D10

OSCAR PRO TITANIUM CEMENT PIERCER
L270MM D10MM STERILE

3

99-04PSFB8

OSCAR PRO TITANIUM FLAT CEMENT PIERCER
L90MM STERILE

Cementless Prosthesis Removal

4

99-04OL09B10

OSCAR PRO TITANIUM FLAT BLADE SMOOTH
L90MM W10MM STERILE

5

99-04CL09B10

OSCAR PRO TITANIUM CURVED BLADE
SMOOTH L90MM W10MM STERILE

6

99-04FL08B06

OSCAR PRO TITANIUM FLAT BLADE SMOOTH
L80MM W6MM STERILE

7

99-04SCL08B04


OSCAR PRO TITANIUM CURVED BLADE
SERRATED L80MM W4MM STERILE

Osteotomies

8



99-04EL09B06

OSCAR PRO TITANIUM OSTEOTOME WEDGE TIP
L90MM W6MM STERILE



OSCAR PRO

6. Others



7. Sales Configuration

SINGLE USE STERILE ITEMS

The following items codes and description of probes are provided in single, sterile packaging configuration. The packaging consists of double pouch inner packaging with a cardboard FSC outer box.

Part #	Description - Single use probes
99-04SL09D06	OSCAR PRO TITANIUM CEMENT SCRAPER L90MM D6MM STERILE
99-04SL18D06	OSCAR PRO TITANIUM CEMENT SCRAPER L180MM D6MM STERILE
99-04SL09D08	OSCAR PRO TITANIUM CEMENT SCRAPER L90MM D8MM STERILE
99-04SL18D08	OSCAR PRO TITANIUM CEMENT SCRAPER L180MM D8MM STERILE
99-04SL27D08	OSCAR PRO TITANIUM CEMENT SCRAPER L270MM D8MM STERILE
99-04SL09D10	OSCAR PRO TITANIUM CEMENT SCRAPER L90MM D10MM STERILE
99-04SL18D10	OSCAR PRO TITANIUM CEMENT SCRAPER L180MM D10MM STERILE
99-04PL09D06	OSCAR PRO TITANIUM CEMENT PIERCER L90MM D6MM STERILE
99-04PL18D06	OSCAR PRO TITANIUM CEMENT PIERCER L180MM D6MM STERILE
99-04PL09D08	OSCAR PRO TITANIUM CEMENT PIERCER L90MM D8MM STERILE
99-04PL18D08	OSCAR PRO TITANIUM CEMENT PIERCER L180MM D8MM STERILE
99-04PL27D08	OSCAR PRO TITANIUM CEMENT PIERCER L270MM D8MM STERILE
99-04PL09D10	OSCAR PRO TITANIUM CEMENT PIERCER L90MM D10MM STERILE
99-04PL18D10	OSCAR PRO TITANIUM CEMENT PIERCER L180MM D10MM STERILE
99-04HL18D06	OSCAR PRO TITANIUM CEMENT HOE L180MM D6MM STERILE
99-04HL18D08	OSCAR PRO TITANIUM CEMENT HOE L180MM D8MM STERILE
99-04EXP	OSCAR PRO TITANIUM CEMENT PLUG EXTRACTOR L180MM STERILE
99-04AL09DNA	OSCAR PRO TITANIUM CURVED CEMENT PIERCER L90MM STERILE
99-04PSFB8	OSCAR PRO TITANIUM FLAT CEMENT PIERCER L90MM STERILE
99-04OL09B10	OSCAR PRO TITANIUM FLAT BLADE SMOOTH L90MM W10MM STERILE
99-04CL09B10	OSCAR PRO TITANIUM CURVED BLADE SMOOTH L90MM W10MM STERILE
99-04FL08B06	OSCAR PRO TITANIUM FLAT BLADE SMOOTH L80MM W6MM STERILE
99-04SCL08B04	OSCAR PRO TITANIUM CURVED BLADE SERRATED L80MM W4MM STERILE
99-04EL09B06	OSCAR PRO TITANIUM OSTEOTOME WEDGE TIP L90MM W6MM STERILE
99-04PL27D10	OSCAR PRO TITANIUM CEMENT PIERCER L 270MM D 10MM STERILE
99-04SL27D10	OSCAR PRO TITANIUM CEMENT SCRAPER L 270MM D 10MM STERILE
99-04EXPL27	OSCAR PRO TITANIUM CEMENT PLUG EXTRACTOR L 270MM STERILE

IRRIGATION KIT

Part #	Description
2000972	STERILE KIT FOR PERISTALTIC PUMP
8000001	STERILE KIT FOR PERISTALTIC PUMP US

The above Part # (not for sale in Canada) are manufactured by:

MANUFACTURER BY:



EMOTEC S.R.L. - VIA MAESTRI DEL LAVORO, 5 - 41036 MEDOLLA (MO) ITALY
C.F. E P. IVA 02578710366 - TEL. : +39 0535 / 53060 - FAX: +39 0535 / 52629 -
E-MAIL: EMOTEC@EMOTEC.IT

DISTRIBUTED BY:



Orthofix Srl - Via delle Nazioni 9 -
37012 Bussolengo (VR) Italy
Telephone +39 045 671 9000 - Fax +39 045 671 9380

The Irrigation Kit is available in U.S. and Europe. For the availability in other countries, please contact the manufacturer.

* Component of the device: OSCARPRO – ULTRASOUND MEDICAL SYSTEM

MULTIPLE USE, NON STERILE ITEMS

The following items are sold individually packed.

Part #	Description
O4HAND*	OSCAR PRO UNIVERSAL HANDSET
O4UHH	OSCAR PRO UNIVERSAL HANDSET HANDLE
O4HOLD	OSCAR PRO HANDSET HOLDER
O4IPL200	OSCAR PRO SLAP HAMMER
SPAN7MM	SPANNER 7 MM
SPAN8/9MM	SPANNER 8/9MM
O4CABLE*	OSCAR PRO UNIVERSAL HANDSET CABLE
O4FOOT-U*	OSCAR PRO FOOTSWITCH FOR ULTRASOUND
O4FOOT-W*	OSCAR PRO FOOTSWITCH FOR IRRIGATION
OSCAR PRO Main lead:	
O4POWEU	OSCAR PRO MAIN LEAD EU (not for sale in Canada)
or	
O4POWUS	OSCAR PRO MAIN LEAD US
or	
O4POWUK	OSCAR PRO MAIN LEAD UK (not for sale in Canada)

A sterilization tray is available to sterilize multiple use.

O4TRAY - OSCAR PRO INSTRUMENT TRAY EMPTY

The tray can accommodate reusable and re-processable instruments, as per the following configuration:

Q.ty	Part #	Description
2	O4UHH	OSCAR PRO UNIVERSAL HANDSET HANDLE
2	O4HOLD	OSCAR PRO HANDSET HOLDER
1	O4IPL200	OSCAR PRO SLAP HAMMER
2	SPAN7MM	SPANNER 7MM
2	SPAN8/9MM	SPANNER 8/9MM
2	O4CABLE	OSCAR PRO UNIVERSAL HANDSET CABLE

OSCAR PRO Instrument Tray spare part:

O4TRAY-01	OSCAR PRO Base tray
O4TRAY-02	OSCAR PRO Lid tray

A dedicated flight case for the instrument tray is available and sold separately; the flight case includes a custom foam in a non-sterile configuration.

OINSTFLY1 - OSCAR INSTRUMENT FLIGHT CASE

The OSCAR PRO generator is provided within a dedicated flight case, including also one water bag holder.

OS4000 - OSCAR PRO GENERATOR

Q.ty	Part #	Description
1		OSCAR PRO GENERATOR
1	OGENFLY2	OSCAR PRO GENERATOR FLIGHT CASE
1	O4BEAM	OSCAR PRO HOLDER FOR WATER BAG L459MM


The OGENFLY2 flight case is designed to accommodate also the following items, as per below configuration:

OGENFLY2 - OSCAR PRO GENERATOR FLIGHT CASE

Q.ty	Part #	Description
1		OSCAR PRO main lead:
	O4POWEU	OSCAR PRO MAIN LEAD EU (not for sale in Canada)
	O4POWUS	OSCAR PRO MAIN LEAD US
	O4POWUK	OSCAR PRO MAIN LEAD UK (not for sale in Canada)
1	O4FOOT-U	OSCAR PRO FOOTSWITCH FOR ULTRASOUND
1	O4FOOT-W	OSCAR PRO FOOTSWITCH FOR IRRIGATION

OSCAR PRO Generator spare parts:

Part #	Description
OGENFLY2	OSCAR PRO GENERATOR FLIGHT CASE
O4BEAM	OSCAR PRO HOLDER FOR WATER BAG L459MM



OSCAR PRO