

K. Contours Proximal Humeral Plate (PHP) + LSP Titanium Humeral Plate



1

1 Contours PHP


See Operative Technique PH-1201-OPT
"Contours PHP Proximal Humeral Plate".



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



1. Contours PHP Components



1
Contours PHP (Sterile)
99-TP3001 Left
99-TP3000 Right


2
Diaphyseal Screw \varnothing 4.5mm (Sterile)
99-TP1422 Length 22mm
99-TP1424 Length 24mm
99-TP1426 Length 26mm
99-TP1428 Length 28mm
99-TP1430 Length 30mm
99-TP1432 Length 32mm
99-TP1434 Length 34mm
99-TP1436 Length 36mm

3
Revision Screw \varnothing 6.5mm (Sterile)
99-TP1526 Length 26mm
99-TP1528 Length 28mm
99-TP1530 Length 30mm
99-TP1532 Length 32mm
99-TP1534 Length 34mm

4


Main Locking Screw \varnothing 6.0mm (Sterile)
99-TP1640 Length 40mm
99-TP1645 Length 45mm
99-TP1650 Length 50mm
99-TP1655 Length 55mm
99-TP1660 Length 60mm

5
Fine Threaded Screw \varnothing 2.2mm, pack of 2 (Sterile)
99-TP1720 Length 120/20mm
99-TP1725 Length 120/25mm
99-TP1730 Length 120/30mm
99-TP1735 Length 120/35mm
99-TP1740 Length 120/40mm
99-TP1745 Length 120/45mm
99-TP1750 Length 120/50mm
99-TP1755 Length 120/55mm



2. Contours PHP Instrumentation (Sterile)

(Out of Tray)




1
99-TP1111K-Wire (Ø 2mm) L. 200mm
(pack of 4)

2
99-1-TP3107 Cannulated Drill Bit (Ø 4.8mm)
quick connect (sterile)

3
99-1-TP1104 Drill Bit (Ø 3.2mm) quick
connect (sterile)



3. Sterilization Box Contours PHP



1

1 TP3995 Contours PHP Sterilization Box, empty

Can accommodate:

- 1xTP3106 Quick connect Screwdriver for hexagonal head
- 1xTP1145 Ruler
- 1xTP3102 3.5mm quick connect Screwdriver
- 1x1-TP3107 Cannulated Drill Bit (\varnothing 4.8mm) quick connect
- 2x1-TP1104 Drill Bit (\varnothing 3.2mm) quick connect
- 1xTP3150 Handle (soft touch) quick connect
- 1xTP3108 Plate Holder
- 1xTP1100 Drill Guide (\varnothing 3.2mm)
- 1xTP1101 Locking T-Wrench

