


D. Prefix2 Fixator



1

1
PreFix² Fixator

See Operative Technique PF-0902-OPT
"The PreFix² Fixator" for application instructions.




CE 0123

Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380



1. Components/Replacements



1



2



3



4



5



6



7

1
Rod
92125 Length 125mm
92175 Length 175mm
92275 Length 275mm
92350 Length 350mm

2
92012 Blue Rod Clip Clamp

3
92016 Blue Screw Clip Clamp


4
92026 Blue Transfixing Pin Clamp

5
92031 Blue Multi Screw Clamp

6
92080 Transfixing Pin
Thread length 80mm,
thread Ø 5mm, shaft Ø 4mm.

7
92050 Transfixing Pin
Thread length 50mm,
thread Ø 5mm, shaft Ø 4mm.

Note: All Transfixing Pins are also available packaged sterile. They can be ordered using the above code numbers preceded by 99- (e.g. 99-92080).




2. Template and Instrumentation



1



2




3

1
19930 Multi Screw Clamp Guide


2
30017 Allen Wrench 5mm

3
92017 Prefix² Clamp Holder


Note: See also instrumentation for bone screw insertion (section "Bone Screws").



3. Kits



1



2

1 99-92501 Prefix² Diaphyseal Complete Kit Sterile Short Screws


Consisting of:

- 2x92125 Rod 125mm long
- 1x92275 Rod 275mm long
- 2x92012 Blue Rod Clip Clamp
- 4x92016 Blue Screw Clip Clamp
- 1x11137 Screw Guide, 80mm
- 1x11138 Drill Bit Guide, Ø 4.8mm, length 60mm
- 1x1100101 Drill Bit 4.8mm, length 180mm
- 1x91150 Universal T-Wrench
- 2x911560 Self-Drilling XCaliber Screw, L. 150mm, thread length 60mm
- 2x911540 Self-Drilling XCaliber Screw, L. 150mm, thread length 40mm


2 99-92502 Prefix² Diaphyseal Complete Kit Sterile No Screws

Consisting of:

- 2x92125 Rod 125mm long
- 1x92275 Rod 275mm long
- 2x92012 Blue Rod Clip Clamp
- 4x92016 Blue Screw Clip Clamp
- 1x11137 Screw Guide, 80mm
- 1x11138 Drill Bit Guide, Ø 4.8mm, length 60mm
- 1x1100101 Drill Bit 4.8mm, length 180mm
- 1x91150 Universal T-Wrench



Kits




3

3 99-92503 Prefix² Diaphyseal Complete Kit Sterile Long Screws

Consisting of:

- 2x92125 Rod 125mm long
- 1x92275 Rod 275mm long
- 2x92012 Blue Rod Clip Clamp
- 4x92016 Blue Screw Clip Clamp
- 1x11137 Screw Guide, 80mm
- 1x11138 Drill Bit Guide, Ø 4.8mm, length 60mm
- 1x1100101 Drill Bit 4.8mm, length 180mm
- 1x91150 Universal T-Wrench
- 2x911550 Self-Drilling XCaliber Screw, L. 150mm, thread length 50mm
- 2x912650 Self-Drilling XCaliber Screw, L. 260mm, thread length 50mm




4


4 99-92504 Prefix² Ankle Complete Kit Sterile

Consisting of:

- 1x92125 Rod 125mm long
- 1x92275 Rod 275mm long
- 1x92175 Rod 175mm long
- 2x92012 Blue Rod Clip Clamp
- 2x92016 Blue Screw Clip Clamp
- 2x92026 Blue Transfixing Pin Clamp
- 1x11137 Screw Guide, 80mm
- 1x11138 Drill Bit Guide, Ø 4.8mm, length 60mm
- 1x1100101 Drill Bit 4.8mm, length 180mm
- 1x91150 Universal T-Wrench
- 2x911540 Self-Drilling XCaliber Screw, L. 150mm, thread length 40mm
- 1x92080 Transfixing Pin, thread length 80mm, thread Ø 5mm, shaft Ø 4mm



4. Sterilization Box



1

1 92995 Prefix² Sterilization Box, empty

Can accommodate:

- 2x92125 Rod 125mm long
- 2x92175 Rod 175mm long
- 2x92275 Rod 275mm long
- 2x92350 Rod 350mm long
- 2x19930 Multi Screw Clamp Guide
- 2x11138 Drill Bit Guide, Ø 4.8mm, length 60mm
- 1X30017 Allen Wrench 5mm
- 1x91150 Universal T-Wrench
- 4x92012 Blue Rod Clip Clamp
- 8x92016 Blue Screw Clip Clamp
- 2x92026 Blue Transfixing Pin Clamp
- 2x1100201 Drill Bit, Ø 4.8mm, length 240mm
- 2x11137 Screw Guide, 80mm
- 1x92080 Transfixing Pin, thread length 80mm, thread Ø 5mm, shaft Ø 4mm

- 1x92050 Transfixing Pin, thread length 50mm, thread Ø 5mm, shaft Ø 4mm
- 2x92031 Blue Multi Screw Clamp
- 1x92017 Prefix² Clamp Holder

