

U. TrueLok™ Elevate

1



1 TrueLok Elevate Frame

See Operative Technique
"TrueLok Elevate Transverse Bone Transport
System" TL-2403-OPT".



CE 0123

Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

EF-TL-Elevate EO A 06/25



1. TrueLok™ Elevate Components

1



2



3




1
50-2505 TrueLok Elevate External Fixator

2
3/8 Rings
50-2511 Half Pin L120mm
D4.0mm thread D4.0xL18mm

3
54-11250 TL+ Half Pin Quick Connect
Ø 6mm x 180mm



2. TrueLok™ Elevate Instruments



1
M210 T-WRENCH FOR BONE SCREWS

2
50-2510 TL STRAIGHT TUBE WRENCH

3
54-1154 TL, WRENCH, COMBO, 10MM

4
93162 T-WRENCH HEXAGON 5-5 QC

5
54-2226 TELESCOPIC NAIL DRILL BIT
D6.0MM STERILE


6
50-2509 TRUELOK ELEVATE DRILL
POSITIONING GUIDE

7
W1003 WIRE CUTTER


8
54-2227 NEEDLE NOSE PLIERS, STAINLESS
STEEL

9
50-2515 DRILL GUIDE D3.2MM L60MM

10
11138 DRILL GUIDE D4.8MM L60MM




3. TrueLok™ Elevate Set Sterile



TL-Elevate is available in a sterile kit including the items below:

99-50-2513 TRUELOK ELEVATE SET STERILE		
Part Number	Description	Q.ty
1-1355001	Drill Bit D.2,7mm L.127mm TIN	1
50-25071	Drill Insert	2
50-25070	Wire Insert	2
50-2505	TrueLok Elevate External Fixator	1
50-2506	TrueLok Elevate Template	1
50-2511	Half Pin L120mm D4.0mm Thread D4.0xL18mm	2
50-2508	K-Wire L180mm D1.8mm	3




4. TrueLok™ Elevate Instruments Tray Complete



The instruments below are not included in the sterile kit and are necessary to implant TL-Elevate

50-2590C TRUELOK ELEVATE INSTRUMENT TRAY COMPLETE

Part Number	Description	Q.ty
M210	T-Wrench For Bone Screws	1
50-2510	TL Straight Tube Wrench	1
54-1154	TL, Wrench, Combo, 10mm	1
93162	T-Wrench Hexagon 5-5 QC	1
54-2226	TL, 90 Degree Tubular Wrench	1
50-2509	TrueLok Elevate Drill Positioning Guide	1
W1003	Wire Cutter	1
54-2227	Needle Nose Pliers, Stainless Steel	1
50-2515	Drill Guide D3.2mm L60mm	1
11138	Drill Guide D4.8mm L60mm	1



TrueLok™ Elevate