


# C. Contours VPS Volar Plating System



## Contours VPS

1

### Long

VPL0310CE Left

VPR0322CE Right

2

### Standard

VPL0309CE Left

VPR0321CE Right


See Manual CV-0702-OPT "Volar Plating with Anatomically Designed Plate and Fixed-Angle Screws" for application instructions



Manufactured by:  
**ORTHOFIX Srl**  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380



# 1. Instrumentation




1

## 1 VP601CE Contours VPS Steribox, empty

Can accommodate:

- 1xDH0411CE 3.0mm Hex Driver
- 2xDH0412CE Bender, Straight Plate
- 1xDH0413CE Bender, Threaded Plate
- 1xDH0416CE Drill Guide, 2.5mm
- 5xDH0421CE K-Wire, 1.1x150mm
- 1xDH0422CE 1.5mm Hex Key, L Pattern
- 1xDH0426CE Depth Gauge, Distal
- 1xDH0427CE Depth Gauge, Proximal
- 1xDH0437CE 1.5 Hex Driver (Tip Only)
- 1xDH0439CE 1.5 Hex Driver
- 2xDH0432CE Drill Guide, 1.6mm
- 2xDH0433CE Drill Guide, 2.0mm
- 2xDH0434CE Drill Bit, 1.6mm
- 2xDH0435CE Drill Bit, 2.0mm
- 2xDH0436CE Drill Bit, 2.5mm

## VP600CE Contours VPS Steribox, complete



## 2. Screws

1



| Total Length (mm) |            | Total Length (mm) |            |
|-------------------|------------|-------------------|------------|
| 14                | PSC2014TCE | 22                | PSC2022TCE |
| 16                | PSC2016TCE | 24                | PSC2024TCE |
| 18                | PSC2018TCE | 26                | PSC2026TCE |
| 20                | PSC2020TCE | 28                | PSC2028TCE |

2




| Total Length (mm) |             | Total Length (mm) |             |
|-------------------|-------------|-------------------|-------------|
| 12                | PSC24P12TCE | 20                | PSC24P20TCE |
| 14                | PSC24P14TCE | 22                | PSC24P22TCE |
| 16                | PSC24P16TCE | 24                | PSC24P24TCE |
| 18                | PSC24P18TCE |                   |             |

3



| Total Length (mm) |            | Total Length (mm) |            |
|-------------------|------------|-------------------|------------|
| 10                | PSC2410TCE | 18                | PSC2418TCE |
| 12                | PSC2412TCE | 20                | PSC2420TCE |
| 14                | PSC2414TCE | 22                | PSC2422TCE |
| 16                | PSC2416TCE |                   |            |

4



| Total Length (mm) |            | Total Length (mm) |            |
|-------------------|------------|-------------------|------------|
| 14                | PSC2714TCE | 22                | PSC2722TCE |
| 16                | PSC2716TCE | 24                | PSC2724TCE |
| 18                | PSC2718TCE | 26                | PSC2726TCE |
| 20                | PSC2720TCE | 28                | PSC2728TCE |

5



| Total Length (mm) |            | Total Length (mm) |            |
|-------------------|------------|-------------------|------------|
| 12                | PSC3512TCE | 18                | PSC3518TCE |
| 14                | PSC3514TCE | 20                | PSC3520TCE |
| 16                | PSC3516TCE | 22                | PSC3522TCE |

1  
2.0mm Locking MicroThreaded Screws

2  
2.4mm Locking Screws

3  
2.4mm Non-Locking Screws

4  
2.7mm Locking Screws

5  
3.5mm Non Locking Cortical Screws

Note: All screws are also available packaged sterile. They can be ordered using the above code numbers preceded by 99- (e.g. 99-PSC2014TCE)

