



RIVAL™

○ REDUCE

□ VIEW

▷ BITE



BITE

Headed Cannulated and
Headless Compression Screws

Self-tapping thread design

Designed to simplify the surgical technique
leading to shorter surgery time

Self-drilling thread design

Designed to reduce the need for pre-drilling

Variety of thread lengths

Optimal implant selection



HEADED SCREWS



Unique, low profile head

To reduce screw prominence



Cruciform driver

In sizes from 2.0 to 8.0mm

2.0mm / 3.0mm / 4.5mm

Hex-head connection

HEADLESS SCREWS

Screw compression

Compression is provided by the difference between thread pitches of the head and main threads





Variety of thread lengths and diameters

From 2.5 to 7.5mm

Universal instrumentation

Designed

- to be used across all RIVAL™ product lines
- to help with inventory management
- to provide a comprehensive range of instruments



Distributed by:



Manufactured by:
 ORTHOFIX Srl
 Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
 Telephone +39 045 6719000, Fax +39 045 6719380



www.orthofix.com

RIBI-1906-PL-E0 A 10/19

