


# TrueLok™

Ring Fixation System

## New TL-Rapid Struts Provide Quick Solutions

- Quick and easy frame assembly
- Acute adjustment with micro-compression or distraction
- Lockable universal hinge on both ends of strut



| Short Range | Medium Range | Long Range  |
|-------------|--------------|-------------|
| 80 - 107mm  | 108 - 191mm  | 170 - 307mm |




Bolt allows  
for quick frame  
adjustments

Acute adjustment  
with micro-compression or  
distraction

Struts are available in three sizes:

| Length | Part #   |
|--------|----------|
| SHORT  | 50-10170 |
| MEDIUM | 50-10180 |
| LONG   | 50-10190 |



Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9, 37012 Bussolengo  
(Verona), Italy  
Telephone +39 045 6719000,  
Fax +39 045 6719380  
www.orthofix.com



Distributed by:

