

TrueLok™ EVO

Ring Fixation System



Evolution All Around

Simple, Stable
and Versatile

Radiolucent

For a clear visualization
of the bone



MR
Conditional

Tested frames

For additional info please refer to
the Instruction for Use (PQEVO)

Quick Application

Thanks to the pre-assembled
sterile sets

Dynamizable

TL-EVO struts have integrated
dynamization module to be
unlocked when needed*




*See the latest Voice of Literature
about dynamization TL-2109-VOL-E0 A 06/21

READY-TO-USE IN PRE-ASSEMBLED STERILE FRAME SETS

TL-EVO is provided in pre-assembled sterile sets for diaphyseal, peri-articular or foot application

Pre-Assembled Sterile Frame Sets


For double-ring or foot & ankle frames, available in different sizes



Different Types of Connection Element Sterile Sets

With all the necessary fixation components to apply the frame to the bone

All TL-EVO components are also provided in sterile single or multiple pack



Please see Orthofix Catalogue or the TL-EVO Operative Technique for more details about the sales configuration. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.

TrueLok™ Family




TrueLok™

Deformity Correction and Trauma

Simple/multiple deformity correction.

Bone transport and best solution for one plane angular deformity.



TL-HEX Trauma

Trauma and Correction

A circular trauma frame that allows minor and acute corrections.

It can be converted into a software-assisted hexapod system (TL-HEX) to treat more complex and multiplanar deformities.



TL-HEX TrueLok Hexapod System™

Deformity Correction and Trauma

Simple/multiple deformity correction (angular/rotational deformity, limb lengthening and bone transport– best solution for multi-plane deformity)

CHAOS – temporary intraoperative reduction system with Unyco Screws.



TrueLok™ EVO

Trauma and Correction

A full sterile, radiolucent system ready for use with integrated dynamization.

It is available in pre-assembled folded configurations, also combined with TL-HEX Rings.



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo
(Verona), Italy
Telephone +39 045 6719000, Fax +39 045
6719380

CE 0123

Rx Only

www.orthofix.com



TL-2124-PL-E0 B 01/22

Distributed by:



 **ORTHOFIX®**
ORTHOPEDICS