


TrueLok™

Ring Fixation System

Achieve Precise Alignment



**TrueLok™
Phantom Hinge**



Designed to maintain:

- Axis of rotation during joint contracture correction
- Joint mobility in arthrodiasis procedures

The Phantom Hinge represents the latest addition to the TrueLok™ Portfolio by providing a stabilizing constrained hinge solution when using circular fixation – whether in cases of trauma or elective deformity surgery.


Radiolucency

The radiolucent material allows for clear monitoring of joint congruence during correction.



Quick and Easy


The Phantom Hinge is a single-use, sterile component that can be quickly and easily added to a TrueLok™ device.



Two Size Offerings

ANKLE HINGE


The smaller size Phantom Hinge is specifically designed for use with the ankle frame.



Cannulation allows to precisely align with reference wire.

KNEE HINGE

The larger size Phantom Hinge is specifically designed for use with the knee-spanning frame.




TrueLok™ Phantom Hinges

Part #	Description	Notes
99-50-2201M	TL PHANTOM HINGE KNEE PACK OF 2	STERILE, 2 PIECES PER PACK
99-50-2202M	TL PHANTOM HINGE ANKLE PACK OF 2	STERILE, 2 PIECES PER PACK


The Phantom Hinge is compatible with the following TrueLok™ Ring Fixation Systems:



TrueLok™
Ring Fixation System



TL-HEX
TrueLok™ Hexapod System



TrueLok™ EVO
Ring Fixation System*

And is part of an extensive portfolio of TrueLok™ Hinges:



TrueLok™
Phantom Hinge



TrueLok™
Tornado Hinge



TrueLok™
Non Lockable
Universal Joint



TrueLok™
Universal Hinge



TrueLok™
Inline Hinge



TrueLok™
Outboard Hinge

TrueLok™, TL-HEX™, and TrueLok™ EVO are trademarks or registered trademarks exclusively licensed to Orthofix Medical Inc., and its affiliate companies. All rights reserved.




Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

Rx Only


orthofix.com
TL-2204-PL-E0 C 04/24

Distributed by:



 **ORTHOFIX®**
ORTHOPEDICS