



# TrueLok™

Ring Fixation System

## Restore Natural Joint Motion



### TrueLok™ Tornado Hinge



Designed to:

- Provide flexible motion reproducing the joint anatomy
- Maintain joint mobility and the natural anatomic axis



The Tornado Hinge represents the latest addition to the TrueLok™ portfolio as a stabilizing unconstrained hinge solution when using circular fixation to carry out arthrodiastasis or during correction of joint contracture in trauma or elective deformity correction procedures.




Quick and Easy

The Tornado Hinge is a single-use, sterile component that can be quickly and easily added to a TrueLok™ frame.

Radiolucency

The radiolucent material allows for clear monitoring of joint congruence during correction.



Versatile Use

It provides flexible and natural motion reproducing the anatomy of the joint.



No need to use a reference wire to target the joint center of rotation.

The central spring mimics the natural motion of the knee cruciate ligament.




TrueLok™ Tornado Hinges		
Part #	Description	Notes
99-50-2203M	TRUELOK TL TORNADO HINGE PACK OF 2	STERILE, 2 PIECES PER PACK


The Tornado Hinge is compatible with the following TrueLok™ Ring Fixation systems:



TrueLok™  
Ring Fixation System



TL-HEX™  
TrueLok™ Hexapod System



TrueLok™ EVO  
Ring Fixation System\*

And is part of an extensive portfolio of TrueLok™ Hinges:



TrueLok™  
Phantom Hinge



TrueLok™  
Tornado Hinge



TrueLok™  
Non Lockable  
Universal Joint



TrueLok™  
Universal Hinge



TrueLok™  
Inline Hinge



TrueLok™  
Outboard Hinge

TrueLok™, TL-HEX™, and TrueLok™ EVO are trademarks or registered trademarks exclusively licensed to Orthofix Medical Inc., and its affiliate companies. All rights reserved.



**Manufactured by:**

ORTHOFIX Srl

Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy

Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

**Rx Only**

  
**orthofix.com**  
TL-2209-PL-E0 C 07/24

**Distributed by:**



 **ORTHOFIX®**  
ORTHOPEDICS